

Product datasheet for **SC318077**

MBTD1 (NM_017643) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MBTD1 (NM_017643) Human Untagged Clone
Tag:	Tag Free
Symbol:	MBTD1
Synonyms:	SA49P01
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318077 representing NM_017643.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTTTGACGGTTATGATAGCTGCAGTGAGGACACAAGCAGCAGCTCCAGCTCCGAAGAGAGTGAGGAA
GAAGTCGCTCCTTACCTTCTAATCTCCGATTATCAAAAACAATGGCAAGCTACACATACCCAGAT
GGTAAATCTGGCATGGCTACCTGTGAGATGTGTGGGATGGTTGGCGTCCGAGATGCTTTTACTCTAAA
ACAAAGCGTTTCTGTAGCGTTTCATGTTCAAGAAGTTACTCGTCAAACTCCAAGAAGGCAAGCATTTTG
GCCAGACTTCAGGGTAAGCCTCCAACAAAGAAAGCAAAAGTTCTTCAGAAACAACCTTTAGTTGCTAAG
CTAGCCGCATATGCTCAGTATCAAGCTACCTTGCAAAATCAAGCAAGACAAAAGCAGCAGTCTCCATG
GAAGGTTTCAGCTGGGTAAGTACATCAATAGCAATAGCTTTATAGCAGCTCCGTTACCTGTTTTAAA
CATGCACCTATGGGACCTGCTGGGTGATATCTCAGAAAATGTGAGAGTAGAAGTTCCCAATACAGAC
TGCAGCCTACCTACCAAGTCTTCTGGATTGCTGGAATTGAAAATTAGCAGGTTACAATGCCCTTTTA
AGATATGAAGGATTTGAAAATGACTCTGGTCTGGACTTCTGGTGAATATATGTGTTCTGATATCCAT
CCAGTTGGTTGGTGTGCAGCCAGCGGAAAACCTCTTGTTCTCCTAGAACTATTACGATAAATATACA
AACTGGAAGCTTTTCTAGTGAAACGACTTACTGGTGCCAAAACACTGCCTCCTGATTTCTCCAAAAG
GTTTCAGAGAGTATGCAGTATCCTTTCAAACCTTGCATGAGAGTAGAAGTGGTTGACAAGAGGCATTTG
TGTCGAACACGAGTAGCAGTGGTGGAAAGTGTAATTGGAGGAAGATTAAGACTAGTGATGAAGAAAGC
GAAGATAGAACAGATGACTTCTGGTGCCATATGCACAGCCCATTAATACATCATATTGGTTGGTCTCGA
AGCATAGGTATCGATTCAAAAGATCTGATATTACAAAGAAACAGGATGGACATTTTGATACACCACCA
CATTTATTTGCTAAGGTAAAAGAAGTAGACCAGAGTGGGAATGGTTCAAGGAAGGAATGAAATTGGAA
GCTATAGACCCATTAAATCTTTCTACAATATGTGTCGCAACCATTAGAAAGGTGCTAGCTGACGGATTC
CTGATGATTGGGATCGATGGCTCAGAAGCAGCAGACGGATCTGACTGGTTCTGTTACCATGCAACCTCT
CCTTCTATTTTCCCTGTCGGTTTCTGTGAAATTAACATGATTGAAGTTACTCCACCCAGAGGTTACACA
AACTTCTTTTAAATGGTTTGACTACCTCAGGAAAAGTGGTCCATTGCAGCACCAGTAAAATATTT
AATAAGGATGTTCCAAATCACGGATTTCTGTAGGAATGAAATAGAAGCAGTAGATCTCATGGAGCCA
CGTTTAATATGTGTAGCCACAGTAAGTCAATATTCATCGTCTCTTGAGGATACATTTTGATGGATGG
GAAGAAGAGTATGATCAGTGGGTAGACTGTGAGTCACTGACCTCTATCTGTAGGGTGGTGTGAGTTA
ACTGGATATCACTACAGCCTCCAGCATCACAGTCATCAAGAGAAAACCAATCAGCTTCATCAAAACAG
AAGAAAAAGGCTAAGTCCAGCAATACAAAGGACATAAGAAAATGACAACACTGCAGCTGAAGGAGGAG
TTGCTGGATGGAGAGGATTATAATTTCTTCAAGGAGCGTCTGATCAGGAAAGCAATGGCTGCCAAC
TTCTACATCAACAAGAGCCATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_017643

Insert Size: 1887 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_017643.2</u>
RefSeq Size:	5397 bp
RefSeq ORF:	1887 bp
Locus ID:	54799
UniProt ID:	<u>Q05BQ5</u>
Cytogenetics:	17q21.33
MW:	70.5 kDa
Gene Summary:	Putative Polycomb group (PcG) protein. PcG proteins maintain the transcriptionally repressive state of genes, probably via a modification of chromatin, rendering it heritably changed in its expressibility (By similarity). Specifically binds to monomethylated and dimethylated 'Lys-20' on histone H4.[UniProtKB/Swiss-Prot Function]