

Product datasheet for **SC318083**

GSTCD (NM_001031720) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GSTCD (NM_001031720) Human Untagged Clone
Tag:	Tag Free
Symbol:	GSTCD
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC318083 representing NM_001031720.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAAAGCCATAAAGAAAAGTCTTACAGAAGAAGAATACCTGTACCTGGACTTTTCTACCAAAACAGAA
GGATGCATCTTTCTCTCATAACATCTGTAACCTTTATTTCTGTTATCTTACTGTGACTGTAAAATCTTT
AAAATTTGCTTAGTTGTACCAAAGAGGTGAGTAGAGATAGTTCACTACTAAGAGATGACCTGATCCAG
GATGTTGAAATACAGATTATTTCAAGGCAGGAGCTCCCACCAATAGTCCAAAATTGCTGTTTGCCTGCA
GTAGTAGAACGATCAGACAATTTTGTAGAGCAGGACTTGCTGTTGTATTGAGACACATAATCCAGAAA
TCCTATGAAGCAGACCCCTTAAAGAAGGAACTTTTGAACCTCTGGGCTTTAAAAAGACTTGCTTGAAA
GCCTGTGCTGAAGTTAGTCAGTGGACCAGGCTATGTGAACCTCACCATCCCTTTGGCTATTGAGAATTTT
CTCAGAGAATCTTCTGACCAGCCCCAACTATACCTGTAGAAATACTACAGCTAGAGAAAAGCTTAGT
GAGCCTGTTAGAGTGCATAATGATGATAAACTCCGCAGGCAGAAGCTCAAGCAACAGAAGGCTGATGGA
GTTGGGCTCCCCTTACTAAGGGAAAGGCAAGAGCAAGGTCCACACACAGGAAACATCTGAAGGGTTG
GATTCTTCATCCAAGAGTCTGGAACGAAAGTGGCATTCTCAAAGCTCACAGTACAGGAAGAACCAGCT
ACTACCAACAGAGAGCCTTCTCACATCAGAAAAGCAAAAGCCTCCGACCTTCCACCTCTGGAGCATGTG
TTTGACAGAAAGGGCTTTACTTCACTCTGGCAGATATTGTGCTCTTGCCCTGTATCCATCATTTCTGGTA
ATTATCAGCAGGAAATTTCTGAGAAGCTGGTAGAATTTCCATTGCTAGCCTCTTGGTACCAGAGGATT
CAGGAAGTGCCAGGAGTAAAAACAGCAGCTTCTAAGTGTGGGATCCAATTTCTCCATTTACCAAAGTTG
TTGACAACCTCAACTGAACAGCATCCAACTTATGTGAAGTCCCAGGTGTAGAAGAGCAAAGCGATCCT
TTATTTATAGGAGGACCAAGACCAACCATGGCCAAGTTAATGAAAAGGGCATTGAAGTGTGTTTTCT
CCCCACCCTTGCCTACTTGGACTCTTGATTGGAATGTTCTCCCTGCAGCAGTCAGCCCAAAGGAAGT
AAAATGTCCAGTGATCGAGCTTTGAGGAAGCAGCAACAGTTGAACAACCTTGCTATGTGTTAACAAAT
CAGGCCAAACCTGGTGACAGAATTGTGGATTTCTGCAGCGGTGGGGCCATGTTGGAATTGCTCTTGCT
CACATGCTGCCATCATGTCAGGTTACATTAATAGAAAACAAGGAATTATCATTAAATTCGTGCTAAGAAG
AGAAGTGATGAACGGGTTAAGCAACATTTGGTTCATTCAAGCAAAATATGGAATATTTTACTGGGATG
TTAATATTGGAGTAGCATTGCATGCTTGTGGAGTGGCAACAGACATGGTATTGAGCACTGTATCAAA
ACACGGGCTTCTTCGTACATGCCCTTGCTGTTATGGTTTCATTGAGAACACCTCAAAGTTCAATTTT
CCAAAAAGTGAACAATTCAGAAAACCTTTATCATAACAAGGAACACATGATTCTGTGCAGATTTGCAGAC
CAGACAGCTGTCCAGCTCCCACCCCAACGAAGGCTCATAGGAAAACAGTGCATGTGCTTGGTGGATCTG
GATCGAGCAAGAGCTGCAGAAGAATGTGGATACTCCGTTCAAGTGATATCCATGGAGCCAGAGAGCTGC
TCTCCAAAAAATAACATGATTGTGGGAGTCCCCATTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_001031720

Insert Size: 1902 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001031720.3](#)

RefSeq Size: 4122 bp

RefSeq ORF: 1902 bp

Locus ID: 79807

UniProt ID: [Q8NEC7](#)

Cytogenetics: 4q24
Protein Families: Druggable Genome
MW: 71.1 kDa