

Product datasheet for **SC318093**

NUMB (NM_003744) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NUMB (NM_003744) Human Untagged Clone
Tag:	Tag Free
Symbol:	NUMB
Synonyms:	C14orf41; c14_5527; S171
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC318093 representing NM_003744.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGAACAAATTACGGCAAAGTTTTAGGAGAAAGAAGGATGTTTATGTTCCAGAGGCCAGTCGTCCACAT
CAGTGGCAGACAGATGAAGAAGCGTTCCGACCCGGAAAATGTAGCTTCCCGTTAAGTACCTTGGCCAT
GTAGAAGTTGATGAATCAAGAGGAATGCACATCTGTGAAGATGCTGTAAAAAGATTGAAAGCTACTGGA
AAGAAAGCAGTTAAAGCAGTTCTGTGGTCTCAGCAGATGGACTCAGAGTTGTGGATGAAAAAACTAAG
GACCTCATAGTTGACCAGACGATAGAGAAAGTTTCTTTCTGTGCCCCAGACAGGAACCTTTGATAGAGCC
TTTTCTTACATATGCCGTGATGGCACCCTCGTCGCTGGATCTGTCACTGCTTCAAGGCTGTCAAGGAC
ACAGGTGAAAGGTTGAGCCATGCAGTAGGCTGTGCTTTTGCAGCCTGTTTAGAGCGCAAGCAGAAGCGG
GAGAAGGAATGTGGAGTGACTGCTACTTTTGTAGTGTAGTGGACCACTTTTACAAGAGAAGGATCATTC
CGTGTACAACAGCCACTGAACAAGCAGAAAGAGAGGAGATCATGAAACAAATGCAAGATGCCAAGAAA
GCTGAAACAGATAAGATAGTCGTTGGTTCATCAGTTGCCCTGGCAACTGCCCATCCCACCTCTCT
CCACCTCTCCTACTTCTGATGCCACGACTCTCTGGAGATGAACAATCCCTCATGCCATCCCACGCCGG
CATGCTCCAATTGAACAGCTTGCTCGCCAAGGCTCTTCCGAGGTTTTCTGCTCTTAGCCAGAAGATG
TCACCTTTAAACGCCAACTATCCCTACGCATCAATGAGTTGCCTTCCACTATGCAGAGGAAGACTGAT
TTCCCCATTAATAATGCAGTGCCAGAAGTAGAAGGGGAGGCAGAGAGCATCAGCTCCCTGTGCTCACAG
ATCACCAATGCCTTCAGCACACCTGAGGACCCCTTCTCATCTGCTCCGATGACCAAACAGTGACAGTG
GTGGCACCACAATCTCCTACCTTCCAAGCTAATGGCACTGACTCAGCCTTCCATGTGCTTGTAAAGCCA
GCCATACTGCTCTAGCACCCGTAGCAATGCCTGTGCGTGAACCAACCCTTGGGCCATGCCCTGAT
GCTGCTAACAAGGAAATTGCAGCCACATGTTCCGGGACCAGTGGGGTCAATCTTCTGGTGTGCTGCTCT
CCAGGTCTCTTCCAGGCCGGTCATAGACGTACTCCCTCTGAGGCCGACCGATGGTTAGAAGAGGTGTCT
AAGAGCGTCCGGGCTCAGCAGCCCCAGGCTCAGCTGCTCCTCTGCAGCCAGTTCTCCAGCCTCCTCCA
CCCAGTGCATCTCCCAGCCAGCATCACCTTTCCAAGGGAATGCATTCTCACCTCTCAGCCTGTGCCA
GTGGGTGTGGTCCCAGCCCTGCAACCAGCCTTTGTCCCTGCCAGTCTATCCTGTGGCCAATGGAATG
CCCTATCCAGCCCTAATGTGCTGTGGTGGGCATCACTCCCTCCAGATGGTGGCCAACGTATTTGGC
ACTGCAGGCCACCCTCAGGCTGCCATCCCCTCAGTCAACCAGCCTGGTCAAGCAGCAGACATTCCT
CACTACGAGGCAAGCAGTGCTACCACAGTCCCTTCTTAAAGCCTCCTGCTCAGCACCTCAACGGTTCT
GCAGCTTCAATGGTGTAGATGATGGCAGTTGGCCTCAGCAGACAGGCATACAGAGGTTCTACAGGC
ACCTGCCAGTGGATCCTTTTGAAGCCAGTGGGCTGCATTAGAAAATAAGTCCAAGCAGCGTACTAAT
CCCTCCCTACCAACCCTTTCTCCAGTACTTACAGAAGACGTTTGAATTTAACTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_003744

Insert Size: 1923 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_003744.5](#)

RefSeq Size: 3614 bp

RefSeq ORF: 1923 bp

Locus ID: 8650

UniProt ID: [P49757](#)

Cytogenetics: 14q24.2-q24.3

Domains: PID

Protein Pathways: Notch signaling pathway

MW: 69.4 kDa

Gene Summary: The protein encoded by this gene plays a role in the determination of cell fates during development. The encoded protein, whose degradation is induced in a proteasome-dependent manner by MDM2, is a membrane-bound protein that has been shown to associate with EPS15, LNX1, and NOTCH1. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2016]
Transcript Variant: This variant (3) lacks an alternate in-frame exon compared to variant 1. The resulting isoform (3) has the same N- and C-termini but is shorter compared to isoform 1.