

Product datasheet for **SC318096**

XKR6 (NM_173683) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	XKR6 (NM_173683) Human Untagged Clone
Tag:	Tag Free
Symbol:	XKR6
Synonyms:	C8orf5; C8orf7; C8orf21; XRG6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318096 representing NM_173683.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGGCGAAATCCGATGGCGGTGGCGTGGGGTGGGCTTCGCTCAGCTGCACAACCTGGACGAGGCG
GTGGGCGCGCGCGGAGGAGGACGGGAGCCCGGGGAGCGCGCTGCGCGCGCGCGCGGACGGCAGC
GAGCCCGCGAGAGCAGCTCGATGCACATCTGCCACTGCTGCAACACCTCCTCGTGCTACTGGGGCTGC
CGCTCCGCTGCTGCGCTCCCTCCTGGGCAGGAAGCCGCGCGCAGCGCGCGCGCGGACGGGGGGAC
CAGCCGCTGCAGCCTCCCGCGGCCCGCGCGCGGCCCAACCCCGACGCCCTCGGCCGCGCGGCCG
GAGCCCGCGCGCGCAGGTGGAGCGGCCGTGGCTCGACTGCCTGTGGATCGTGCTGGCGCTGCTGGTG
TTCTTCGGGGACGTGGGCACCGACTGTGGCTGGCCCTCGACTACTACCGCAAGGGGACTACGTCTAC
TTCGGGCTGACCCTCTTCTCGTGCTGGTGCCGTGCTGCTGGTGAGAGCCTGAGCTTCCGCTGGTTC
GTGCAGGACTACACGGGCGCGGGCTGGCGCCGTGGAGGGCTCACCAGCCGGGCCCCCATGATG
GGGGCCGCTACGTCCACGGCGCGGCCGCGGTGGCCAGGCGTGAGGGTCTCCCCACGCCGGGGGCG
CAGCGCCTGTGTGCGCTCTCCGTGTGGATCTGGCAGTCGTCATCCACCTGCTGCAGATGGGGCAGGTG
TGGAGGTATATCCGCACCATGTACCTGGGATTAGAGCCAGCGCGGAAGGAACACCAGCGACGTTTC
TACTGGGCTATGATGATGAATATGCAGACGTCAACATGCTGCGCCTCCTGGAGACCTTCTGGAGAGC
GCGCCCCAACTGGTGCTACAGCTCTACATCATGCTCCAGAAGAACAGCGCGGAGACCTGCCCTGTGTC
TCCTCTGTGACTTCCCTGATGTCCTGGCTTGGTGCTAGCCTCCTATCACAAGCTGTCGCGGACTCC
AGGGACGACAAGAAGAGCATGAGCTACAGAGGGCCATCATCCAGGTCTTCTGGCGCCTCTTACCATC
TCATCCCGAGTGATCTCTTTGCCCTCTTGTCTCCATCTTCCAGCTCTATTTGGGATCTTCGTGGTG
GTTCACTGGTGCGCCATGGCCTTCTGGATCATCCATGGCGGAACAGACTTCTGCATGTCCAAGTGGGAG
GAGATCCTCTTCAACATGGTGGTAGGGATCGTGACATTTTCTGCTGTTTAAACGTCAGGAAGGGCGG
ACTCGATATCGAATGTTTGCATATTATACGATAGTCTTGACCGAGAATGTGCCTTGACGTTCTTTGG
TATTTTTACAGAGACCGGAGACCACTGACTCCTATGCGGTGCCAGCACTGTGTTGTGCTTTATTAGC
TTTGTGGCTGGGATCGCAATGATGCTCTTATACTATGGCGTGCTGCATCCACAGGACCACGAGCTAAG
ATCCTTGCCAGCTCCTGTTGTGCCGAGCTGCTCTGGGGCATCCCTTTGCCCCCGATGTTGAGCCCATG
GCGCCTGAGATCCCTGGGTACCGGGGGACCCAGGTTACGCCACCAGAGCCGTACGGAACAACAGGAG
GATCTCACGGCTGACACTTGCTTGCTGCTGTTTCCAAGTGAGACCATGGGGCCCCCTACCCGTTGGGG
CGTCCTTACCTCCAGAAGGGCCCTCATTAAAGATTGACATGCCAAGAAAGCGATACCCAGCTTGGGAT
GCTCATTTTGTAGACAGGAGGCTGAGAAGGACTATTAACATTCTACAATATGTCACCCCCACCGCAGTA
GGCATTGATATCGAGACGGACCACTCCTCTATGAGTTGCTACAGTATGAGTCTTCACTCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_173683

Insert Size: 1926 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_173683.3](#)

RefSeq Size: 3761 bp

RefSeq ORF: 1926 bp

Locus ID: 286046

UniProt ID: [Q5GH73](#)

Cytogenetics: 8p23.1

Protein Families: Transmembrane

MW: 71.6 kDa