

Product datasheet for **SC318100**

ZNF74 (NM_003426) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF74 (NM_003426) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF74
Synonyms:	COS52; hZNF7; ZFP520; ZNF520
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318100 representing NM_003426.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGAGATCCCTGCCCGGAGCCCGAGAAGACAGCTCTTCTCTCAGGATCCTGCTCTTTCCCTGAAA
GAGAATCTCGAGGATATATCGGGTTGGGTCTTCCCGAAGCCAGGTCCAAGGAATCGGTGAGTTTCAAG
GATGTGGCTGTGGACTTCACCCAGGAGGAGTGGGGTCAACTAGACTCCCCTCAGAGGGCCTGTACCGG
GATGTGATGTTGGAGAATACCAGAACCTTCTTGCCCTAGGACCTCCACTGCACAAGCCAGATGTGATC
TCTCATCTGGAACGAGGCGAGGAGCCATGGAGCATGCAGAGGGAAGTCCCCAGAGGGCCCTGTCCAGAA
TGGGAGCTGAAGGCGGTGCCCTCTCAACAGCAGGGCATTGCAAAGAAGAACGGGCCAGGAGCCCATC
ATGGAGCGGCCCTCGCGGGGCGCAGGCGTGGGGGCCAGGCGAGTGTCTGCAGAGGAGTCAAGCT
GCGCCCTGGGCGCCGACCTGCCATGGTCTGGGACGTCCCTGTAGAGGAATCCCCCTCAGGTGTCCC
CTCTTCGCCCAGCAACGCGTTCCCGAGGGGGGACCTTGTGTGACACACGCAAGAAGCTCCAGGCCACT
GAGGGCAGAACCAAGGCCCGCGAGACTGTGTGACGGGAAAACGCCTCCACGCCAAGTGAAGCCAGAA
AAGTTCCCCAGGTGCGCCGGCAGCGCGGGGCGGCGCCGGGAGGGCGAGTTTCTGTGCGGCGAGTGC
GGGAAGGCGTTCCGCGAGAGCTCCTCCCTCACGCTGCACCGGCGCTGGCAGCCGGGAGAAGGCTTAC
AAGTGCATGAATGCGGCAAGGCCTTACCTGGAGCACCACCTTCTGGAGCACCAGGCGCATCCACACC
GGCGAGAAGCCCTTCTTCTGCGGCGAGTGGGGAAGGCCTTCACTGCCACTCGTCCCTCAACGTGCAC
CAGCGCATCCACAGGGCGAGCGGCCCTACAAGTGCAGCGCTGCGAGAAGGCCTTCACTGCAGCTCG
CTGCTCAGCATGCACCTGCGGGTGCACACCGGCGAGAAGCCCTACCGGTGCGGCGAGTGGGCAAGGCC
TTCAACACAGCGTACACACCTCACACGCCACCACCGCATCCACAGGGCGAGAAGCCCTACAGTGGGC
TCCTTCGCGCAAGGCCTTCACTGCCACTCATCCCTCACCCTGCATGAGAAGATCCACAGCGGGGACAAG
CCGTTTCAAGTGCAGGACTGCGAGAAGGCCTTCAACAGCCGCTCGCGCCTCACCCTCCACAGAGGACG
CACACGGGCGAGAAGCCCTTCAAGTGCAGGCGACTGCGGGAAGGGCTTCACTGCCACGCGTACCTGCTC
GTGCACCGGCGCATCCACAGCGGCGAGAAGCCCTTCAAGTGAACGAGTGGGCAAGGCCTTCACTGCC
CACGCTACCTCATCGTGCACCGGCGCATCCACAGGCGAGAAGCCCTTCACTGCAGCCAGTGTGG
AAGGCCTTCACTGCCACTCGTCCCTCATCGTGCACCGGCGCATCCACAGCGGTGAGAAGCCCTACAAG
TGCAGCGAGTGGGCGAGGCCTTCAAGCAGAACCACTGTCTCATTAAACATCAGAAAATCACTCCGGG
GAGAAGTCGTTTAAAGTGTGAGAAATGTGGGAGATGTTCAACTGGAGCTCGACCTCACTGAGCACCAG
AGGCTGCACAGCGAGGGGAAGCCCTTGGCCATCCAGTTCAACAAACACCTGCTCAGCACATACTACGTG
CCTGGCAGCCTGCTGGGTGCAGGGGATGCTGGACTGAGGGACGTGGATCCCATCGACGCGCTGGATGTG
GCAAAGCTCTTGTGCGTGGTTCCCCCAGAGCTGGCAGGAATTTCTCCCTGGGGAGCAAACCTCGAAAC
TAA
AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

Restriction Sites: SgfI-RsrII

ACCN: NM_003426

Insert Size: 1935 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_003426.3</u>
RefSeq Size:	3954 bp
RefSeq ORF:	1935 bp
Locus ID:	7625
UniProt ID:	<u>Q16587</u>
Cytogenetics:	22q11.21
Domains:	zf-C2H2
Protein Families:	Transcription Factors
MW:	72.2 kDa
Gene Summary:	May play a role in RNA metabolism.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longest isoform (a). Both variants 1 and 3 encode isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.