

Product datasheet for **SC318103**

SYT16 (NM_031914) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SYT16 (NM_031914) Human Untagged Clone
Tag:	Tag Free
Symbol:	SYT16
Synonyms:	CHR14SYT; Strep14; SYT14L; syt14r; yt14r
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318103 representing NM_031914.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGTGTGGCCATGGCGTCTCAGGATGTTTCAAGAACTTCTCCAGCCTTTCTCTTCTGGATATCTCGG
GTTTATGAAGCTCTCCAGCAAGCAGGAGATATGTTATCTGCTTCGCTGGTTAACATAAGCAAACAAGAC
TCTAAATTGAGTGACAACTAGATCAGGACTTAGATAAATTTAGATTTCAGGAAACGTACTTTGAAGAT
GAAGAACAAGACAATGATTGGAGTCAAGAGGATGCAAAATTCCTTGTCTTGAAGTGATCATTTCTCA
TGTTGTAATAGTGATTTGCAGGACTCTGCCAAAATTCAGCCCAAGCCTTAGCCAACATGCAAAGGAC
TCATGTTCCACAATGTCCAGTGGCCCAATTGGGCCAGTGATGACCGCAAGTTACCACATGTGCTTTCT
TCTATTGCGGAGGAAGAGCATCACCTTGAAAAGCAAAGAAGTGGCCTTCAACATGGCTTTGACAGCCAG
CTCCCTGGTACTTTAGAACTGTTAATGGAAAAAGCAAGTCAACAGCTTTGGGGATGACGAAGAGCTG
TCCACATCTTCTGACAGTGACGAGGAGTGATCAAAACAATTTGAGATTTCCGTGTCCCGTCCCAGAGT
TTCCGTTCAGTGACATCTGAGAAAGGAAAGCAGACAGGATTGGAGCAGAAACCAAAATTCAGCCGTTCTG
TTGTTGACACACGGAGAAGATGGCACAGAAGTATCTGCCTGCGAAGATTTGGATGGAGCCAGCCAACGG
CGTTATTCTGAGAATCTCTCCTACGGTGAAGATGACCACATCCCTGCTCACTCACAGTCCCCATGTGAA
AGAGGGGATGCCAAACACCACGGCACATCTACCAAGAGTCCAGTGTTGGTCCAAAGCCTCAGGCGCCAA
TCCACAGAGGGCAGCTTGGAGATGGAGACAGCTTTTAATAGCCGGGGATTTGAAGATTCCTATGCCACT
GACAGCTCCTCCATGTGGAGTCCAGAGGAACAGGACAGGACCAATTTGCAGGTGCCATCCGGGGTCTCA
GAGCCCATCTCAAAGTGTGGTACCTAGATGTCATCTTTGAATATAGAGCCGCCAGCCAGAAGCTCACA
GTGACCATTTGTGAGGGCACAGGGCCTCCAGATAAGGACCGAAGTGGTGTCAACTCCTGGCAAGTTTCA
GTAGTGCTGCTGCCTGGTAAGAAACACAGGGGCAGGACGAACATACAGAGAGGGCCCAACCCCGTCTTC
AGGGAGAAGGTACCTTTGCCAAGCTGGAGCCAGAGATGTGGCTGCCTGTGCTGTCCGCTTCCGCTG
TACGCTGCCCGGAAGATGACCCGAGAGAGAATGATGGGAGAGAACTATTCTATCTCAGCCACCTGCAC
CCAGAAGGGGAAATGAAAGTGACTCTGGTTCTGGAGCCAAGAAGTAATATAAGCAGTGGAGGGTCTCCG
CTCAGCCCATCTGCGTTTCTCAGAGTATAGTACTTATCCACGCAGTCGCTGTCTCATGGAGGGGCG
CCAGAGCTGTTGGTGGGGCTCTCGTACAATGCCACAACGGGGCGATTATCTGTGGAATGATCAAAGGC
AGCCATTTCCGAAACCTCGCTGTTAACCGAGCACCTGATACATATGGAAGTCTTTCTCCTCAATTCT
GTGGGTCAAGAGATGTCCGTTGCAAGACGTCCATTTCGGCGTGGTCAGCCCAATCCTGTCTATAAGGAG
ACCTTTGTTTTCCAGGTGGCCCTCTTTAGCTGTCTGATGTACGTTGATGATTTCCGTTTATAACAGG
CGTACTATGAAGCGTAAAGAGATGATTGGCTGGATTGCCCTGGGCCAGAACAGCAGTGGAGAGGAGGAA
CAAGATCACTGGGAGGAGATGAAGGAAACCAAGGCCAGCAGATCTGCAGATGGCACACTTTGCTGGAA
TCCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_031914

Insert Size: 1938 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_031914.2</u>
RefSeq Size:	3137 bp
RefSeq ORF:	1938 bp
Locus ID:	83851
UniProt ID:	<u>Q17RD7</u>
Cytogenetics:	14q23.2
MW:	72.3 kDa
Gene Summary:	May be involved in the trafficking and exocytosis of secretory vesicles in non-neuronal tissues. Is Ca(2+)-independent.[UniProtKB/Swiss-Prot Function]