

Product datasheet for **SC318117**

ZNF12 (NM_006956) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF12 (NM_006956) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF12
Synonyms:	GIOT-3; HZF11; KOX3; ZNF325
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318117 representing NM_006956.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGAATAAATCCCTGGGGCCAGTGTCTTCAAGGACGTGGCTGTGGACTTCACCCAGGAGGAATGGCAG
CAGCTGGATCCTGAGCAGAAGATACTTACAGGGATGTGATGCTGGAGAACTACAGCAATCTAGTTTCT
GTGGGGTATCACATTATCAAACCGGATGTTATCAGCAAGTTGGAGCAAGGAGAAGAGCCATGGATAGTA
GAAGGAGAATTCCTACTTCAGAGCTATCCAGATGAAGTCTGGCAAATGATGACCTAATAGAGAGAATC
CAGGAAGAGGAAAATAAACCTTCAAGGCAAATGTGTTTATTGAGACCCTGATTGAAGAGAGAGAATAT
ATTAGTAGTGATGGAAGCTATGCAAGAATGAAAGCTGATGAATGTAGTGGATGTGGGAAATCACTCTC
CATATTAAGCTTGAGAAAATCATCCAGGAGATCAAGCTTATGAATTTAATCAAAATGGGAACTTAT
ACTCTAAATGAAGAAAGTCTTTATCAGAAAATTCGTATTTTGGAGAAACCTTTTGAATATATTGAATGC
CAGAAAGCCTTCCAAAAGGACACTGTTTTTGTAAATCACATGGAAGAAAAGCCCTATAAGTGAATGGA
TCTGAAATAGCCTTTCTCCAGATGTCGGACCTCACTGTACATCAGACATCTCATATGGAATGAAGCCC
TATGAATGCAGTGAATGTGGGAAATCCTTCTGTAAAAAGTCAAAATTTATTATACATCAGAGGACTCAC
ACAGGAGAGAAAACCTTACGAATGTAATCAGTGTGGGAAATCCTTCTGCCAGAAGGAAACCTTACTGTG
CATCAGAGAACACACACAGGGGAGAAGCCCTATGAATGTAATGAATGTGGGAAGAACTTTTACCAGAAG
TTACACCTCATTACAGATCAGAGAACTCACTCAGGAGAGAAGCCCTATGAATGTAGTTATTGTGAAAA
TCCTTTTGCCAGAAGACACACCTCACACAACATCAGAGAACACATTACAGGAGAGAGACCTTATGTTTGT
CATGACTGTGGGAAAACCTTCTCGCAGAAGTCAGCACTTAATGACCATCAGAAAATTCACACAGGTGTG
AACTCTACAAGTGTAGTGAATGTGGGAAATGCTTCTGCCGAAGTCTACTCTCACGACCCACCTGAGG
ACCCACACAGGAGAGAAAACCGTATGAATGTAATGAGTGTGGAAAATCTTCTCTCGTTGTCATATCTC
ACTGTACATTATAGAATCATTACAGGAGAGAAAACCTTATGAATGTAATGAATGTGGAAAACCTTCTAC
CTGAATTCAGCCCTCATGAGACATCAGAGAGTGCACACAGGAGAGAAAACCTTACGAATGTAATGAATGT
GGAAAGTTATTCTCCAGTTGTACATCTCACTATCCATCATAGAACTCATTACAGGAGTAAAACCTTAT
GAATGTAGTGAATGTGGGAAAACCTTCTACCAGAATCAGCCCTTTGTAGACATCGGAGAATACACAAA
GGAGAGAAGCCCTATGAATGCTATATATGTGGAAAATCTTCTCTCAGATGTCATACCTCACTATACAT
CATAGAATTCATTACAGGAGAGAAGCCCTATGAATGTAGTGAATGTGGGAAAACCTTCTGCCAGAATTC
GCCCTTAATCGACATCAGAGAACACACACAGGAGAGAAAAGCCTACGAATGTTATGAATGTGGGAAGTGC
TTCTCTCAGATGTCCTATCTCACTATACATCATGAATTCATTACAGGAGAGAAAACCTTTGAATGTAAT
GAGTGTGGAAAAGCCTTCTCTCGGATGTACATCTCACTGTACATTATAGAATCATTACAGGAGAGAAA
CCCTATGAGTGTACTGAATGTGGAAAAAATCTACCACAAATCAGCATTCAACGCCATCAGAGAATT
CATAGGAGAGGCAATATGAATGTAATAGATGTGGGAAGGCTTCTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-MluI

ACCN: NM_006956

Insert Size: 1980 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_006956.2</u>
RefSeq Size:	4960 bp
RefSeq ORF:	1980 bp
Locus ID:	7559
UniProt ID:	<u>P17014</u>
Cytogenetics:	7p22.1
Protein Families:	Transcription Factors
MW:	77.1 kDa
Gene Summary:	This gene is a member of the krueppel C2H2-type zinc-finger protein family and encodes a protein with eight C2H2-type zinc fingers and a KRAB domain. This nuclear protein is involved in developmental control of gene expression. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008] Transcript Variant: This variant (2) lacks an alternate in-frame segment, compared to variant 1, resulting in a shorter protein (isoform b), compared to isoform a.