

Product datasheet for **SC318124**

KY (NM_178554) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KY (NM_178554) Human Untagged Clone
Tag:	Tag Free
Symbol:	KY
Synonyms:	MFM7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318124 representing NM_178554.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGAGCTGAAGAAGGACATCAACGCTGTATCTATCGACATGCTGCTGATCGTCACTCGGAGAAGCGG
CGCGCCGCACAGGTACGCTCTCAGACCAGCAGCGCAACCCGAGCTCGCTGCTGCAGCGCGGAGGAGGA
TTCCAAGGTGTTGGAAATGGAGTCCGAAGATGGCAGAAATTAGAAGGAAATGACTTTCACGAAAACCTTG
GTGGAGAAGCAGCACCCCTCAGCAGCCCCAGGTCATCACTTCTACAACAGCCAAGGGACACAACCTGACT
GTGGAAGTCCACCCTCGAGATGCCATGCCCCAGCTGCTCAAGAAGTTCTCCTTGGCTAAACGTTTACAA
GGTGATAAAAAATGGAAACACAAGACCCCGCAACCTGGAGGGAAGATGCCCATGCCTACCCCTGGGAT
CGATCTAGCCTGAAATCCATGTCCTGGACCTGCAGCAGTTTGAGAAATTGGATATCTACGCTCACAG
GTGACAGCCAAGAGTGGCTAGACGAACGGTGAGTGACCTGCTCCAGGAGGCCACACTGACCTGGAA
AGGGTCCGCGCCATCTGGATCTGGATCTGCCATCACATAGAGTATGACATTGCAGCTGCTCAGGAGAAG
GACCGCCAAGCCTTCAAACCCACTGACATCCTGCGGACCCAGAAGACCAACTGTGATGGCTATGCTGGC
CTCTTCGAGAGAATGTGAGGCTCGCCGAGTGCAAGTGTATGACCGTGCCTGGCTACTCCAAGGGTTTC
GGCTACCAGACAGGGCAGAGCTTCTCGGGGAGTTTGACCATGCCTGGAATGCTGTGTACCTGGAGGGA
AGATGGCACCTGGTGGACAGCACCTGGGGCAGCGGCCTGGTGGACACCATCACCTCCAAATTCACCTTC
CTCTACATGAGTTCTACTTCTCACCCACCCTGCACTGTTTATCGAGGACCACTTCCAGACAACAAG
AACTGGCAGCTGCTAAAACCTCTCAATCTCTGAGGCAGTTTGAGAACACATGTATCACAAGAGTGAA
TTCTACAACAAAGGGATGCTGAGTGCCACCCAGAGACTTCCATGATCAGAACAGTGAATGGGAAGGCC
ACGGTCACCATTGAGAGCTGCGCCCCGACGCTGTTTCTGTTTCAATGGCAAGCAAGAGCATGGG
CTGCTGAGCCTAAGGAAGAATGGGATGAAGTTGGAGGTGTACCCTCCAACCATGGGCACTCACAAGCTG
CAGATCTTTGCCAAGGGCAACTCCGACATCTACAGCTCAGTGCTGGAGTACACGCTCAAGTGCAATTAT
GTGGACATGGGTGTCCAGCTGCCTGCTGAGCTTACCAGCCCGTGGGCCCCAGCTGGTTCTCGGAGCAG
ATGGGCATCATGAAGCCCTCCACCCTGACCCTATCATCCACACCAGCGACGGGCGCTGCTCCATCAGC
TTCAGCGTGGAGGAGGCATTAATGTCTGGTTCCCTCCAGGGGATGATGGCCCCATCACTGAGGAG
ACACAGCGGCGCTACATCTTCCAGCTGCACCGGGAAGCAGACCGAGCTGAAAGTCCAGCTGCCCCAT
GCAGGCAAGTTTGCCCTCAAGATCTTTGTCAAGAAGAGGCAGGAACCAGGAACTACATCTTCGCTTT
AATTACCTTGATGCTGTGCCAACCAAGGTGAAGTGGCCCATGTTCCCTGAGAGCTTTGGCAACTGG
GGACAGGACAATGAGCTGCTGGAACCTCTGTCAGGTGTGCTTCTGCCAACCGGAATGTCCCATTTAAA
TTGAAGCTGCATGGTATTGCCAAGTCCCTGGTGAAGGGGCAGGACACATGGCCCCGACCTGAACCAC
GAGGGCTACTGGGAGGGAAGCTGCAGCAGAGCTGGCTGCCAGGAAGTCTATGTCATGGTGTGGAGAAT
GCCAACCACAATTTCTACTCCTACATCCTGAAATACAAAGTGAATGCCAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_178554

Insert Size: 1986 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_178554.4</u>
RefSeq Size:	5718 bp
RefSeq ORF:	1986 bp
Locus ID:	339855
UniProt ID:	<u>Q8NBH2</u>
Cytogenetics:	3q22.2
MW:	75.2 kDa
Gene Summary:	The protein encoded by this gene belongs to the transglutaminase-like superfamily. The protein is involved in the function, maturation and stabilization of the neuromuscular junction and may be required for normal muscle growth. Mutations in this gene are associated with myopathy, myofibrillar, 7. [provided by RefSeq, Apr 2017] Transcript Variant: This variant (1) encodes the longest isoform (1). Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.