

Product datasheet for **SC318147**

REXO1L1 (REXO1L1P) (NM_172239) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	REXO1L1 (REXO1L1P) (NM_172239) Human Untagged Clone
Tag:	Tag Free
Symbol:	REXO1L1
Synonyms:	GOR; REXO1L1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC318147 representing NM_172239.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTTGCGAGCCACAGCTCCCTGCTGGTTCCCCCTGGATACCCAGAAGCTAAGAAGGTGGCCGAGGAG
GCGGCCCTGGAGCAAGCCGCCATTTGGGAGGGGAGCAGAGCCAGGCCGGTGCTCCCGAAGGCAGCAAG
ATGTTGCGAGCCACAGCTCCCTGCTGGTTCCGCCCTGGATACCCAGAAGCTAAGAAGGTGGCCAAGGAG
GCGGCCCCCGAGGCAAGCCGCCATTTGGGAGCGGAGCAGAGCCCGGCCGGTGCTCCCGAAGGCAGCAAG
ATGTTGCGAGCCACAGCTCCCTGCTGGTTCCACCTGGATACCCAGAAGCTAAGAAGGTGGCCGAGGAG
GCGGCCCTCGAGGCTCCAGAATCCCACTGCCCTCTCATCAGCCTGCCAGAGCTTCGGGCTCTGGGTG
CCCCAGATGCACAAGCAGGCTCAGCATTTGTGGACATCCAGGCGGAGCCCCAAAACAGGGGTCCGGCG
GTGCCCCAGCGTGGCCCAAGATGGTGACGGAGTCGTGCTACTTCCCTGCGCAGAGGGGATCGGCCCTGC
CGCTTGCCAGCCGCCCAAGGCTGACAGAGAGGCCCTCGGGAGTCCGCATCTCAGCCCCAGGAAGAGG
AAGACGATCGCCACTTCCAGCCCTTGCTTGGTCCACAGTTACACAGATGCCAAGAGAACCCGGGTG
GCCAGCAGCAGCAACGCTCCCGTGGCTCCAAGGTCGGCAGACAGCCAGGGAAGACGCGCAACAGGTCA
GGGATGGCATGCAAGACCACCGCCACCACCAGCTCTAAGCGAATCGTCCGTGCTGCATCCTTACCGAGT
TTGAGTTTGAAGAAACCCATTATCCTCCGAAGCTCTGGGTGCCAAGTCCCACCCTCCTCCGCCGAGGC
TATCTCCAAGTGTACCGAAGAGTGTCTCAAGTTCTGCGCCTCCAAGCAGGAGGCCGAGGAGAAGGCG
CTGAACGAGGAGAAGGTGGCCTACGACTGCAGCCCCAACAAAGAACAGGTACCTGAACGTGGTCTGAAAC
ACCCTCAAGAGACTGAAGGGCTGACCCCAAGCTCCATGCCCGCCCTCAGCAGGGCCGCCCTGTACAGC
CGCCTCCAGGAGTTCCTGCTCACCCAGGACCAGCTCAAGGAGAACGGCTACCCCTTCCCGCACCCCGAG
CGGCCCGAGGCGCCGCTCTTCACTGGCCAGGGGAAGGGGCCCGGCGACTCCTCTGTAGGGTCTGC
TGCCGTTGTGGCACCGAGTACCTGGTGTCTCCTCGGGCCGCTGTGTACGCGACCAGTTGTGTTATTAT
CACTGGGGGCGGGTCCGCTCGAGCCAGGTGGTGGAGGCCGGTTAGCCAGTACACCTGCTGTGCAGCT
GCTCCTGGCTCTGTGGGCTGCCAGGTGGCAAAGCAGCAGCTGCGGGACGGCCGAAGGAGAGCCTCGAT
GGCTTCGTGGAGACCTTCAAGAAAGAGTTGTCCAGAGACGCTTATCCAGGAATCTACGCCTTGGACTGT
GAGATGTGCTACACCACGCATGGCCTAGAGCTGACCCGCGTCACCGTGGTGGACGCCGACATGCCAGTG
GTGTACGACACCTTCGTCAAGCCCGACAACGAGATCGTGGACTACAACACCAGGTTTTCCGGAGTACC
GAGGCCGACGTCGCAAGACGAGCATCACGTTGCCCAAGTCCAAGCCATCCTGCTGAGCTTTTTCAGC
GCCAAACCATCCTCATCGGGCACAGCCTGGAGAGCGACCTGCTGGCCCTGAAGCTCATCCACAGCACC
GTGGTGGACACGCGCTGCTCTTCCCGCACTACCTGGGTTTTCCCTACAAGCGCTCCCTCAGGAATCTC
GCGGCCGACTACCTGGCACAGATCATCCAGGACAGCCAGGACGGCCACAACCTCCAGCGAGGACGCAAGC
GCCTGCCTGCAGCTGGTGTGTGGAAGGTCGACAGCGCGCCAGATCCAGCCACGCCACCAGTCCCGC
TCTCCCGCCGCTGGCCTGTCTTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_172239

Insert Size: 2028 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_172239.4](#)

RefSeq Size: 7032 bp

RefSeq ORF: 2028 bp

Locus ID: 254958

Cytogenetics: 8q21.2

MW: 73.9 kDa