

## Product datasheet for **SC318151**

### **RBM35A (ESRP1) (NM\_017697) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	RBM35A (ESRP1) (NM_017697) Human Untagged Clone
Tag:	Tag Free
Symbol:	RBM35A
Synonyms:	DFNB109; RBM35A; RMB35A
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC318151 representing NM\_017697.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGACGGCCTCTCCGGATTACTTGGTGGTGCTTTTTGGGATCACTGCTGGGGCCACCGGGCCAAGCTA
GGCTCGGATGAGAAGGAGTTGATCCTGCTGTTCTGGAAAGTCGTGGATCTGGCCAACAAGAAGTGGGA
CAGTTGCACGAAGTGCTAGTTAGACCGGATCAGTTGAACTGACGGAGGACTGCAAAGAAGAACTAAA
ATAGACGTCGAAAGCCTGTCTCGGCGTCGACGCTGGACCAAGCCCTCCGACAGTTTAACAGTCAGTG
AGCAATGAACTGAATATTGGAGTAGGGACTTCCTTCTGTCTGTACTGATGGGCAGCTTCATGTCAGG
CAAATCTGCATCCTGAGGCTTCCAAGAAGAATGTACTATTACCTGAATGCTTCTATTCTTTTTGAT
CTTCGAAAAGAATTCAAGAAATGTTGCCCTGGTTCACCTGATATTGACAACTGGACGTTGCCACAATG
ACAGAGTATTTAAATTTGAGAAGAGTAGTTCAGTCTCTCGATATGGAGCCTCTCAAGTTGAAGATATG
GGGAATAAATTTAGCAATGATTTAGAGCCTTATAATCACAGGTTTTAGATCCAGAGAGAGTGAAT
TACAAGTTTGAAGTGGAACTTGCAGCAAGATGGAAGTTATTGATGATAACACCGTAGTCAGGGCACGA
GGTTTACCATGGCAGTCTCAGATCAAGATATTGCAAGATTCTTCAAAGGACTCAATATTGCCAAGGGA
GGTGACGCACTTTGTCTGAATGCTCAGGGTCAAGGAACGGAGAAGCTCTGGTTAGGTTTGAAGTGAG
GAGCACCGAGACCTAGCACTACAGAGGCACAAACATCACATGGGGACCCGGTATATTGAGTTTACAAA
GCAACAGGTGAAGATTTCTTAAATGCTGGTGGTACTTCCAATGAGGTAGCCAGTTTCTCTCCAAG
GAAATCAAGTCATTGTCGCATGCGGGGGCTCCCTTTCACGGCCACAGCTGAAGAAGTGGTGGCCTTC
TTTGGACAGCATTGCCCTATTACTGGGGGAAAGGAAGGCATCCTTTTGTACCTACCCAGATGGTAGG
CCAACAGGGGACGCTTTTGTCTCTTTGCTGTGAGGAATATGCACAGAATGCGTTGAGGAAGCATAAA
GACTTGTGGGTAAAAGATACATTGAACCTTTCAGGAGCACAGCAGCTGAAGTTCAAGTTCAGAGTGTGAAT
CGATTCTCTCGGCCCTCTCATTCCACTTCCAACCCCTCCATTATTCCAGTACTACCTCAGCAATTT
GTGCCCTTACAATGTTAGAGACTGTATACGCCCTTCGAGGTCTTCCCTATGCAGCCACAATTGAGGAC
ATCCTGGATTTCTGGGGGAGTTCGCCACAGATATTCGTAATCATGGGGTTCACATGGTTTTGAATCAC
CAGGGCCGCCATCAGGAGATGCCTTTATCCAGATGAAGTCTGCGGACAGAGCATTATGGCTGCACAG
AAGTGTCATAAAAAAACATGAAGGACAGATATGTTGAAGTCTTTCAGTGTTCAGCTGAGGAGATGAAC
TTTGTGTTAATGGGGGCACTTTAAATCGAAATGGCTTATCCCCACCGCCATGTAAGTTACCATGCCTG
TCTCTCCCTCTACACATTTCCAGCTCCTGCTGCAGTATTCTACAGAAGCTGCCATTTACCAGCCC
TCTGTGATTTTGAATCCACGAGCACTGCAGCCCTCCACAGCGTACTACCCAGCAGGCACTCAGCTCTTC
ATGAATACACAGCGTACTATCCAGCCCCAGGTTGCGCTAATAGTCTTGGCTACTTCCCTACAGCT
GCTAATCTTAGCGGTGTCCCTCCACAGCCTGGCAGGTGGTCAAGTGCAGGGCCTGGCCTACAATACT
GGAGTTAAGGAAATTTCTAATCTTCCAAAGTTACCAGTATGCAACCGAGGATGGACTTATACACACA
AATGACCAGGCCAGGACTCTACCAAAGAATGGTTTTGATTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**ACCN:** NM\_017697

**Insert Size:** 2046 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_017697.3</u>
<b>RefSeq Size:</b>	3806 bp
<b>RefSeq ORF:</b>	2046 bp
<b>Locus ID:</b>	54845
<b>UniProt ID:</b>	<u>Q6NXG1</u>
<b>Cytogenetics:</b>	8q22.1
<b>Domains:</b>	RRM
<b>MW:</b>	75.6 kDa
<b>Gene Summary:</b>	<p>ESPR1 is an epithelial cell-type-specific splicing regulator (Warzecha et al., 2009 [PubMed 19285943]).[supplied by OMIM, Aug 2009]</p> <p>Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).</p>