

Product datasheet for **SC318155**

GPSM2 (NM_013296) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GPSM2 (NM_013296) Human Untagged Clone
Tag:	Tag Free
Symbol:	GPSM2
Synonyms:	CMCS; DFNB82; LGN; PINS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC318155 representing NM_013296.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGGAAAAATTTGATAAGCATGAGAGAAGACCATTCTTTTCATGTTTCGTTACAGAATGGAAGCTTCT
TGCCTAGAGCTGGCCTTGAAGGGGAACGTCATGTAAATCAGGAGACTGCCGCGCTGGCGTGTCATTC
TTTGAAGCTGCAGTTCAAGTTGGAACCTGAAGACCTAAAAACACTTAGCGCTATTTACAGCCAGTTGGGC
AATGCTTATTTCTATTTGCATGATTATGCCAAAGCATTAGAATATCACCATCATGATTTAACCCCTTGCA
AGGACTATTGGAGACCAGCTGGGGGAAGCGAAAGCTAGTGGTAATCTGGGAAACACCTTAAAAAGTTCTT
GGGAATTTTGACGAAGCCATAGTTTGTGTGACGACACCTAGATATTTCCAGAGAGCTTAATGACAAG
GTGGGAGAAGCAAGAGCACTTTACAATCTTGGGAATGTGTATCATGCCAAAGGGAAAAAGTTTGTGTC
CCTGGTCCCAGGATGTAGGAGAATTTCCAGAAGAAGTGAGAGATGCTCTGCAGGCAGCCGTGGATTTT
TATGAGGAAAACCTATCATTAGTACTGCTTTGGGTGACCGAGCGGCACAAGGACGTGCCTTTGGAAT
CTTGGAACACACATTACCTCCTTGGCACTTCAGGGATGCAGTTATAGCTCATGAGCAGCGTCTCCTT
ATTGCAAAAAGAAATTTGGAGATAAAGCAGCTGAAAGAAGAGCATATAGCAACCTTGGAATGCATATATA
TTTCTTGGTGAATTTGAAACTGCCTCGGAATACTACAAGAAGACACTACTGTTGGCCCGACAGCTTAAA
GACCGAGCTGTAGAAGCACAGTCTTGTACAGTCTTGGAAATACATATACTTTACTTCAAGACTATGAA
AAGGCCATTGATTATCATCTGAAGCACTTAGCAATTGCTCAAGAGCTGAATGATAGAATGGTGAAGGA
AGAGCATGTTGGAGCTTAGGAAATGCATACACAGCACTAGGAAATCATGATCAAGCAATGCATTTTGTCT
GAAAAGCACTTGGAAATTTCAAGAGAGGTTGGGGATAAAAAGTGGTGAACCTAACAGCAGCACTTAATCTC
TCAGACCTTCAAATGGTCTTGGTCTGAGCTACAGCACAATAACTCCATAATGTCTGAAAATACTGAA
ATTGATAGCAGTTTGAATGGTGTACGCCCAAGTTGGGACGCCGCATAGTATGGAAAATATGGAACCT
ATGAAGTTAACACCAGAAAAGGTACAGAACTGGAACAGTGAAATTTCTTGCTAAGCAAAAACCTCTTATT
GCCAAAACCTTCTGCAAAGCTACTCTTGTCAACAGACTGAAGGGGAAAAAATACAAAACGAATTCCTCC
ACTAAAGTTCTCCAAGATGCCAGTAATTCTATTGACCACCGAATTCCAAATTCAGAGGAAAATCAGT
GCAGATACTATTGGAGATGAAGGTTCTTTGACTTATTAAGCCGATTTCAAAGCAATAGGATGGATGAT
CAGAGATGTTGCTTACAAGAAAAGAACTGCCATACAGCTTCAACAACAACCTTCTTCCACTCCCCCTAAA
ATGATGCTAAAAACATCATCTGTTCTGTGGTATCCCCAACACGGATGAGTTTTTAGATCTTCTTGCC
AGCTCACAGAGTCGCCGCTGGATGACCAGAGGGCTAGTTTCAGTAATTTGCCAGGGCTTCGTCTAACA
CAAAACAGCCAGTCGGTACTTAGCCACCTGATGACTAATGACAACAAGAGGCTGATGAAGATTTCTTT
GACATCCTTGAAAATGTCAAGGATCCAGATTAGATGATCAAAGATGTGCTCCACCACCTGCTACCACA
AAGGGTCCGACAGTACCAGATGAAGACTTTTTTCAGCCTTATTTTACGGTCCCAGGGAAAAGAGAATGGAT
GAACAGAGAGTTCTTTTACAAGAGATCAAAACAGAGACACTGACTTTGGGCTAAAGGACTTTTTGCAA
AATAATGCTTTGTTGGAGTTTAAAAATTCAGGGAAAAAATCGGCAGACCATTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_013296

Insert Size: 2055 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_013296.4](#)

RefSeq Size: 3039 bp

RefSeq ORF: 2055 bp

Locus ID: 29899

UniProt ID: [P81274](#)

Cytogenetics: 1p13.3

Domains: TPR, GoLoco

Protein Families: Druggable Genome

MW: 76.7 kDa

Gene Summary: The protein encoded by this gene belongs to a family of proteins that modulate activation of G proteins, which transduce extracellular signals received by cell surface receptors into integrated cellular responses. The N-terminal half of this protein contains 10 copies of leu-gly-asn (LGN) repeat, and the C-terminal half contains 4 GoLoco motifs, which are involved in guanine nucleotide exchange. This protein may play a role in neuroblast division and in the development of normal hearing. Mutations in this gene are associated with autosomal recessive nonsyndromic deafness (DFNB82). Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2016]