

Product datasheet for **SC318157**

MAMDC2 (NM_153267) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MAMDC2 (NM_153267) Human Untagged Clone
Tag:	Tag Free
Symbol:	MAMDC2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC318157 representing NM_153267.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGCTGTTAAGGGGCGTCCTCCTGGCGTTGCAAGCCCTGCAGCTCGCCGGTGCCCTCGACCTGCCCGCT
GGTCCTGTGCCCTTTGAAGAGAGCACTTGGCGCTTTGACTCCGTTGGCTCTCTGCCGTGGATTTTA
AATGAGGAAGGCCATTACATTTATGTGGATACCTCCTTTGGCAAGCAGGGGGAGAAAGCTGTGCTGCTA
AGTCTGACTTACAGGCTGAGGAATGGAGCTGCCTCCGTTTGGTCTACCAGATAACCACATCTTCGGAG
TCTCTGCAGATCCCAGCCAGCTGAACCTCTACATGAGATTTGAAGATGAAAGCTTTGATCGCTTCTT
TGGTCAGCTAAGGAACCTTCAGACAGCTGGCTCATAGCCAGCTTGGATTTGCAAAACAGTTCCAAGAAA
TTCAAGATTTTAAAGAAGGTGTACTAGGACAGGGAAACACAGCCAGCATCGCACTATTTGAAATCAAG
ATGACAACCGGCTACTGTATTGAATGTGACTTTGAAGAAAATCATCTCTGTGGCTTTGTGAACCGCTGG
AATCCCAATGTGAACGGTTTGGTGGAGGAGGAAGTATCGGAATGTCCACTCCATTCTCCACAGGAT
CACACCTCAAGAGTGAACCTGGGCCACTACATGTACGTGGACTCAGTTTATGTGAAGCACTTCCAGGAG
GTGGCACAGCTCATCTCCCCGTTGACCAGGCCCCCATGGCTGGCTGCCTGTCATTTTATTACCAGATC
CAGCAGGGGAATGACAATGTCTTTCCCTTACACTCGGGATGTGGCTGGCCTTACGAGGAAATCTGG
AAAGCAGACAGGCCAGGGAAATGCTGCCTGGAACCTTGCGGAGGTCGAGTTTCAGTGCTCCTTACCCCATG
GAGGTTATTTTGAAGTTGCTTCAATGGTCCCAAGGGAGGTTATGTTGCCCTGGATGATATTTCAATC
TCTCCTGTTCACTGCCAGAATCAGACAGAACTTCTGTTCACTGCCGTGGAAGCCAGCTGCAATTTTGA
CAAGATCTCTGCAACTTTTACCAAGATAAAGAAGGTCCAGGTTGGACCCGAGTGAAGTAAAACCAAC
ATGTATCGGGCTGGAGACCACACTACAGGCTTAGGGTATTACCTGCTAGCCAACACAAAGTTCACATCT
CAGCCTGGCTACATTGGAAGGCTCTATGGGCCCTCCCTACCAGGAACTTGCAGTATTGTCTGCGTTTT
CATTATGCCATCTATGGATTTTTAAAAATGAGTGACACCCTAGCAGTTTACATCTTTGAAGAGAACCAT
GTGGTTCAAGAGAAGATCTGGTCTGTGTTGGAGTCCCCAAGGGGTGTTTGGATGCAAGCTGAAATCACC
TTAAGAAGCCCATGCCTACCAAGGTGGTTTTCATGAGCCTATGCAAAAGTTTCTGGGACTGTGGGCTT
GTAGCCCTGGATGACATTACAATACAATTGGGAAGCTGCTCATCTTACAGAGAACTTCCACCTCCACCT
GGAGAGTGTACTTTGAGCAAGATGAATGTACATTTACTCAGGAGAAAAGAAACCGGAGCAGCTGGCAC
AGGAGGAGGGGAGAACTCCCACTTCTACACAGGACCAAGGGAGATCACACTACTGGGGTAGGCTAC
TACATGTACATTGAGGCCTCCCATATGGTGTATGGACAAAAGCACGCCTCTTGCCAGGCCTCTGCGA
GGAGTCTCTGAAAACACTGCTTGACCTTTTTCTACCACATGTATGGAGGGGCACTGGCCTGCTGAGT
GTTTATCTGAAAAGGAAGAAGACAGTGAAGAGTCCCTCTTATGGAGGAGAAGAGGTGAACAGAGCATT
TCTGGCTACGAGCACTGATTGAATACAGCTGTGAGAGGCAACACCAGATAATTTTTGAAGCCATTCGA
GGAGTATCAATAAGAAGTGATATTGCCATTGATGATGTTAAATTTAGGACAGGACCTGTGGAGAAATG
GAAGATACAACCTCAACAATCATCAGGATATTCTGAGGACTTAAATGAAATTGAGTATTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_153267

Insert Size: 2061 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_153267.4](#)

RefSeq Size: 3618 bp

RefSeq ORF: 2061 bp

Locus ID: 256691

UniProt ID: [Q7Z304](#)

Cytogenetics: 9q21.12

Domains: MAM
Protein Families: Secreted Protein
MW: 77.6 kDa