

Product datasheet for **SC318160**

AGAP6 (NM_001077665) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AGAP6 (NM_001077665) Human Untagged Clone
Tag:	Tag Free
Symbol:	AGAP6
Synonyms:	AGAP-6; CTGLF3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318160 representing NM_001077665.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGGAACATACTGACCTGTCGTGTGCACCCTAGCGTCAGCCTCGAGTTTGACCAGCAGCAGGGGTGCG
GTGTGTCCCTCTGAATCTGAGACCTATGAGGCAGGAGCTAGGGACAGGATGGCAGGAGCGCCCATGGCT
GCTGCTGTACAGCCTGCTGAGGTGACTGTTGAAGTTGGTGAGGACCTCCACATGCACCACGTTTCGTGAC
CGGGAGATGCCTGAAGCTTTGGAGTTTAACTTTCTGCCAATCCAGAGTCAAGCACAATATTCCAGAGG
AACTCTCAAACAGAAGCTTTGGAGTTTAACTTTCTGCCAATCCAGAGGCAAGCACAATATTCCAGAGG
AACTCTCAAACAGATGTTGTAGAAATAAGAAGAAGCAACTGTACAAACCATGTATCTGCTGTGCGTTTC
AGTCAACAATACAGCTTGTGTTGACAATATTCCTTGATGACAGCACAGCCATCCAGCATTATCTTACA
ATGACAATAATATCTGTGACCTTGAGATACCTCATCATATCACACAAAGAGATGCAGATAGAACTTTG
AGCATACCTGATGAACAGTTACACTCATTTCGGGTTTCCACCGTGCACATTATGAAGAAAAGAAATGGA
GGTGGGAGTTTAAATAACTATTCTCTCCATTCCATCGACTCCAGCACCAGCCAGGAGGACCCTCAG
TTCACTGTTCTCCCACTGCCAACACACCCACGCCGTTTGAAGCGGTCCATGCCCTGGTCCAACCTG
TTTACATCTGAGAAAGGAGTGACCCAGACAAAGAGAGGAAAGCCCCGGAGAATCATGCTGACACCATC
GGGAGCGGTAGAGCCATCCCCATTAAACAGGGCATGCTCTTAAAGCGAAGTGGGAAATGGCTGAAGACA
TGGAAAAAGAAATACGTACCCCTGTGTTCCAATGGCATGCTCACCTATTATTCAAGCTTAGGTGATTAT
ATGAAGAATATTCATAAAAAAGAGATTGACCTTCAGACATCTACCATCAAAGTCCAGGAAAGTGGCCA
TCCCTAGCCACATCGGCTGCACACCCATCTCCAGCTCTAAAAGCAATGGCCTATCCAAGGACATGGAC
ACCGGGCTGGGTGACTCCATATGCTTCAGCCCCAGTATCTCCAGCACCACCAGCCCCAAGCTCAACCCG
CCCCCTCTCCTCATGCTAATAAAAAAGAAACACCTAAAGAAGAAAAGACCAACAACCTTTATGATTGTG
TCTGCCACTGGCCAAACGTGGCACTTTGAAGCCACGACGTATGAGGAGCGGGATGCCTGGGTCCAAGCC
ATCCAGAGCCAGATCCTGGCCAGCTGCAGTCATGCGAGAGCAGTAAAAGCAAGTCCAGCTGACCAGC
CAGAGCGAGGCCATGGCCCTGCAGTCGATCCAAAACATGCGTGGGAACGCCCACTGTGTGGACTGTGAG
ACCCAGAATCCTAAGTGGCCAGTTTGAAGTGGGAGTCTCATGTATTGAATGCTCAGGTATCCAC
CGCAGTCTTGGCCCCACCTTTCCCGTGTGCGATCTCTGGAGCTGGATGACTGGCCAGTTGAGCTCAGG
AAGTTATGTATCTATTGTCAATGACCTAGCCAAACAGCATCTGGGAAGGGAGCAGCCAGGGGCAGACA
AAACCCTCAGAAAAGTCCACGAGGGAAGAGAAGGAACGGTGGATCCGTTCCAAATATGAGGAGAAGCTC
TTTCTGGCCCCACTACCCTGCAGTGTCCCTGGGCCAGCAGCTGCTGCGGGCCACCGCTGATGAG
GACCTGCAGACAGCCATCCTGCTGCTGGCACATGGCTCCTGTGAGGAGGTGAACGAGACCTGTGGGAG
GGAGACGGCTGCACGGCGCTCCATCTGGCCTGCCGAAGGGGAATGTGGTCTGGCGCAGCTCCTGATC
TGGTACGGGGTGGACGTATGGCCCGAGATGCCACGGGAACACAGCGCTGACCTACGCCCGGCAGGCC
TCCAGCCAGGAGTGCATCAACGTGCTTCTGCAGTACGGCTGCCCGCAGCAGTGTGTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001077665

Insert Size: 2061 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001077665.2</u>
RefSeq Size:	2772 bp
RefSeq ORF:	2061 bp
Locus ID:	414189
UniProt ID:	<u>Q5VW22</u>
Cytogenetics:	10q11.23
MW:	75.7 kDa
Gene Summary:	Putative GTPase-activating protein.[UniProtKB/Swiss-Prot Function]