

## Product datasheet for **SC318164**

### PTCD3 (NM\_017952) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PTCD3 (NM_017952) Human Untagged Clone
Tag:	Tag Free
Symbol:	PTCD3
Synonyms:	MRP-S39
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC318164 representing NM\_017952.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGCGGTTGTATCTGCTGTTGCTGGCTGGCCCTCCGCAGCAGGCTTGGCCAGCCGCTGACGGGTCCG
CGGGCGGGTTTGTGTGAACAGGCACGCAGCTGCAGATTTTATTCTGGTAGTGCAACCCTCTCAAAGTT
GAAGGAAGTGAATGAACAGGGATTGAAGAAGTAGTAATTCAAAAAAGAAAACCTGGGATAAAGTAGCC
GTTCTTACGGCACTTGCATCCACAGTAAACAGGGATACCACAGCTGTGCCTTATGTGTTTCAAGATGAT
CCTTACCTTATGCCAGCATCATCTTTGGAATCTCGTTCATTTTTACTGGCAAAGAAATCCGGGGAGAAT
GTGGCCAAGTTTATTATTAATTCATACCCAAATATTTTCAGAAGGACATAGCTGAACCTCATATACCG
TGTTTAAATGCCTGAGTACTTTGAACCTCAGATCAAAGACATAAGTGAAGCCGCCCTGAAGGAACGAAT
GAGCTCAGAAAAGTCAAAGCCTCTGTGGACATGTTTGATCAGTTTTGCAAGCAGGAACCACTGTGTCT
CTTGAAACAACAAATAGTCTCTGGATTATTGTGTTACTATGGTGACCAGGAGCCCTCAACTGATTAC
CATTTTCAACAACTGGACAGTCAAGCATTGGAAGAGGAAAAATGATGAGACATCTAGGAGGAAAGCT
GGTCATCAGTTTGGAGTTACATGGCGAGCAAAAAACAACGCTGAGAGAATCTTTCTCTAATGCCAGAG
AAAAATGAACATTCTATTGCACAATGATCCGAGGAATGGTGAAGCACCGAGCTTATGAGCAGGCATTA
AACTTGTACACTGAGTTACTAAACAACAGACTCCATGCTGATGTATACACATTTAATGCATTGATTGAA
GCAACAGTATGTGCGATAAATGAGAAAATTTGAGGAAAAATGGAGTAAAATACTGGAGCTGCTAAGACAC
ATGTTTGACAGAAAGGTGAAACCAATCTTCAGACTTTTAAATACCATTCTGAAATGTCTCCGAAGATTT
CATGTGTTTGAAGATCGCCAGCCTTACAGTTTTACGTGAAATGAAAGCCATTGGAATAGAACCCTCG
CTTGCAACATATCACCATATTATTCGCCTGTTTGATCAACCTGGAGACCCTTTAAAGAGATCATCCTTC
ATCATTTATGATATAATGAATGAATTAATGGAAAGAGATTTTCTCAAAGGACCCGGATGATGATAAG
TTTTTTCAGTCAGCCATGAGCATATGCTCATCTCTCAGAGATCTAGAACTTGCCTACCAAGTACATGGC
CTTTTAAAAACCGGAGACAACCTGGAATTCATTGGACCTGATCAACATCGTAATTTCTATTATTCCAAG
TTCTTCGATTTGATTTGTCTAATGGAACAAATTGATGTTACCTTGAAGTGGTATGAGGACCTGATACCT
TCAGCCTACTTTCCCACTCCCAACAATGATACATCTTCTCAAAGCATTGGATGTGGCCAAATCGGCTA
GAAGTGATTCTAAAATTTGAAAGATAGTAAAGAATATGGTCATACTTTCCGCAGTGACCTGAGAGAA
GAGATCCTGATGCTCATGGCAAGGGACAAGCACCCACCAGAGCTTCAGGTGGCATTGCTGACTGTGCT
GCTGATATCAAACTGCGTATGAAAGCCAACCCATCAGACAGACTGCTCAGGATTGGCCAGCCACCTCT
CTCAACTGTATAGCTATCCTCTTTTAAAGGGCTGGGAGAACTCAGGAAGCCTGAAAAATGTTGGGGCTT
TTCAGGAAGCATAATAAGATTCTAGAAGTGAGTTGCTGAATGAGCTTATGGACAGTGCAAAAGTGCT
AACAGCCCTTCCAGGCCATTGAAGTAGTAGAGCTGGCAAGTGCCTTCAGCTTACCTATTTGTGAGGGC
CTCACCCAGAGAGTAATGAGTGATTTTGAATCAACCAGGAACAAAAGGAAGCCCTAAGTAATCTAACT
GCATTGACCAGTGACAGTGATACTGACAGCAGCAGTGACAGCGCAGTGACACCAGTGAAGGCAAATGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_017952  
**Insert Size:** 2070 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_017952.5</a>
<b>RefSeq Size:</b>	6752 bp
<b>RefSeq ORF:</b>	2070 bp
<b>Locus ID:</b>	55037
<b>UniProt ID:</b>	<a href="#">Q96EY7</a>
<b>Cytogenetics:</b>	2p11.2
<b>Domains:</b>	PPR
<b>MW:</b>	78.5 kDa
<b>Gene Summary:</b>	Mitochondrial RNA-binding protein that has a role in mitochondrial translation. [UniProtKB/Swiss-Prot Function]