

Product datasheet for **SC318166**

ARMH3 (NM_024541) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ARMH3 (NM_024541) Human Untagged Clone
Tag:	Tag Free
Symbol:	ARMH3
Synonyms:	C10orf76
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_024541, the custom clone sequence may differ by one or more nucleotides

```

ATGGCACAGGTAGAGAAACGTGGAGGTTTGCTCCGGAAATCTTCAGCCTCCAAAAACCA
CTGAAGGAAAAAGTGGTGCTGATGTATGATGAGATCTTCATGACAGAGGACCCAGTAAG
TGCAGTCCTCGGTTTTGGGAGGAGCTCTTCTCATGAAGTGAATTTAGAGTACCTAGAA
GGCAAGCTGGAATCTCTTGATGGTGAGGAGTTAATGAAGATCAAGGACAATATTAATTGC
TTATTCCAACACTGCATCCAGGCTCTGGGAGAGGAGCATCCAATTCGGGTTGTCAATGCA
TTGCAGACCCCTGTGCGCACTCATTGAGGAGTCCATCAAAAGAATAAGTCTACCTCTGGG
TTTGACATTATCAACATGCTGATGGGCTTTGACAAGGCGGAGCTGTGCATGAAGAACTTG
ATGGAGAGTTTGGATTCAATTGCTTTGTGAGAAGGTTCTGAAAGTCTGAAGAGTTTATGT
CTGAAACTTCTCCTTTGCTTAGTGACCGTGACAGATAACATCAGCCAGAACACTATTCTC
GAGTATGTAATGATCAACAGCATATTTGAAGCAATTTTACAGATACTTTCCCATCCCCCA
AGTCGTAGGGAGCATGGGTATGATGCTGTCGTCCTCTTGCTTTGCTGGTGAAGTATAGA
AAATATGAGTCTGTGAATCCTTATATTGTGAAGCTGTCTATCGTGGATGATGAGGCCACA
CTCAATGGAATGGGACTTGTAATTGCTCAGGCTTTATCTGAGTACAACAGGCAGTATAAA
GACAAGGAAGAAGAACACCAAGTGGTTTTTCTCTGCTTTAACAATATGGTGGGCAGC
ATGTTTCATAGCAGATGCCATGAGAAAATCTCAGTACAAACCAATGAGGCCATTCTCTG
GCACTTTATGAAGCTGTTCAATTTAAATCGCAACTTCATCAGATATTAGCTCAGAGCCAT
CCAGAAATGGGCTTGGTGACGACCCCTGTGAGTCTGCTCTACAACCCAGTCACACCA
CTTGGGACCACACCGCCTTCTCTGATGTTATAAGTCTGTGGAAGTCCCACTGGATGCA
GATGTACAGACAGTAATCTACTGATAACCTTTTTAAAGTATAGCTCAATTGTCTATGCA
GACACCAAGATGAACACAGGCTTCACAGTGGCAAACTCTGTCTGATTATCCTTACATGT
ATTGCAGAGGATCAATATGCCAATGCATTTTGCATGATGACAACATGAATTTTCAGTA
AACTTACATAGAATGCCTATGAGACACAGGAAGAAGGCAGCAGACAAGAATCTTCCCTGC
CGTCCTTTAGTATGTGCACTACTGGACCTGATGGTAGAGTTTATTGTAACACACATGATG
AAGGAGTTTCTATGGATCTCTATATACGCTGCATCCAGGTAGTACACAACTGCTCTGC
TACCAGAAGAAGTGTGGGTACGCTGCATTACACCTGGCGGGAGCTCTGGTCAGCCTTG
ATAAATTTGCTGAAGTTCCTTATGTCAAATGAGACTGTACTTTTGGCCAAACACAACATT
TTTACATTAGCCCTTATGATTGTGAACCTATTTAATATGTTTATCACATATGGCGACACA
TTTCTGCCAACCCCGCAGCAGCTATGATGAACCTTACTATGAGATTATCCGCATGCACCAG
AGCTTTGACAACCTCTACTCCATGGTCCTGAGGCTTTCTACCAATGCAGGCCAGTGAAG
GAAGCAGCTAGCAAGGTGACCCATGCATTGGTTAATATCAGAGCCATCATCAACCACTTT
AACCCCAAAATTGAGTCCTACGCTGCTGTGAATCACATATCCCAACTGTGAGGAGCAG
GTGCTGGAGGTGGTGAGAGCCAACTATGACACGCTCACGCTGAAGCTGCAGGATGGCCTG
GACCAGTATGAGCGCTACTCAGAGCAGCACAAGGAAGCTGCCTTCTTCAAAGAGCTGGTT
CGATCCATTAGCACCAACGTCCGGAGAAACCTGGCCTTCCACACACTCAGCCAAGAAGTC
CTGCTCAAGGAGTTCTCCACTATCTCC

```

Restriction Sites: Please inquire

ACCN: NM_024541

OTI Disclaimer:	Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery. The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_024541.2 , NP_078817.2
RefSeq Size:	4101 bp
RefSeq ORF:	2070 bp
Locus ID:	79591
UniProt ID:	Q5T2E6
Cytogenetics:	10q24.32
Gene Summary:	Involved in GBF1 recruitment, Golgi maintenance and protein secretion.[UniProtKB/Swiss-Prot Function]