

Product datasheet for **SC318167**

ZNF185 (NM_007150) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF185 (NM_007150) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF185
Synonyms:	SCCELL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC318167 representing NM_007150.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGTATCTCAGCTCTTGGAGGCCGCACCAAAGGGAAGCCTCTGCCACCAGCGAGGAGGAGCGCAAT
AACGTTCTCAAGCAGATGAAAGTGCGAACCACGCTGAAGGGGGACAAGAGCTGGATTACCAAGCAGGAT
GAATCGGAGGGTCGCACCATAGAGCTGCCCTCAGGCCGGAGTCGCGCCACATCCTTTTCATCAGCTGGG
GAGGTTCCGAAGCCTAGGCCTCCGAGCACAAGGGCTCCCACTGGCTACATCATCCGGGGAGTGTTACC
AAGCCATAGACAGCTCTCCAGCCCCAGCAGAGTTCCCAAGGCCAACGGGACTCCAAAAAGTGCT
GCCAGTCTGGTGAGAACAGCTAACGCTGGTCTCCCGCCCCCTCCTCTGGCTACAAGATGACCACT
GAGGATTACAAGAAGCTGGCACCTACAATATCAGGCGCAGCTCTACATCAGGGGACACCGAGGAGGAG
GAGGAGGAGGAGTGGTCCATTCTCTCAGATGAACAGAAACGGAGGTGAGAGGCTGCAAGCGGTGTT
CTGAGGAGGACAGCTCCCGGGAGCACTCTACGTCTGTGAGCGCCAGAAGAGCACTGGTCTACC
CAGGAGACACAGGCACCGTTTATCGCGAAGAGGGTGGAGGTGGTGAAGAGGACGGGCCCTCTGAGAAG
AGCCAGGACCCACCTGCTCTGGCAAGATCCACTCTGGCTCAAACAGTGGCGATGGAGGCAGGACCAAAA
GCGTCTCGGGCAATTTGGATCGAGTGCCTGCCAAGTATGCCTAGCCCCGCTGGGAGCCAGGAGCTCAGC
TCAAGAGGTGAGGAAATTTGTCGGCTGCAGATCCTGACACCCAGGGCAGGACTCCGCCTGGTGGCCCA
GACGTGGAAGGCATGAGGTCTTCCCGAGCAACAAGACAAGGAGGCCCTGCTCCAGAGAGCTCCAG
AGGGACTTGGCTGGTGGAGAGGCTTTCAGGGCCCCAACACAGATGCTGCAAGGTCAAGTGCACAGTTG
AGTGATGGCAATGTGGGATCCGGAGCCACGGGCTCCCGCCCTGAAGGCTTGGCTGCAGTAGACATCGGC
TCCGAGAGAGGAAGCTCCAGTGCCACTTCACTCTGCTGTCCCTGCTGATAGGAAGAGCAACAGCACA
GCGCCCAAGGAGGATGCAAAGGCAGACCCAAAGGGGGCCTTGGCTGATTATGAGGGGAAGGATGTGGCC
ACCAGGTCGGAGAGGCCTGGCAGGAGAGGCCTGGAGCTCCAAGAGGTGGCCAAGGAGACCCAGCTGTA
CCCGCTCAGCAACCTGCAGATCCAGCACCCAGAGCGGCAGAGCAGCCCAAGCGGATCTGAGCAACTT
GTCAGACGAGAGAGTTGTGGCAGCAGCGTGTGACTGATTTTGGGGGAAGGATGTGGCCACCAAGGTC
GGAGAGGCCTGGCAGGACAGGCCTGGAGCCCCAAGAGGTGGCCAAGGAGACCCAGCTGTACCCACTCAG
CAACCTGCAGATCCAGTACCCAGAACAGCAGAACAGCCCAAGCGGATCTGAGCAATTCGTGAGACGA
GAGAGCTGCACCAGCAGGGTGGAGAGCCCTCGAGCTGCATGGTCACTGTTACTGCACTGCCACATCT
GAGCAGCCTCACATTTATATTCCAGCCCCGCAAGTGAATTGGACTCCAGCTCTACCACCAAGGGATT
CTCTTCGTGAAGGAGTACGTGAATGCTAGTGAAGTGTCTTCTGGGAAGCCAGTATCTGCAGCTATAGC
AACGTCAGCAGCATTGAGGACTCATTCCGCATGGAGAAGAAGCCTCCATGTGGCAGCACTCCATACTCT
GAGAGGACAACCTGGAGGGATCTGTACTTACTGCAACCGTGAGATCCGAGACTGTCCAAAGATTACCTA
GAACATCTTGGTATCTGCTGCCATGAATATTGCTTTAAGTGTGGGATTTGCAGTAAACCGATGGCGGAT
CTCCTGGATCAGATCTTCATTACCGTGACACCATTCACTGTGGGAAATGCTATGAGAAGCTCTTTAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_007150

Insert Size: 2070 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_007150.3](#)

RefSeq Size: 4344 bp

RefSeq ORF: 2070 bp

Locus ID: 7739

UniProt ID: [O15231](#)

Cytogenetics: Xq28

Domains: LIM

MW: 73.5 kDa

Gene Summary: Zinc-finger proteins bind nucleic acids and play important roles in various cellular functions, including cell proliferation, differentiation, and apoptosis. This gene encodes a LIM-domain zinc finger protein. The LIM domain is composed of two contiguous zinc finger domains, separated by a two-amino acid residue hydrophobic linker. The LIM domain mediates protein:protein interactions. Multiple alternatively spliced transcript variants encoding different isoforms have been identified.[provided by RefSeq, May 2010]
Transcript Variant: This variant (4) lacks an in-frame exon in the CDS, as compared to variant 1. The resulting isoform (4) lacks an internal segment, as compared to isoform 1.