

Product datasheet for **SC318175**

MICALCL (NM_032867) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MICALCL (NM_032867) Human Untagged Clone
Tag:	Tag Free
Symbol:	MICALCL
Synonyms:	Ebitein1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318175 representing NM_032867.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGTCACCTCCTAAGGACCCTTCTCCTTCTTCTCCTTCTCCTTCCATTTCATCTTCCCCA
CCATCTTCTTCTCAACCAGTGTCTGGGAATGCTCCAGATGGTTCTCCCGCCTCAGATGACAGCT
TCTGAGCCCTCTCGCAAGTCTCGAGAGGTCATCCAAGTCTCCACCCAAACTTCGGAGGCGAGCC
GTAGCCCAAGGAGCACCCAGGAAATCCCTGTATCTGCCTATCACCCAAAGCCAGAGTGGGAGAG
TACTGCCTGGTGAACCGTGGTGAAGATGGCCTCTCAGACCCTGCAGAGATGACTTCTGATGAGTGCCAG
CCAGCAGAGGCCCTCTGGGGACATCGGAAGCAACCACAGAGACCCACCCCACTGGGGGAAGGAC
AGGAGCTGGACAGGGCAAGAGCTATCTCCCTGGCTGGAGAAGACCGGAAAAAGGGAGTACTGGAGCC
AGGAAGGAAGAAGAGGGAGGGCCAGTGTGGTAAAGGAGAAGTTGGCCTGAAGAAGTTAGTCTCACT
CAGGAGCAGAAGACCATGTTGTTGGATTGGAATGACTCCATCCCTGAGAGTGTGCACCTCAAAGCTGGG
GAGCGAATTTCCAGAAAAGTGTGAGAATGGTAGAGGAGCCGTGTGCTAAAACAGTCCGCCCTCTG
CTGCTCCCTAGGGCAGCAGGAGAGCCCTGCCAACCCAGAGAGGGGCCAGGAGAAGATGGGGACCCT
GCGGAACAAGCTCAAGGGGAGCGAAACGTGCCTCCACCCAAGTCCCACTGCGGCTCATAGCCAATGCC
ATCCGAAGGTCTCTAGAGCCCTCTTTCCAACCTGAAGGCGGGGAAGAAGGCTGGGCCAAGCAAGAA
TCCAAAATTTGCCCGCACAGGCCTGCACTCGCTCATTAGCCTTCGAAAACCAATTCATAAAGAC
GGGACCAGCATTCCCTGGGAGAAACAGTCTCAGCCTTTAGCCCTCTGACCCTGCCCTCCGCACC
CACAGTTTGCCAATCGGCCATCCAAGTCTTCTGCACTTAGGTCCCAACCTGCAGCAAGATTGAA
GATGTCCCCACACTCTCGAAAAGTGAAGTTGCAAGAGAAGTCCAGATGCTTCAAGCCTCAAAG
AAAAGAATCTCACTTTTTCTCCTCCTCAGACTCAAAGCAAATCTTTGAGAGTTTCTCAAGAATCC
AGACAAAAGAAAGGACATCAGGGACCTCTTTGGCAGCCCAAGAGGAAGGTGCTGCCTGAAGATAGTGCG
CAGGCCCTGGAGAAGCTGCTGCAGCCTTTCAAAGCACCTCCCTGCGCCAGGCAGCTCCTCCTCCT
CCTCCTCCTCCTCCTCCTCCTCCTACAGCGGGAGGTGCAGACTCCAAGAAGTTCCTCCTCAGA
GCACAGGTAACAGAGGCTTCTCTTCTGCTCTTCAACCTCCTCCTCCTGAGATGAAGAATTTGAT
CCCCAGCTTCTTGCAGTTAAAGGAGAAGAAGCACTTAGAAGAAGAAAGAGCTAGAAAAGCAATG
AAGCAGTTGGTTAAGCAAGAAGAATTGAAAAGACTCTATAAGGCTCAGGCCATCCAGAGGCAGCTGGAG
GAGGTGGAGGAGCGGCAGAGGGCTTCTGAGATCCAGGGTGTGAGGCTGGAGAAGGCGTTGCGAGGAGAA
GCAGATTCAGGCACACAGGATGAAGCACAGCTTTTGCAGGAATGGTTAAGCTGTTCTGGAGAAGAT
AAATTAATGCGATATGAGTCGGAGCTCCTAATCATGGCCAGGAAGTGAATAGAAAGATCATCAAAGC
AGACTGGAGCAGAACTGAGAGAGAAAATGCTCAAGGAGGAGAGCCAGAAAGATGAGAAGGATCTAAAC
GAAGAGCAAGAAGTATCACCGAGCTGATGCAAGTATTGAGCAAAGGGACAAACTCGTCGATTCCTTA
GAGGAACAACGCATCAGAGAAAAAGCCGAGGACCAGCACTTTGAAAGCTTCGTATTCTCCAGAGGCTGT
CAGCTGAGCAGGACTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_032867

Insert Size: 2088 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_032867.2
RefSeq Size:	3044 bp
RefSeq ORF:	2088 bp
Locus ID:	84953
UniProt ID:	Q6ZW33
Cytogenetics:	11p15.3
Protein Families:	Protease
MW:	77.3 kDa
Gene Summary:	May cooperate with MAPK1/ERK2 via an intracellular signal transduction pathway in the morphogenetic development of round spermatids to spermatozoa. May act as Rab effector protein and play a role in vesicle trafficking.[UniProtKB/Swiss-Prot Function]