

Product datasheet for **SC318177**

TMEM168 (NM_022484) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM168 (NM_022484) Human Untagged Clone
Tag:	Tag Free
Symbol:	TMEM168
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318177 representing NM_022484.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGTGTAATCACTGCGTTATTGCTTTAGTCATTGTCTCTATTTAGCAATGACAAGACTGGAAGAAGTA
AATAGAGAAGTGAACATGCATTCTTCAGTCCGGTATCTTGCTATTTAGCCAGAATCAATTTATTGGTT
GCTATATGCTTAGGTCTATACGTAAGATGGGAAAAACAGCAAATTCCTTAATTTTGGTAATTTTATT
CTTGGTCTTTTGTCTTGAATCGCCAGCATACTCTATTACTATTTTCAATGGAAGCAGCAAGTTTA
AGTCTCTCCAATCTTTGGTTTGGATTCTTGCTTGGCTCCTATGTTTTCTTGATAATTCATCCTTTAAA
AATGATGTAAAAGAAGATCAACCAAATATTTGCTTCTAACATCCATAGTGTTAAGGATATTGTGCTCT
CTGGTGGAGAGAATTTCTGGTTATGTCGTCATCGGCCCACTTTACTAACCACAGTTGAATTTCTGGAG
CTTGTTGGATTGGCATTGCCAGCACAACTATGTTGGTGGAGAAGTCTCTGAGTGCATTTTGTGTTT
GTAGCTCTGGCTATGCTGATTATTGATCTGAGAATGAAATCTTCTTAGCTATTCCAAACCTAGTTATT
TTTGCAGTTTTGTATTTTTTCTCATTGGAACCTCCAAAAATCCGATTGCTTTTGCCTGTTTTTTT
ATTTGCCTGATAACTGATCCTTTCTTGACATTTATTTTAGTGGACTTTCACTAAGTGAAGATGAAAA
CCCTTTTTGTACCGTGAAGAATTTGCAGAAGACTTTCAGTCGTTTTTGTGGAATGATTGAGCTTACA
TTTTTTATTCTTTCCGCATTCAAACCTAGAGACACTCACCTCTGGTATTTTGTAAACCTGGCTTTTCC
ATTTTTGGAATTTCTGGATGATTTGTCATATTATTTTCTTTAACTCTTTGGGATTCATACCAAAA
TTAAATGACTGCCATAAAGTATATTTTACTCACAGGACAGATTACAATAGCCTTGATAGAATCATGGCA
TCCAAAGGGATGCGCCATTTTGTCTGATTTCAGAGCAGTTGGTGTCTTTAGTCTTCTTGAACAGCG
ATTTTGGGAGCAGTTTCTGGCAGCCAACAAATGGAATTTTCTTGAGCATGTTTCTAATCGTTTTGCCA
TTGGAATCCATGGCTCATGGGCTCTTCCATGAATTGGTAACTGTTTAGGAGGAACATCTGTTGGATAT
GCTATTGTGATTCCACCACTTCTGCAGTCCTGATGGTCAGCCAACACTGCTTCCCCCAGAACATGTA
CAGGAGTTAAATTTGAGGTCTACTGGCATGCTCAATGCTATCCAAAGATTTTTTGCATATCATATGATT
GAGACCTATGGATGTGACTATTCACAAGTGGACTGTCATTTGATACTCTGCATTCCAACTAAAAGCT
TTCCTCGAACTTCGGACAGTGGATGGACCCAGACATGATACGTATATTTTGTATTACAGTGGGCACACC
CATGGTACAGGAGAGTGGGCTCTAGCAGGTGGAGATACACTACGCCTTGACACACTTATAGAATGGTGG
AGAGAAAAGAATGGTTCTTTTGTCTCCCGCTTATTATCGTATTAGACAGCGAAAATCAACCCCTTGG
GTGAAAGAAGTGAGGAAAATTAATGACCAGTATATTGCAGTGCAAGGAGCAGAGTTGATAAAACAGTA
GATATTGAAGAAGCTGACCCGCCACAGTAGGTGACTTTACAAAAGACTGGGTAGAATATAACTGCAAC
TCCAGTAATAACATCTGCTGGACTGAAAAGGGACGCACAGTGAAGCAGTATATGGTGTGTCAAAACGG
TGGAGTGACTACACTCTGCATTTGCCAACGGGAAGCGATGTGGCCAAGCACTGGATGTTACACTTTCCT
CGTATTACATATCCCCTAGTGCATTTGGCAAATTTGGTTATGCGGTCTGAACCTTTTTTGGATCTGCAAA
ACTTGTTTTAGGTGCTTGAAAAGATTAAAAATGAGTTGGTTTCTTCTACTGTGCTGGACACAGGACAA
GGCTTCAAACCTTGTCAAATCTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_022484

Insert Size: 2094 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_022484.5</u>
RefSeq Size:	7317 bp
RefSeq ORF:	2094 bp
Locus ID:	64418
UniProt ID:	<u>Q9H0V1</u>
Cytogenetics:	7q31.1
Protein Families:	Transmembrane
MW:	79.8 kDa