

Product datasheet for **SC318189**

RUFY1 (NM_025158) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RUFY1 (NM_025158) Human Untagged Clone
Tag:	Tag Free
Symbol:	RUFY1
Synonyms:	RABIP4; ZFYVE12
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318189 representing NM_025158.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGCCGACCCGGGAAGGCGGCTGCGCTGCTGGCGGGGGCGGAGCTGGAGCCGGAGCTGGAGCCGGGG
CCGGGGCCCGGGTACGCGCTTGAGCCGGGAGAAGAGTTTGAGATCGTGGACCGAAGCCAGCTGCCGGC
CCAGGGCAGCTGCGGAGCGCAACGAGGCCGGCGGGCGGAGGGCTGGTGGCGCCCATCCTGACCCCTG
GCACGCAGGGCCACCGGGAACCTGTCGGCGAGCTGCGGGAGCGCCTGCGCGCGGGCCGGGGCTGGGC
GGCGGGGACAGCGGGGACGGCACGGCGCGCAGCTTCTAAGTGCCAGATGATGGAGGAGCGTGCCAAC
CTGATGCACATGATGAACTCAGCATCAAGGTGTTGCTCCAGTCGGCTCTGAGCCTGGGCCCGAGCCTG
GATGCGGACCATGCCCTTGCAGCAGTCTTTGTAGTGATGGAGCACTGCCTCAAACATGGGCTGAAA
GTAAAGAAGATTTTATTGGCCAAAATAAATCATTCTTTGGTCTTTGGAGCTGGTGAGAACTTTGT
CCAGAAGCATCAGATATAGCGACTAGTGCAGAAATCTCCAGAATTAAGACAGCTGTGGGAAGAGGC
CGAGCGTGGCTTTATCTGCACTCATGAAAAGAACTGGCAGATTATCTGAAAGTGCTTATAGACAAT
AAACATCTCTTAAGCGAGTTCTATGAGCCTGAGGCTTTAATGATGGAGGAAGAAGGGATGGTGATTGT
GGTCTGCTGGTGGGACTCAATGTTCTCGATGCCAATCTCTGCTTGAAAGGAGAAGACTTGGATTCTCAG
GTTGGAGTAATAGATTTTTCCCTCTACCTTAAGGATGTGCAGGATCTTGATGGTGGAAGGAGCATGAA
AGAATTACTGATGTCCTTGATCAAAAAAATATGTGGAAGAAGCTTAACCGGCACTTGAGCTGCACAGTT
GGGGATCTTCAAACCAAGATAGATGGCTTGAAAAAGACTAACTCAAAGCTTCAAGAAGAGCTTTCAGCT
GCAACAGACCGAATTTGCTCACTTCAAGAAGAAGCAGCAGTTAAGAGAACAAAATGAATTAATTGAA
GAAAGAAGTAAAAGAGTGTAGAGATAACAAAACAGGATACCAAAGTTGAGCTGGAGACTTACAAGCAA
ACTCGGCAAGTCTGGATGAAATGTACAGTGATGTGGAAGCAGCTAAAAGAGGAGAAGAAAGTCCGG
TTGGAACCTGAAAAAGAAGTGGAGTTACAATTTGGAATGAAAACCGAAAATGGAAAATTGCAATGAAGTTA
CTGAAAAAGGACACCCACGAGAAGCAGGACACACTAGTTGCCCTCCGCCAGCAGCTGGAAGAAGTAAA
GCGATTAATTTACAGATGTTTCAAAAGCTCAGAATGCAGAGAGCAGTTTGCAGCAGAAGATGAAGCC
ATCACATCTTTGAAGGAAAAACCAACCAAGTTATGTCCAGCATGAAACAAATGGAAGAAAGGTTGCAG
CACTCGGAGCGGGCAGGAGGGGGCTGAGGAGCGGAGCCACAAGCTGCAGCAGGAGCTGGCGGGAGG
ATCGGCGCCCTGCAGCTGCAGCTCTCCAGCTGCACGAGCAATGCTCAAGCCTGGAGAAAGAATTGAAA
TCAGAAAAAGAGCAAAGACAGGCTCTTCAGCGCAATTACAGCAGAGAAAGACACTTCTCTCTACTC
AGGATGGAGCTGCAACAAGTGAAGGACTGAAAAAGGAGTTGCGGAGCTTCAAGACGAGAAAGGCAGAG
CTGAGAAAGATCTGCGAGGAGCAGGAACAGCCCTCCAGGAAATGGGCCTGCACCTCAGCCAGTCCAAG
CTGAAGATGGAAGATATAAAGAAGTGAACCAGGCACTGAAGGGCCACGCCTGGTGAAAGATGACGAA
GCGACACACTGTAGGCAGTGTGAGAAGGAGTTCTCCATTTCCCGGAGAAAGCACCCTGCCGGAAGTGT
GGCCACATCTTCTGCAACACCTGCTCCAGCAACGAGCTGGCCCTGCCCTCTACCCCAAGCCGGTGCGA
GTGTGCGACAGCTGCCACACCCTGCTCTGCAGCGCTGCTCTCCACGGCCTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_025158
Insert Size: 2127 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_025158.4</u>
RefSeq Size:	2660 bp
RefSeq ORF:	2127 bp
Locus ID:	80230
UniProt ID:	<u>Q96T51</u>
Cytogenetics:	5q35.3
Domains:	FYVE, RUN
Protein Pathways:	Endocytosis
MW:	79.8 kDa
Gene Summary:	<p>This gene encodes a protein that contains a RUN domain and a FYVE-type zinc finger domain. The encoded protein binds to phosphatidylinositol-3-phosphate (PI3P) and plays a role in early endosomal trafficking, tethering and fusion through interactions with small GTPases including Rab4, Rab5 and Rab14. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Dec 2011]</p> <p>Transcript Variant: This variant (1) encodes the longer isoform (a, also known as rabip4').</p>