

Product datasheet for **SC318190**

FBXO43 (NM_001029860) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FBXO43 (NM_001029860) Human Untagged Clone
Tag:	Tag Free
Symbol:	FBXO43
Synonyms:	EMI2; ERP1; FBX43
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC318190 representing NM_001029860.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGTTTTAAAGACAAAGATGAGAGAATTTCTGTTTGAAGCCTACGTAACCTTTGACATCTAAGAGC
TCAAGATTTACTGATGAAACAGAGATTTGAAGATGTCGCAAAGGCACTCAGGTCAAGCTGGCACTGAA
GCAGGAAATGGGGCGGACTCTCCTCCAATTGTCAACTCCAAGTACTCCACCTTCAGAGATTTTGTTC
ACATCTTCATTTCAAGATAGTGGCTACAATGAGTTAAAATCTTGTAGCTTTGATAATATAGATAAAGAA
TATCTTGAAAGAAAGAAAAAGGCCAACATTACTCTATGAGCACCTGAAACTCAGGCCTGGGCTTA
ACACATCCTTTAGAATCTCCCACTCAAAAAAGAAATGTATCTTGCTAGAAAGGAAAAGGATAAAACC
CCAGAATTTGTAAACACCTAAAATCAGTGGGAAAAATGTTTACCTCGCAGAAGGTTGAATGTATCT
TTCGCTCTTCTAAAAGGGGACTTTGAATCACAAAATAGTTCTTTAGAAAGTAGTATAAGCCAAGTTATC
AACTTAGAAAAAATATCCAAGCAGTGCTTCAGGTTTTCCAGGGCAAATAATTTAGCCCTTTAGTT
ACTAGCACTTTAAAACAGAAGAAGTACTTCATGCAGTCAAAAATGAGGCTTAATTTTTCTCAGCAA
AAGACTTCCACAATTGATGATTCCAAAGATGATTGTAGCCTATTTGAAGTTGAATGTATATCTCCAATT
CAGGGCAATAATTTAAAGACTCTATCACACATGACTTTAGTGATAGCAGTTTATGCATTATGATGAG
AATGCATGTCCAGAGCTCCTGGGCTCCTGTTAGTGGAACTTGTGGAAACAGATGAGGACATATTT
GTGACTCCGATAAGTAATCTTGTGGCAACATTAGATTTAACGCAAGTCAAATACCTTTCTCCTTCACCT
GAAGTGAGAGGCAGTATTTCAACGCCTGAAGACAGTGGTTTTAACTCACTTAGCTTGAGAAATCAGAA
GATTCCTGTCTGACCAGGAGGTTCTTTCAAGAACTACTGCAGAAACATAAGGGGACTCCCAAAGTT
GGGACACCATAAGAAAGACAAGACATCTTGAAGGTGAGAAAGACTGTCCACCTTCGGAACAAAAGC
TCGCAGTCAGAGACAGAAGAGGAAAAGCAGATTTGCCACCTGACTCTGAAAAAGAGCAGCAGCTGCT
TCTGCCATCTCAGAGGGTCAGCTGAGTAGTGATGAGAGTGGGGATTTGACCTTTAGCTTAAAGAATTTA
TCAAAGACCCAGCCTTGAATTGGTACATGAGCTGTTTCAAAAAGCAAGAGGAAAAGATTACAGGAA
AATAGTGGACATGAATTCTTAGAGCAAGGGGATGGGGAGAAAATAGCTGTACTGCAGTGTACTTGCA
GGACTGATCGGCAAGAAAATGGGTATAGAAAACTGGACATCTTAACAGAATTAATAATAGAAATTTA
AAGCATATCTTGCTATGGTTTTAGAGTCTTGACCGCAGAGAGCCTATGCAGTGTGGAAAGTAAAGC
AGAAATGGCGTGAATTTGTTGTTCAAGATAAAAATGCAATCGGAGGAGGAAATTTATATCACACAA
CTGAAAACAGATTCTGAGGGGCTGTATTAATGTGAGGATGCTGCCACTCGGCTCCAGCTTTTAAAT
CGCTCAGCTTTAAGATCTGTGCAGGCACAGGCTAGGATACCTGGTTCTCAGAGAGAGCAAGGGTCAACA
TTATCTCCCTGGGAGAAGTTTTGACACCTTAGCAAGCTCTTCTGTTACTCACTTAAGTAGTAAACAG
GAAGAATATGTTAAGGTTGCCAAAACACTTTTTACTGATGAAGCATTAAAACCTTGCCCAAGGTGCCAG
TCCCTGCTAAGTACCAGCCATATAAGAAAAGGGGACTGTGTAGCCGAACAGCCTGTGGTTTTGACTTT
TGTGTGTATGTCTGTGTGCTTATCATGGGTCTGAAGAAATGTAGTAGAGGAGCAGCAAAGCCAAGAAAT
AGAAAAGATGCTCTCCAGGAAGTGCCAGAGTAAGCGGAATTTAAAACGCCTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_001029860

Insert Size: 2127 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001029860.3](#)

RefSeq Size: 3305 bp

RefSeq ORF: 2127 bp

Locus ID: 286151

UniProt ID: [Q4G163](#)

Cytogenetics: 8q22.2

Protein Pathways: Oocyte meiosis

MW: 78.4 kDa

Gene Summary: Members of the F-box protein family, such as FBXO43, are characterized by an approximately 40-amino acid F-box motif. SCF complexes, formed by SKP1 (MIM 601434), cullin (see CUL1; MIM 603134), and F-box proteins, act as protein-ubiquitin ligases. F-box proteins interact with SKP1 through the F box, and they interact with ubiquitination targets through other protein interaction domains (Jin et al., 2004 [PubMed 15520277]).[supplied by OMIM, Mar 2008]
Transcript Variant: This variant (2) represents the shorter transcript but encodes the supported protein.