

Product datasheet for **SC318192**

KRI1 (NM_023008) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KRI1 (NM_023008) Human Untagged Clone
Tag:	Tag Free
Symbol:	KRI1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318192 representing NM_023008.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGGCCACAGAACCGCCATGCCGAACCGCGCGGTTCGTGCGAGCTGCGGGTGAACCGCGCGTTTGCC
GCGCGGTACAACCGCTACCGGGAGCGCGAGGAAGTGCAGCGGCTGAAGGATCGCTACGGGACCGAGAC
AGCAGCAGCGACTCCAGCTCCGAGTCCGACTCAAGCGACGAGCGGTGAATTTGATCCCCAGCAGGAG
CGGGACTTTTACAAAACGCTCTCCTTGTGAAGAAGAAGGACCCCGCATTTATCAGAAAGATGCCACC
TTCTATAACAGAACAGCATCGTCATCAGACAGTGAAGGAGACCCAGAAGCCTTGGAAGAAGCAGAGAAA
GTGCGGCCCATGTACCTGAAGGACTACGAGAGGAAGGTTATCTTGGAAGAAGCAGGCAAAATATGTTGAT
GAGGAGAACTCAGACGGGAGACTTCCAATCACAGACTCCAGGAGACATCGTCGCAAGTTATGTGGAG
GAACAGAAACAGCTCAAGGAAAGCTTCCGGGCATTTGTGAGGACAGTGAAGACGAGGACGGCGCTGGG
GAGGGCGGCTCCAGTTTGTGCGAGAACTGCCAAAACAGGACGAGAGAGGCCAGGAGGAGGCCGAC
TACATCGAGTGGCTGAAGGGACAGAAAGAGATTCGGAACCCAGATTCCCTGAAGGAACTGACGCATCTC
AAGGAATACTGGAACGACCTGAGTTGGATGAAGGGGAGCGGTTCTGCGGGATTACATCCTCAACAAA
CGCTATGAGGAGGAGGAAGAGGAGGAGGAAGATGAAGAGGAAATGGAGGAAGAGGAGGGGTCCACGGT
CCCCCAGTCCAGCTGGCTGTGGACGACTCCTCAGACGAAGGGGAGCTGTTTCTGAAGAAACAGGAGGAC
TTTGAACAGAAGTACAATTTCCGTTTCGAGGAGCCGGACTCAGCATCGGTCAAGACCTACCCACGCGAC
ATCGCGTCTCCGTGCGCCGTAAAGATGAGCGCAGAAAGGAGAAGAGGGAAGAGACTCGGGAGCGAAAG
AAGAGGGAGAAAGCAAAGAAGCAGGAAGAGCTCAAGCAGCTGAAGAACCTGAAGAGGAAGGAGATTCTG
GCCAAGCTGGAGAAGCTGCGGAAAGTAACAGGCAACGAGATGCTGGGCCTCGAGGAGGGGACCTTGAA
GACGACTTCGACCCTGCCAGCACGACGACTCATGCAGAAAGTCTTTGGGACGAGTACTACGGGGCC
GTGGAGGAGGAGAAGCCACAATTTGAGGAAGAAGAAGGGCTTGAAGACGACTGGAAGTGGGACACGTGG
GACGGGCCTGAGCAGGAGGAGACTGGAGCCAGCAGGAGCTGCACTGTGAGGACCCCACTTCAACATG
GACGCCGACTACGACCCAGCCAGCCGAGGAAGAAAAAGCGGAGGCCCTTGACGGGCAAGAAGAAG
CGCAAGTCGCCCTTCGCCGCGCCGTGGGGCAGGAGAAGCCCGTGTGTAACCCGGGACAAAGACGTTT
GAGGAGTACCTGGATGAGTATTACCGGCTGGACTACGAGGACATCATCGACGACCTGCCCTGTGCTTC
AAGTACCGCACAGTGGTGCCCTGTGACTTTGGCCTCAGCACTGAGGAGATCCTCGCTGCTGACGATAAG
GAGCTGAACCGGTGGTGCTCCCTAAAGAAGACCTGCATGTACAGGTGAGGACGAGGAGCTGCGGGAC
AAGCGGGCGTACAGCCAGAAGGCCAGAACTCATGAAAAAGCGGAGGTCTTCAAGTCACTCTGCCGA
GAAGAGGCAGAGACACCTGCGGAAGCCACAGGGAAGCCACAGAGAGATGAAGCCGCCACAGAGGCAG
CTGCCAGCCCTTGATGGCAGCTTGATGGGGCCGAGAGTCCCCCAGCACAGGAAGAGGAAGCCCTGTA
TCACCCCAACAAGAAGCCAGCCCCCAGAAAGCGGAGGAGGCCAAGAAGGCACGGCTGCTGGGCCCACT
GTGATGCTTGGTGGATGCGAGTTTCAGCCGCCAGAGACTGCAGGCTTTGGCCTCAACCCAAACGGCTG
CACTTCCGCCAGCTGGGCCGCGAGCGGAGGAAACAACAGGGGCCCAAGAACAGCTCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_023008

Insert Size: 2130 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_023008.3</u>
RefSeq Size:	3023 bp
RefSeq ORF:	2130 bp
Locus ID:	65095
UniProt ID:	<u>Q8N9T8</u>
Cytogenetics:	19p13.2
Domains:	Krr1
MW:	83.3 kDa
Gene Summary:	This gene overlaps with the gene for cysteine endopeptidase AUT-like 4 in a head-to-tail orientation. [provided by RefSeq, Jul 2008]