

Product datasheet for **SC318194**

INTS10 (NM_018142) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	INTS10 (NM_018142) Human Untagged Clone
Tag:	Tag Free
Symbol:	INTS10
Synonyms:	C8orf35; INT10
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318194 representing NM_018142.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGCTCTGCCAGGGGACTGCGAGTTCTGGTGCAGCGAGCCCGGAGTTGGTGCCGAAGACCTGTGG
GCAGCCAAGGCGTGGCTGATCACGGCCCGCAGCCTCTACCCGGCAGACTTTAACATCCAGTATGAGATG
TACACCATCGAGCGGAATGCAGAGCGGACCGCCACCGCCGGGAGGCTGCTGTACGACATGTTTGTGAAT
TTCCAGACCAGCCGGTGGTGTGAGAGAAATCAGCATTATTACATCAGCATTAAAGAACGATTCACAG
GACAAACAAACCAATTTTAAAGAAGTTTATTTGAAACTCTTCTGGTCCGGTCCAGTGTGAAATGTTA
CTAAAGGTCACGAACAATGCTTCAACAGTTAGAACGATCAGAAATGTTGCTTCTACTTTTGGGCGC
TTCCCTGAAACGGTGGTGCAGCATGGGTTGGCCTTGGGAGGCACTATTAGAGGCTGAACTATTGAA
GAACAAGAATCTCCAGTGAAGTCTTTAGAAAATTATTTGTTGTGATGCTTCTCTAATAATTAAC
AACCATGATGTTGATTACCTGCCAATTTATTGTATAAGTACTTGAACAAAGCAGCTGAATTTTATATC
AATTATGTCACTAGGTCTACTCAATAGAAAATCAGCATCAAGGCGCCAGGATACATCTGATTTAATG
TCACCTAGCAAACGTAGCTCTCAGAAGTACATAATAGAAGGGCTGACGGAATAATCATCCCAGATCGTG
GACCTTGGGAGAGGTTGTTAAGATTTGAATGTTGTTGGAATGAGATGTGAATGGCAGATGGATAAA
GGAAGACGAAGCTATGGAGATATTTGCATAGAATGAAGGATCTCTGCAGATACATGAACAACTTGGAT
AGTGAAGCACATGCAAAATATAAAACCAAGTGGTGTATTCACCATGCTGGTCTTCTTTAAGAATGCA
TTCCAGTATGTCAACAGCATACAGCCATCTCTTCCAAGGTCTAATGCCCCGAGCCAAGTTCCTACTG
GTTCTTCTGAAGATGTATCGAATGTGTATGGTGTATAGAAATTGATCGTAATAAACACATCCATAAA
AAGAGGAACTAGCTGAAGGAAGAGAAAAACCATGAGTTCAGACGATGAAGACTGTTCCGCGAAAGGA
AGAAATCGTCACATTGTAGTCAATAAAGCCGAACCTGCTAACTCCACTGAAGTGTAGAAAGCTTTAAA
TTGGCCAGGGAGAGCTGGGAGTTGCTCTATTCCCTAGAATTCTTGACAAAGAATTTACAAGGATTTGC
TTGGCCTGGAAGACGGATACTTGCTTTGGTTAAGAATCTTCTCACTGATATGATCATCTATCAGGGT
CAATATAAAAAGGCGATAGCCAGCCTGCATCACTAGCAGCTCTCCAGGGATCCATTTCTCAGCCACAG
ATCACAGGGCAGGGGACCCTGGAGCATCAGAGGGCGCTCATCCAGCTGGCGACGTGCCACTTTGCGCTA
GGGAGTACAGAATGACATGTGAAAAAGTCTTGAATTTGATGTGCTACATGGTACTCCCCATTCAAGAT
GGAGGCAATCCAGGAGGAACCTCGAAAGTAAAGCCCAATTTAGAAAAGTTTCGGATCTGAAGCTC
CTGCTTGTACCAGCAAGGCTATCATGCCATACTGCCTCCATTTAATGTTAGCCTGTTTTAAGCTTAGA
GCTTTACAGACAACAGAGACGACATGGCATTGGGCGATGTGATTGTGTTGCTTCAGCAAGAGTGGCCA
CGGGGCGAGAATCTTTCTGAAAGCTGTCAATAAAATTTGCCAACAAGGAAATTTCCAATATGAGAAT
TTTTTCAATTACGTTACAAATATTGATATGCTGGAGGAATTTGCCCTACTTGAGAACTCAGGAAGTGGG
AAAATTCATCTGGAATTACTACCAATCAAGGAATGCTGATCAAGCACCACACTGTAACTCGAGGCATC
ACCAAAGGCGTGAAGGAGGACTTTCGCCTGGCCATGGAGCGCCAGGTCTCCCGCTGTGGAGAGAATCTG
ATGGTGGTTCTGCACAGGTTCTGCATTAATGAGAAGATCTTGCTCCTTCAGACTCTGACCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI
ACCN: NM_018142
Insert Size: 2133 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_018142.2</u>
RefSeq Size:	2539 bp
RefSeq ORF:	2133 bp
Locus ID:	55174
UniProt ID:	<u>Q9NVR2</u>
Cytogenetics:	8p21.3
MW:	82.2 kDa
Gene Summary:	INTS10 is a subunit of the Integrator complex, which associates with the C-terminal domain of RNA polymerase II large subunit (POLR2A; MIM 180660) and mediates 3-prime end processing of small nuclear RNAs U1 (RNU1; MIM 180680) and U2 (RNU2; MIM 180690) (Baillat et al., 2005 [PubMed 16239144]).[supplied by OMIM, Mar 2008]