

Product datasheet for **SC318195**

PHACTR4 (NM_023923) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PHACTR4 (NM_023923) Human Untagged Clone
Tag:	Tag Free
Symbol:	PHACTR4
Synonyms:	PPP1R124
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318195 representing NM_023923.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGACAAGCTGATGTCTCCAGACCGGTAAATCCAGATGCAGTTGAGGAAGCAGACCAGCCCACTACA
GAGCCAGGCATGGTCTTGACAGTGTGGAAGCAGGAGACACAACACCTCCTACCAAAAGGAAGAGCAAG
TTCTCAGGCTTTGGCAAGATCTTCAAGCCCTGGAAATGGAGGAAAAAAAAAAGTAGTGATAAATTTAAA
GAGACTTCAGAAGTTTGAACGAAAAATATCTATGCGAAAGCCAAGAGAAGAGCTGGTTAAAGAGGG
GTTCTGTTGGAAGACCCTGAGCAAGGTGGTGAGGATCCAGGAAAGCCAAGCGATGCCATGTTAAAGAAT
GGCCATACCAACCCCATAGGGAATGCCAGATCATCTAGTCCAGTCCAAGTAGAGGAAGAGCCAGTAAGA
TTAGCAAGTCTTAGGAAAGCTATTCCAGAAGAGGACCTAAAGAAACGACTAGGCTCAACTGGAAGCCAG
CCTAATTCTGAAGCAGAGTCTGTTCTGAGAATGTACCAAAACACCTTTACTTCTCCCAAAAGACCC
TTGTCCTCTTCTCATGAAGCAAGTGAAGGGCAAGCAAAGGATGCCACTTCTCTGGCGGCACGGCAAGG
TTCATCATCTCCACCTCCATCACCACAGCACCCGCTGCCACCCTGCTGCCACAAGCCTTGCAAAGACT
GTTAATCTCTCTGTCAACCCCTTCCCAAGCACCCAGGACTCTGCCTGCTGCTCCTGCCAGCACTAACACT
ACTGCTACCCCAAGCCTCACTCATATGGTCCCTGCCAAGCAGCCCTTATCCCTCCCTTAAACAGCT
CACAGAAATAGCAACCTGTGATTGCTGAACTGTCCCAAGCAATAACAGTGTTACATTGTTATCAAAA
CCGTCCCCACCCCTTACCACCTAAGAGAGGCATTCCATCAACCTCAGTACCCACCTTGAGTCTGCTGCT
GCCATCACCACAAAAACCAAGTGATGAAAGAGAGAAGAGCAGTGTCTATGGGCTCGGAAGTACTA
CCAATGATCTCACCTCGCTCTCCGTCCCCCCTGCTGCTACTCATATACCTCCAGAGCCTCCACGCACC
CCTCCATTCCCTGCTAAGACTTTTCAAGTTGTGCCAGAAATTGAGTTTCCACCATCCTTAGATCTACAC
CAGGAGATTCCCCAGCAGGAAGATCAGAAAAAGGAAGTCCCAAGAGGATACTGGACCAGAACTTTGGG
GAGCCCCATATACCCCTTAGGCTGCCTCCACTCCCACTGCATATTGGAATCCAGCAGGCCCTCACCAGC
CCACTTCCCATGACTCCTATTCTGGAGGGTCTCACAGAGCTCATTGTTGCTTTTTGAAAAAGTGAC
AGCTTTTCTGAGGACAGCAGTACGCTGGGTCGGACAGGTCTCTTCCCATCACTATTGAAATGCTAAAA
GTTCCAGACGATGAAGAAGAAGAGGAGCAACCTGTCCATCCACATTGAGTGAAGAAATGACACCTACC
TCAGTCATTCTTAAATTACCACAGTGTCTACGGGAGGAAGAAGAGAAGGAGAGCGACTCTGATTGAGAA
GGTCCCATTCAGTACCGAGATGAAGAAGATGAAGATGAAAGCTATCAGAGTGTCTCGCCAAACAAAGTG
AAGAGGAAAGACACACTGGCAATGAAGTTGAACCACAGACCCAGTGAACCAGAGTTGAACCTGAATTCT
TGGCCTTGTAAGCAAGGAGGAGTGAATGAAATACGGCACCAGATTGGAACACACTGATCCGCGCA
CTGAGTCAAAGACCAACACCAGAAGAACTAGAACAACGCAATATATTGCAACCTAAAAATGAAGCTGAT
CGTCAGGCAGAAAAACGAGAAATTAACGTCGGCTCACTAGAAAGCTCAGTCAAAGGCCAACTGTGCT
GAACCTCTTGCCAGGAAGATTCTGAGGTTAATGAATATGTAGAGGTAACAGATGCTCAAGATTATGAC
CGGCGAGCCGACAAACCTTGACCAAACTGACCCTGCTGACAAGGCTGCCATAAGAAAAGAAATTAAT
GAATTTAAAGCTCCGAGATGGAGGTTTCATGAAGAGAGCAACATTTTACACGCTACCATCGTCCATGA
```

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites: SgfI-MluI
ACCN: NM_023923
Insert Size: 2139 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_023923.3</u>
RefSeq Size:	6208 bp
RefSeq ORF:	2139 bp
Locus ID:	65979
UniProt ID:	<u>Q8IZ21</u>
Cytogenetics:	1p35.3
Domains:	RPEL
MW:	79.1 kDa
Gene Summary:	This gene encodes a member of the phosphatase and actin regulator (PHACTR) family. Other PHACTR family members have been shown to inhibit protein phosphatase 1 (PP1) activity, and the homolog of this gene in the mouse has been shown to interact with actin and PP1. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008] Transcript Variant: This variant (2) contains a distinct 5' UTR and initiates translation at an alternative AUG, compared to variant 1. The resulting isoform (2) has a distinct and longer N-terminus when compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.