

Product datasheet for **SC318197**

TNS4 (NM_032865) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TNS4 (NM_032865) Human Untagged Clone
Tag:	Tag Free
Symbol:	TNS4
Synonyms:	CTEN; PP14434
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC318197 representing NM_032865.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGTCCCAGGTGATGTCCAGCCCACTGCTGGCAGGAGGCCATGCTGTGAGCTTGGCCCTTGTGATGAG
CCCAGGAGGACCCTGCACCCAGCACCCAGCCCAAGCCTGCCACCCAGTGTCTTACTACACCACGGAA
GGCTGGGAGCCAGGCCCTGATGGCCCCGTGCCCTGCATGGGGCCCCCTGGCCGACTCCAGCAAGCC
CCACAGGTGGAGGCCAAAGCCACCTGCTTCTGCCGTCCCTGGTGAGAAGGCCCTTGGGGACCCAGAG
GACCTTGACTCCTACATTGACTTCTCACTGGAGAGCCTCAATCAGATGATCCTGGAAGTGGACCCACC
TTCCAGCTGCTTCCCCAGGGACTGGGGCTCCCAGGCTGAGCTGGCCAGAGCACCATGTCAATGAGA
AAGAAGGAGGAATCTGAAGCCTTGACATAAAGTACATCGAGGTGACCTCCGCCAGATCAAGGTGCCAC
GATGGCCCCCAGCACTGCTCCAGCCCTCTGTACCCCGCCCTTCGGCTCCCTTCGCAGTGGTGGCTC
CTCCTTTCCAGAGACGTCCCCGAGAGACGAAGCAGCAGTGAAGCCATCTTCTCTGGGAACCAG
GGCAGGGGGCACCAGCGCCCTCTGCCCCCTCAGAGGGTCTCTCCCTCGACCCCAAATCCCCCAGC
ATCTCAATCCCTTGCATGGGGAGCAAGGCCCTCGAGCCCCATGGTTTGGGCTCCCGCTGGTGGCTTCT
CCAAGACTGGAGAAGCGGCTGGGAGGCCTGGCCCCACAGCGGGGAGCAGGATCTCTGTGCTGTGAGCC
AGCCCACTGTCTGATGTGAGCTATATGTTTGAAGCAGCCAGTCCCTCCTGCACTCCAGCAACTCCAGC
CATCAGTCATCTTCCAGATCCTTGGAAAGTCCAGCCAACCTTCTCCTCCAGCCTCCACAGCCTTGGCTCA
GTGTCCCTGTGTACAAGACCCAGTGACTTCCAGGCTCCAGAAACCCACCCTAACCATGGGCCAACCC
AGAACACCCCACTCTCCACCACTGGCCAAAGAACATGCCAGCAGTGGCCCCATCCATCACCACCTCC
ATGGTGGACATACCCATTGTGCTGATCAACGGCTGCCAGAACCCAGGGTCTTCTCCACCCAGCGGACC
CCAGGACACCAGAACTCCGTTCAACCTGGAGCTGCTTCTCCAGCAACCCCTGTCCAGCCACCAGGAGC
AACAGCCAGACCCTGTCAGATGCCCTTTACCACATGCCCAGAGGGTCCCGCCAGGGACATGCAGCCC
ACCATGAAGTTCTGTGATGGACACATCTAAATACTGGTTTAAAGCCAAACATCACCCGAGAGCAAGCAATC
GAGCTGTGAGGAAGGAGGAGCCAGGGCTTTTGTATAAGGGACAGCTTTCATACCGAGGCTCCTTC
GGCCTGGCCCTGAAGGTGCAGGAGGTTCCCGCTCTGCTCAGAGTCGACCAGGTGAGGACAGCAATGAC
CTCATCCGACACTTCTCATCGAGTCGTCTGCCAAAGGAGTGCATCTCAAAGGAGCAGATGAGGAGCCC
TACTTTGGGAGCCTCTCTGCCTTCGTGTGCCAGCATTCCATCATGGCCCTGGCCCTGCCCTGCAAATC
ACCATCCCACAGAGAACTGGGAGGTGCAGATGGGGCTCGGACTCTACAGACAGCCAGCCTCCTGC
CAGAAGAAATCTGCGGGCTGCCACACCCTGTACCTGAGCTCAGTGAAGCGTGGAGACCCTGACTGGAGCC
CTGGCCGTGCAGAAAGCCATCTCCACCCTTTGAGAGGGACATCCTCCACAGCCACCCGTTGGTCCAC
TTCAAAGTCACAGAGCAGGGCATCACTCTGACTGATGTCCAGAGGAAGGTGTTTTCCGGCGCCATTAC
CCTCACCACCCTCCGCTTCTGTGGTATGGACCTGAGCAACGGAAGTGGCAGAAGTACTGCAAACCC
TCTGGATCTTTGGGTTTGTGGCCAAGAGCCAGACAGAGCCTCAGGAGAAGTATGCCACCTCTTTGCG
GAGTATGACATGGTCCAGCCAGCCTCGCAGGTATCGGCTGGTACTGCTCTGCTGCAGGACGCAGAA
AGGATGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_032865

Insert Size: 2148 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_032865.5](#)

RefSeq Size: 4090 bp

RefSeq ORF: 2148 bp

Locus ID: 84951

UniProt ID: [Q8IZW8](#)

Cytogenetics:	17q21.2
Domains:	SH2, PID
MW:	76.8 kDa
Gene Summary:	May be involved in cell migration, cartilage development and in linking signal transduction pathways to the cytoskeleton (By similarity). May promote apoptosis, via its cleavage by caspase-3.[UniProtKB/Swiss-Prot Function]