

Product datasheet for **SC318225**

RAI16 (FAM160B2) (NM_022749) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RAI16 (FAM160B2) (NM_022749) Human Untagged Clone
Tag:	Tag Free
Symbol:	RAI16
Synonyms:	RAI16
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318225 representing NM_022749.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCTGAGCCGGCTCGGGGCGCTGCTGCAGGAAGCCGTGGGGCGCGGAGCCAGCATTGACCTGCTG
CAGGCCCTTCGTGGAGCACTGGAAGGGCATCAGCACTACTACATCGAGAGCACAGATGAAAGCACCCCC
GCCAAGAAGACAGACATTCCCTGGCGGCTGAAGCAGATGCTGGATATCCTGGTGTATGAAGAGCAGCAG
CAGGCGGCCGCGGGTGAAGCAGGGCCCTGCCTGGAGTACCTGCTGCAGCAAGATCCTGGAGACTCTC
TGCACGCTGGCAAGGCCGAGTACCCCCAGGCATGCGGCAGCAGGTGTTCCAGTCTTCAGCAAGGTT
CTGGCGCAGGTGCAGCACCCCTGCTGCATTACCTCAGCGTCCACAGGCCTGTGCAGAACTCCTCCGA
CTTGGTGGGACTGCTTCCGGATCCGTTACAGAAAAGGAGGAGGTGCAGTTCACCACCGTCTCTGCTCC
AAGATCCAGCAGGACCCAGAGCTGCTCGCTACATCCTGGAAGGTAAAAAGATTGTAGGTAGGAAGAAA
GCATGCGGAGAACCCTGCCCCTGCCTAAGGACACAACCAGCCACGGGGACAAGGACTGCTCCCAGAT
GGTGTCTCTGCCAGGCCAGCTGGACGGGGAGTCTGTGGGGCCAGGCCTTGAACAGCCACATGCCT
GCTGAGACCGAGGAGCTGGACGGTGGGACCACAGAGAGCAACCTGATTACCTCCCTGCTTGGGCTGTGC
CAGAGCAAGAAGAGTCCGGTGGCCTTGAAGGCCAGGAGAACCTGCTGCTCCTGGTGAGCATGGCCTCC
CCAGCAGCTGCCACCTACCTGGTACAGAGCAGCGCCTGCTGCCCTGCGATCGTCCGGCACCTTTGCCAG
TTGTACCGGTCCATGCCTGTCTTCTGGACCCCGAGACATTGCCACCTTAGAGGGCATCAGCTGGAGG
TTACCCAGTGCCCGTCTGATGAGGCTTCTTCCCTGGCAAGGAGGCCTTGGCTGCCTTCTTGGGCTGG
TTTGATTACTGCGACCACCTCATCACAGAGGCACACACGGTGGTTGCGGACGCCTTGGCGAAGGCTGTG
GCTGAGAACTTCTTCGTGGAGACCCTGCAGCCCCAGCTCCTGCAGTGTCCGAGCAGAGCATCTTGACC
TCCACCGCCCTCCTCACAGCCATGCTGCGCCAGCTTCTGCTCCCTGCGCTGCTGCGGGAGGCCGTGGCT
TTCTCCTGGGCACAGACCGGCAGCCTGAAGCCCCCGGGGACAACCCCCACACCCTGTATGCTCATCTC
ATCGGGCATTGTGACCACCTCTCTGATGAGATCAGCATACCACACTCCGGCTGTTTGGAGGCTGCTG
CAGAAGCCCCACGAGGGGATCATCCACAGCCTGGTCTGCGCAACCTTGAAGGCCGCCCTTACGTGGCC
TGGGGCTCACCAGAGCCTGAGAGCTATGAGGACACCCTAGACCTGGAGGAAGACCCTACTTCACCAGC
AGCTTCTGGATTCCGGCTTTCAAACCTCCCGCAAAGCCTGCCTAGCTCCTGCTACCAGTTACGATGGC
AAAACAGCAGTGACCGAGATCGTCAACAGTTTCTGTGCCTGGTCCCGGAGGAAGCCAAGACCTCTGCC
TTCTGGAGGAGACAGGCTATGACACATACGTCCACGATGCTTATGGCCTGTTCCAGGAGTGCAGCTCC
CGCGTCGCTCCTGGGGCTGGCCTTGACCCCCACACCTTGGACCCCCATGAGCCCGAGCGACCTTTC
TTGAGGGCCACTTCTCCGAGTGTGTTGACCGCATGTCCCGATTCTGGATCAGCCATACAGCCTG
AACCTGCAGGTGACCTCGGTCTGTCCCGGCTTGCCCTTCCCCACCCCATATTCATGAGTACCTG
CTGGATCCGTACATCAGCCTGGCCCCGGCTGCAGGAGCCTATTCTCCGTGTTGGTGAGGGTGATCGGG
GACTTGATGCAGAGAATCCAGAGGGTACCCAGTCCCAGGCAAGCTGCTCCTGGTGCAGCAAGCAGTTG
ACGGGCCAGGCTCCTGGGGAGCAGCTGGACCACCAGACCTCCTCCAGGGCGTGGTGGTGTGGAGGAG
TTCTGCAAGGAGCTGGCTGCCATTGCCTTCGTCAAGTTTCCCCACATGATCCTCGCCAGAACGTCTCC
CCAGCCCCGAAGGGCAGGTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_022749

Insert Size: 2232 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_022749.5</u>
RefSeq Size:	3759 bp
RefSeq ORF:	2232 bp
Locus ID:	64760
UniProt ID:	<u>Q86V87</u>
Cytogenetics:	8p21.3
MW:	82.3 kDa
Gene Summary:	<p>Able to activate MAPK/ERK and TGFB signaling pathways (PubMed:22971576). May regulate the activity of genes involved in intestinal barrier function and immunoprotective inflammation (By similarity). May play a role in cell proliferation (PubMed:22971576). [UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (1) encodes the longest isoform (1). Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.</p>