

Product datasheet for **SC318226**

WDPCP (NM_015910) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	WDPCP (NM_015910) Human Untagged Clone
Tag:	Tag Free
Symbol:	WDPCP
Synonyms:	BBS15; C2orf86; CHDTHP; CPLANE5; FRITZ; FRTZ
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC318226 representing NM_015910.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGAGGCGAGAGTTTTGCTGGGACGCCTACTCCAAAGCGGCCGGGAGTCGCGCTTCTCCCACTCCCG
AGACAGGATAGAGATTCCTTCTGCCATCAGATGTCTTTCTGCTTGACTGAACTGCACCTGTGGTCTTTG
AAGAATACCTTACACATTGCGGATAGAGACATTGGGATCTACCAGTATTATGACAAGAAAGATCCACCA
GCGACAGAGCATGGTAACTTAGAAAAGAAGCAGAAGCTGGCAGAGTCACGAGATTATCCTTGACGCTC
AAAAACAGACGCCGAGAAAACTCCGAGACTCGCTCAAAGAGTTGGAGGAGCTGATGCAAAAACAGTCGG
TGTGTGCTGAGCAAAATGGAAGAACAATATGTCTGTGACGCTCCTTTTGGTTTCAGGTGTGCTGGTGTCT
CTAAGCCTTTCTGGGCCGACGCTGGAGAAAGTGGTGATTGACAGAAGCCTGGTGGGAAGCTCATCTCA
GACACCATCAGTGATGCTCTTCTCACAGACAGTTTTATCATCTTATCATTTTTGGCACAAAAACAATA
TGTTTTATTCAGTTTACCAAGAAGATGGAGTCTTCTGATGTAACAAAAGACTGGAAAACTCTCAGCC
TTGGATTATAAGATTTTCTATTATGAAATACCCGGCCCAATAACAAGACAACAGAGCGACATCTAGCT
ATCAACTGTGTTTCATGATAGAGTTGTTTCTGGTGGCCACTGGTCAACGATGATGCTTGGCCTTGGGCC
CCCATTTCTTCTGAGAAGGACAGAGCCAATCTACTCCTCCTGGGTTATGCTCAAGGAAGACTAGAGGTT
CTGAGTTCTGTCCGCACAGAATGGGACCCACTGGATGTTGCTTTGGCACAAACAGCCTTATCAGGTG
TTCACAGTGGAGCACTCCGTAAGTGTAGACAAAGACCCATGGCTGACAGCTGCATCTATGAATGCATT
CGGAATAAAATCCAGTGTGTGTCAGTACCAGAATACCACTAAAGTCAAAGGCCATCAGCTGCTGCAGG
AATGTTACTGAAGACAACTGATTTCTGGGCTGTGAAGATCTTCGCTAATCTTTATGAACTCACCGT
AGAGTACTCTCTTAGCACAGACTGAACTTTTGCCTTCATTAATAAGCTGCCACCAAGTGGTGCCTATT
CTGCTAGTTGGCAGCAACCAAGGGGAGTTGCAAATTTTTGATATGGCTCTATCCCTATTAACATCCAA
CTGTTGGCTGAAGACCGCTTACCAGGGAGACTCTGCAATTCAGTAAATATTTGATGCCTCCAGCAGT
CTTGTTCAAATGCAATGGATAGCTCCTCAGGTTGTTTCTCAGAAGGGTGAAGGTAGTGATATCTATGAT
CTCCTCTTCTCAGGTTTGAAGAGGACCTTTGGGTGTGCTGTTGTTAAACTAGGCGTCTTCACTCGA
GGACAGCTGGGCTGATAGACATCATCTCCAGTACATTCAGTGTGATGAGATCTATGAGGCAATAAAC
ATCCTGAGCAGCATGAACTGGGACACTCTGGGCCACCAAGTGTCTTATCAGCATGAGCGCCATTGTAAC
CATCTTCTTAGACAGAAGCTCACTCCAGAGAGAGAAGCACAGCTTGAGACAAGCCTTGGAACCTTCTAT
GCTCCAACAAGACCACTTCTGGATTCCACTATATTGGAATATAGAGATCAAATCAGCAAATATGCAAGG
AGATTCTCCATCACTTGTCTCAGGTACCAGAGGTTTGAAGGCAATTTCTCCTAGCTGTTGACGTTGGT
GCTCGTGACCTCTTTATGGATATTCATTACCTTGCACTAGATAAAGGTGAATTGGCACTAGCTGAAGTG
GCAAGAAAAAGAGCTAGTGACATTGATGCAAGATCAATAACCTCTGGGTTGAACTCCTGGGACCTTG
GACAGAGGGGATATGCTAAATGAAGCATTTATTGGCCTGTCTTTAGCACCTCAAGGAGAAGACTCATTT
CCAGATAACCTCCCTCCCTCTTGCCCAACCCACAGACATTTTTACAACAAGAACTACTGAATGGCTCT
TCTAACAGACAAATAATTGACAGAAGGAATGAACCTGAAAAAGACATCTGTTCTGGATTTTTGATGACT
AATACCTGTAATGCAGAAGACGGAGAAGTGAAGAGACGGCAGAGAAGCAAGGAAATCAGAGATGGTGGT
TCTCTCAAAATGATTCACCTTTGGTCTGGTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_015910

Insert Size: 2241 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_015910.5](#)

RefSeq Size: 3326 bp

RefSeq ORF: 2241 bp

Locus ID: 51057

UniProt ID: [O95876](#)

Cytogenetics: 2p15

MW: 85.1 kDa

Gene Summary: This gene encodes a cytoplasmic WD40 repeat protein. A similar gene in frogs encodes a planar cell polarity protein that plays a critical role in collective cell movement and ciliogenesis by mediating septin localization. Mutations in this gene are associated with Bardet-Biedl syndrome 15 and may also play a role in Meckel-Gruber syndrome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2014]
Transcript Variant: This variant (2) represents the longest transcript and encodes the longer isoform (2).