

## Product datasheet for **SC318228**

### **NOX5 (NM\_024505) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	NOX5 (NM_024505) Human Untagged Clone
Tag:	Tag Free
Symbol:	NOX5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC318228 representing NM\_024505.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGAACACATCTGGAGACCCAGCCAGACGGCCCTGAAGGCTGTAGAGGCACCATGAGTGCCGAGGAG
GATGCCAGGTGGCTCCGGTGGGTGACTCAGCAGTTTAAGACCATTGCAGGAGAAGATGGGAGATCAGC
CTGCAAGAATTCAAAGCAGCTCTGCATGTGAAAGAGTCCTTCTTTGCAGAGCGATTCTTTGCCCTATTT
GACTCCGATAGAAGTGGCACCATCACCTCCAGGAGCTGCAGGAGGCACTGACCTGCTCATCCATGGC
AGCCCCATGGACAACTCAAATTCCTCTCCAGGTGTATGACATCGATGTGTGTGCACGGCAGGGGGCG
TCTGCAGGTACAGAGTGGGTGCTGGGGCAGGCCCGCACTGGGCTTCATCCCACCTCGGGACAGGCAGT
GGCTCCATTGACCCGGATGAGCTGCGCACTGTGCTGCAGTCGTGTCTGCGCAGAGCGCCATCTCGCTG
CCTGACGAGAAGCTGGACCAGCTGACGCTGGCGCTCTTGAATCGGCCGACGCGGACGGCAACGGGGCC
ATCACCTTCGAGGAGCTCCGGGACGAGCTGCAGCGCTTCCCGGAGTCATGGAGAACCTGACCATCAGC
GCTGCCACTGGCTGACGGCCCCGCCCGCCACGCCGCGCCGGCCGCGCCAGCTGACCCGCGCC
TACTGGCACAACCACCGAGCCAGCTGTCTGCCTGGCCACCTATGCAGGCCTCCACGTGCTGCTTTC
GGGCTGGCGGCCAGCGCCACCGGGACCTCGGGCCAGCGTCATGGTGGCCAAGGGCTGGCGCCAGTGC
CTCAACTTCGACTGCAGCTTCATCGCGGTGCTGATGCTCAGACGCTGCCTCACCTGGCTGCGGGCCACG
TGGCTGGCTCAAGTCTACCACTGGACCAGAACATCCAGTTCCACCAGCTTATGGGCTACGTGGTAGTG
GGGCTGTCCCTCGTGACACCCGTGGCTCACACTGTGAATTTGTACTCCAGGCTCAGGCGGAGGCCAGC
CCTTCCAGTTCTGGGAGCTGCTGCTCACCACGAGGCCGGCATTGGCTGGGTACACGGTTCCGGCTCC
CCGACAGGTGTCGCTCTGCTGCTGCTCCTCCTCATGTTTCTGCTCCAGTTCTGCTCCAGTCCGAGG
AGTGGCCACTTTGAGGTGTTCTATTGGACTCACCTGCTCCTACCTCCTCGTGTGGCTTCTGCTCATCTTT
CATGGGCCCAACTTCTGGAAGTGGCTGCTGGTGCCTGGAATCTTGTTTTTCTGGAGAAGGCCATCGGA
CTGGCAGTGTCCCGCATGGCAGCCGTGTGCATCATGGAAGTCAACCTCCTCCCCTCCAAGTCACTCAT
CTCCTCATCAAGCGGCCCTTTTTTTTCACTATAGACCTGGTACTACTTGTATCTGAACATCCCCACC
ATTGCTCGCTATGAGTGGCACCCCTTACCATCAGCAGTGTCTCCTGAGCAGAAAGACACTATCTGGCTG
CACATTCGGTCCCAAGGCCAGTGGACAAACAGGCTGTATGAGTCCTTCAAGGCATCAGACCCACTGGGC
CGTGGTCTAAGAGGCTGTGAGGAGTGTGACAATGAGAAAGAGTCAAAGGTCGTCGAAGGGCTCTGAG
ATACTTTTGGAGAAACAAAATTCTGTAACATCAAGTGTACATCGATGGGCCTTATGGACCCCCACC
CGCAGGATCTTTGCCTCTGAGCATGCCGTGCTCATCGGGCAGGCATCGGCATCACCCCTTTGCTTCC
ATTCTGCAGAGTATCATGTACAGGCACCAGAAAAGAAAGCATACTTGCCCGAGCTGCCAGCACTCCTGG
ATCGAAGGTGTCCAAGACAACATGAAGTCCATAAGGTGGACTTTATCTGGATCAACAGAGACCAGCGG
TCTTTCGAGTGGTTTGTGAGCCTGCTGACTAACTGGAGATGGACCAGGCCGAGGAGGCTCAATACGGC
CGTTCCTGGAGCTGCATATGTACATGACATCTGCACTGGGAAGAATGACATGAAGGCCATTGGCTG
CAGATGGCCCTTGACCTCCTGGCCAACAAGGAGAAGAAAGACTCCATCACGGGGCTGCAGACGGCACC
CAGCCTGGGCGGCTGACTGGAGCAAGGTGTTCCAGAAAGTGGCTGCTGAGAAAGGGCAAGGTGCAG
GTCTTCTGTGGCTCCCAGCTCTGGCAAGGTGCTGAAGGGCCATTGTGAGAAGTTCGGCTTCAGA
TTTTTCCAAGAGAATTTCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_024505

**Insert Size:** 2298 bp

**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:	<u>NM_024505.3</u>
RefSeq Size:	2809 bp
RefSeq ORF:	2298 bp
Locus ID:	79400
UniProt ID:	<u>Q96PH1</u>
Cytogenetics:	15q23
Protein Families:	Ion Channels: Other, Transmembrane
MW:	86.4 kDa
Gene Summary:	<p>This gene is predominantly expressed in the testis and lymphocyte-rich areas of spleen and lymph nodes. It encodes a calcium-dependen NADPH oxidase that generates superoxide, and functions as a calcium-dependent proton channel that may regulate redox-dependent processes in lymphocytes and spermatozoa. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Oct 2011]</p> <p>Transcript Variant: This variant (1, alternately referred to as NOX5 gamma and v3) encodes the longest isoform (1).</p>