

Product datasheet for **SC318234**

SLC39A6 (NM_012319) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SLC39A6 (NM_012319) Human Untagged Clone
Tag:	Tag Free
Symbol:	SLC39A6
Synonyms:	LIV-1; ZIP6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC318234 representing NM_012319.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGCGAGGAAGTTATCTGTAATCTTGATCTGACCTTGGCCTCTCTGTACAAATCCCCTTCATGAA
CTAAAAGCAGCTGCTTTCCCCAGACCCTGAGAAAATTAGTCCGAATTGGGAATCTGGCATTAAATGTT
GACTTGGCAATTTCCACACGGCAATATCATCTACAACAGCTTTTCTACCGCTATGGAGAAAATAATTCT
TTGTCAAGTTGAAGGGTTCAGAAAATTACTTCAAAATATAGGCATAGATAAGATTAAGAAATCCATATA
CACCATGACCACGACCATCACTCAGACCACGAGCATCACTCAGACCATGAGCGTCACTCAGACCATGAG
CATCACTCAGACCACGAGCATCACTCTGACCATGATCATCACTCTCACCATAATCATGCTGCTTCTGGT
AAAAATAAGCGAAAAGCTCTTGGCCAGACCATGACTCAGATAGTTCAAGTAAAGATCCTAGAAACAGC
CAGGGGAAAGGAGCTCACCGACCAGAACATGCCAGTGGTAGAAGGAATGTCAAGGACAGTGTAGTGT
AGTGAAGTGACCTCAACTGTGTACAACACTGTCTCTGAAGGAACACTTTCTAGAGACAATAGAGACT
CCAAGACCTGAAAACCTTTCCCAAAGATGTAAGCAGCTCCACTCCACCCAGTGTACATCAAAGAGC
CGGGTGAGCCGGCTGGCTGGTAGGAAAACAAATGAATCTGTGAGTGAGCCCCGAAAAGGCTTTATGTAT
TCCAGAAACACAAATGAAAATCCTCAGGAGTGTTCATGCATCAAAGCTACTGACATCTCATGGCATG
GGCATCCAGGTTCCGCTGAATGCAACAGAGTTCACACTATCTGTCCAGCCATCATCAACCAAATGAT
GCTAGATCTTGTCTGATTCATACAAGTAAAAGAGGCTGAAATCCCTCCAAAGACCTATTCATTACAA
ATAGCCTGGGTTGGTGGTTTTATAGCCATTTCCATCATCAGTTTCTGTCTCTGTGGGGTTATCTTA
GTGCCTCATGAATCGGGTGTTTTTCAAATTTCTCCTGAGTTTCTGTGGCACTGGCCGTTGGGACT
TTGAGTGGTGTGCTTTTTACACCTTCTCCACATTCTCATGCAAGTCAACCCATAGTCATAGCCAT
GAAGAACCAGCAATGAAAATGAAAAGAGGACCCTTTTCACTCATCTGTCTTCTCAAACATAGAAGAA
AGTGCCTATTTTGAATCCACGTGGAAGGTCTAACAGCTCTAGGAGGCCTGTATTTTCATGTTTCTGTT
GAACATGTCTCATTGATCAAACAATTTAAAGATAAGAAGAAAAGAATCAGAAGAAACCTGAAAAT
GATGATGATGTGGAGATTAAGAAGCAGTTGTCCAAGTATGAATCTCAACTTTCAACAAATGAGGAGAAA
GTAGATACAGATGATCGAACTGAAGGCTATTTACGAGCAGACTCACAAAGGCCCTCCCACTTTGATTCT
CAGCAGCCTGCAGTCTTGAAGAAGAAGAGGTGATGATAGCTCATGCTCATCCACAGGAAGTCTACAAT
GAATATGTACCCAGAGGGTGCAAGAATAAATGCCATTACATTTCCACGATACACTCGGCCAGTCAGAC
GATCTCATTACCACCATCATGACTACCATCATATTCTCCATCATCACCACCACCAAAACCACCATCT
CACAGTCACAGCCAGCGTACTCTCGGGAGGAGCTGAAAGATGCCGGCGTCGCCACTCTGGCCTGGATG
GTGATAATGGGTGATGGCTGCACAATTTCAAGCATGGCCTAGCAATTGGTGTGCTTTTACTGAAGGC
TTATCAAGTGGTTAAGTACTTCTGTTGCTGTGTTCTGTGATGAGTTGCCTCATGAATTAGTGACTTT
GCTGTTCTACTAAAGGCTGGCATGACCGTTAAGCAGGCTGTCTTTATAATGCATTGTGACCCATGCTG
GCGTATCTTGAATGGCAACAGGAATTTTCATTGGTCATTATGCTGAAAATGTTTCTATGTGGATATTT
GCATTAAGTGTGGCTTATTCATGTATGTTGCTCTGGTTGATATGGTACCTGAAATGCTGCACAATGAT
GCTAGTGACCATGGATGTAGCCGCTGGGGTATTTCTTTTACAGAATGCTGGGATGCTTTTGGGTTTT
GGAATTAAGTACTTATTTCCATATTTGAACATAAAAATCGTGTTCGTATAAATTTCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:



ACCN: NM_012319

Insert Size: 2268 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:	NM_012319.3
RefSeq Size:	3637 bp
RefSeq ORF:	2268 bp
Locus ID:	25800
UniProt ID:	Q13433
Cytogenetics:	18q12.2
Domains:	Zip
Protein Families:	Druggable Genome, Transmembrane
MW:	85 kDa
Gene Summary:	<p>Zinc is an essential cofactor for hundreds of enzymes. It is involved in protein, nucleic acid, carbohydrate, and lipid metabolism, as well as in the control of gene transcription, growth, development, and differentiation. SLC39A6 belongs to a subfamily of proteins that show structural characteristics of zinc transporters (Taylor and Nicholson, 2003 [PubMed 12659941]).[supplied by OMIM, Mar 2008]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>