

## Product datasheet for **SC318237**

### ASCC2 (NM\_032204) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ASCC2 (NM_032204) Human Untagged Clone
Tag:	Tag Free
Symbol:	ASCC2
Synonyms:	ASC1p100; p100
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC318237 representing NM\_032204.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGCCAGCTTGCCCTGGACCAACTCCAGATCACCCACAAGGACCCGAAGACAGGAAAGCTGAGGACT
TCACCAGCGCTGCACCCGAGCAGAAGGCAGACCGGTATTTTGTGTTATACAAACCGCCCTAAAGAC
AACATTCGCGCCCTAGTGGAGGAGTACCTGGAACGCGCCACCTTCGTAGCCAATGACCTCGACTGGCTC
CTGGCCTTGCCCTCACGATAAATTCTGGTGCCAGGTGATCTTTGACGAGACTCTACAGAAGTGCCTGGAC
TCCTACCTGCGCTATGTCCCGCAAATTCGACGAGGGGGTGGCCTCAGCCCTGAGGTTGTTGACATG
CAGAAGCGCCTCCATCGAAGTGTCTTCTCACCTTCTCCGCATGTCCACTCACAAGGAATCCAAAGAT
CACTTCATTTCCCTTCTGCGTTTGGAGAAATCCTCTACAATAACTTCTCTTTGACATTCCAAAGATC
CTGGACCTCTGCGTGCTCTTTGGAAAAGGCAACTCACCCTGCTCCAGAAGATGATAGGAAACATCTTT
ACACAGCAGCCAAGTTACTACAGTGACCTGGATGAAACCCTGCCTACCATCCTTCAGGTCTTCAGCAAT
ATCCTCCAGCACTGTGGTTTGAAGGGGACGGGGCAATACCACACCCAGAAGCTTGAGGAGAGGGGC
CGATTGACCCCAAGTACATGCCTCTCCTGGAATTAAGGACATTGTTCTCTACCTTTGTGATACCTGC
ACCACACTTTGGGCTTTCTGGATATCTCCCTTTGGCTTGCCAGACCTCCAGAAGCAGCACTTTTGT
TACAGACTAGCTTCTTCTACGAAGCAGCAATCCCGAAATGGAGTCTGCAATTAAGAAGAGGAGGCTT
GAAGATAGCAAGCTTCTTGGTGACCTGTGGCAGAGGCTCTCCATTCCAGGAAGAAGCTAATGGAGATT
TTCCACATCATCTGAACCAGATCTGCCTCCTTCCCATCTAGAAAGCAGCTGTGACAACATTAGGGC
TTCATCGAAGAGTTCTTTCAGATCTTTCAGCTCCTTGTGCAGGAGAAGAGGTTCTCCGGGACTATGAT
GCACTCTTCCCGTGGCCGAAGACATCAGCTTGTGCAGCAGGCTCATCAGTCTTGACGAGACGCGG
ACTGCCTACATCCTCCAGGCAGTCGAGAGTGCATGGGAAGGGTGGACAGACGGAAGCCACAGATGCT
AAAGACCCATCGGTGATTGAGGAGCCTAATGGGGAGCCTAACGGGGTACGGGTGACAGCAGAGGCAGTC
AGTCAAGCATCATCACATCCGGAGAACTCGGAGGAAGAGGAGTGCATGGGAGCAGCCGCGCTGTGGGC
CCTGCCATGTGTGGGGTGGAACTGGACTCTCTCATCTCCAAGTGAAGGACCTGCTGCCAGACCTTGGT
GAGGGCTTTCATCTGGCCTGCCTGGAGTACTACCACTACGACCCAGAGCAGGTGATCAACAATATCTG
GAGGAGCGGCTGGCCCCACCCTCAGCCAGCTGGACCGCAACCTAGACAGAGAAATGAAACCAGACCCT
ACACCCCTGCTGACGTCTCGCCACAACGTCTTCCAGAATGACGAGTTTGTGTTTCAGCAGGGACTCA
GTAGACCTGAGCCGGTGCACAAGGGCAAGAGCACCAGGAAGGAGGAAAACACCGGGAGTTTGTGAAC
GACAAGCGTGCAAGTGGCGGCACAGCGGCAGCGCTACGAGCAGTACAGCGTGGTGGTGGAGGAGTGCCA
CTGCAGCCAGGGCAGAGCCTGCCCTACCACAGTGTCTACTACGAGGATGAGTACGATGACACATACGAT
GGCAACCAGGTGGGCGCAATGATGCAGACTCTGATGACGAGCTCATCAGCCGAGGCCATTACCATC
CCTCAGGTGCTGAGAACCAAGTGCCTAGAGAAGGGCAGGAGGAGGATGACGACGATGAGGAAGACGAT
GCTGACGAGGAGGCTCCCAAGCCGACCATTTTGTTCAGGACCCTGCAGTGTGAGAGAGAAGGCAGAA
GCCAGGCGCATGGCCTTTCTCGCCAAGAAAGGGTACCGGCATGACAGCTCAACAGCAGTGGCCGGCAGC
CCCCGAGGCCATGGGCAGAGCCGCGAGACAACCCAGGAACGCAGGAAGAAGGAAGCAACAAGGCGACA
AGAGCCAACCACAACCGGAGAACCATGGCCGACCGCAAGAGGAGCAAAGGCATGATCCCATCTGTA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**


**ACCN:** NM\_032204

**Insert Size:** 2274 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_032204.4](#)

**RefSeq Size:** 2840 bp

**RefSeq ORF:** 2274 bp

**Locus ID:** 84164

**UniProt ID:** [Q9H118](#)

<b>Cytogenetics:</b>	22q12.2
<b>Domains:</b>	CUE
<b>Protein Families:</b>	Druggable Genome, Transcription Factors
<b>MW:</b>	86.4 kDa
<b>Gene Summary:</b>	Plays a role in DNA damage repair as component of the ASCC complex. Recruits ASCC3 and ALKBH3 to sites of DNA damage by binding to polyubiquitinated proteins that have 'Lys-63'-linked polyubiquitin chains (PubMed:29144457). Part of the ASC-1 complex that enhances NF-kappa-B, SRF and AP1 transactivation (PubMed:12077347).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).