

Product datasheet for **SC318238**

CEP68 (NM_015147) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CEP68 (NM_015147) Human Untagged Clone
Tag:	Tag Free
Symbol:	CEP68
Synonyms:	KIAA0582
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC318238 representing NM_015147.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGCCCTGGGTGAAGAAAAGGCAGAAGCGGAAGCATCTGAAGACACAAAGGCCAGTCCTATGGGAGA
GGGAGCTGCAGGGAGCGGGAGCTGGACATCCCAGGGCCATGAGTGGGAGCAGCCCCACGCCTGGAA
GCTGAGGGAGGGCTCATCTCCCCTGTATGGGGGGCAGAAGGGATACCTGCCCTACTTGCTGGATTGGG
ACTGACCCCTGGCGGCCCTCTAGAGCCACCAGCCACAGGCCAGTGATGCCAACAGAGAGCCGTAGCT
GAGAGGTCTGAGCCTGCACTCAGTGGCCTGCCTCCTGCCACCATGGGGTCTGGAGACCTTCTGCTCTCC
GGGAAAGCCAGGTGGAGAAGACCAAGCTTTCTCCTCCGAGGAGTTCCTCAGACTCTGAGCCTTCCC
AGAACAACAATAATTTGCTCAGGACATGATGCTGATACCGAAGATGATCCATCCCTAGCAGATTTGCC
CAGGCACTGGACCTCAGCCAGCAGCCTCACAGCTCAGGTCTCTTGCCTGTACAGTGGAAAGTCCGTG
CTGAGCCAGGTTCCGCAGCTCAGCCTCCAGCTGCAGCATCTCTGCTTCTCCACAGGCAGCAGTCTC
CAGGGTACCAGGAGAGGGCGGAGCCTCGTGGTGGTTCTCTGGCCAAGGTCTCCTCCTCCCTGGAGCCG
GTCGTCCCCCAGGAACCTTCTCTGTGGTGGGGCTAGGACCTCGGCCCCAGTGGTACCACAGCCTGTG
TTCTCTGGGGTATGCTTCTGGGCTAGGCAGGAGACGCCTCTCCTTCCAGGCTGAGTACTGGGCCTGT
GTGCTGCCAGATTCCTGCCTCCATCACCCGACCGCCACTCCCCTCTCTGGAACCCAAATAAAGAGTAT
GAAGATCTGCTTGACTATACTTACCCACTGAGGCCCGGGCCTCAGCTCCCAAAGCACCTTGATAGCCGT
GTGCCAGCTGACCCTGTCTGCAGGACTCCGGGTAGACTGGATAGCTTCTCTGTCTCTCCAGCAAGC
ACCCTCAAATCACCTACTAATGTCTCCCCAACTGCCACCAGCAGAGGCCACTGCCCTGCCATTTTCT
GGGCCCAGAGAGCCAAGCCTTAAGCAGTGGCCCTCCAGAGTACCCAGAAACAGGGTGGCATGGGCTTG
GCATCTTGGAGCCAATTCATCTACCCCCAGAGCCCCAGCAGTAGGGATGCTCGTTGGGAGCGCAGA
GAGCCAGCCCTGAGGGGTGCGAAGGACCGCTGACTATAGGCAAGCACCTTGATATGGGCTCTCCCAG
CTAAGGACACGGGACAGAGGGTGGCCCTCGCCCAGGCCAGAGAGGGAGAAGAGGACCAGCCAGAGTGCC
CGGCGCCCTACCTGCACAGAGTCTAGGTGGAAATCAGAAGAGGAAGTGGAAAGTATGACAGATATCTT
GCCCTCCCCTCGGCTGACACAGGTTTCTAGCCTGGTTTCGATCTAGGATCCATTTCTACCTGGTT
ACCCTGCCACTGGGGATATCAAAGGGCAGAGCCCCTTGGAAAGTGTGACAGAGTATGGGCCAGCTTCC
TTCCCTTCAAGCTCCAGCCAAAGCCAGCTTCCCCTGGAGCTGCCCTCCAAGGATCTGGGGATCCTGAG
GGCCAGAAATCCCTGTTTCTGCGCTCCTTCGTCCGTGCCACGACTCCGCAGGGGAAGGCAGTCTGGGG
AGCAGCCAGGCCCTCGGGTCTCCTCTGGACTGCTGAAAACACGCCCTCCTTCCAGCTAGGTTGGAC
CGGTGGCCATTCTCAGACCCAGATGTTGAAGGGCAGCTTCCCAGGAAAGGAGGAGAACAGGGAAAAGAA
TCACTGGTGAATGTGTGAAGACATTTTGTGTGACAGTGGAAAGAGCTGATCTGCTGGCTGTATAATGTT
GCAGATGTTACTGACCACGGGACTGCAGCCAGGTCCAATCTTACAAGTCTCAAGTCTTCTCTGCAGCTT
TACCGCAATTTAAGAAAGATATAGATGAACATCAGTCTCTGACGGAGAGTGTCTTACAGAAGGGGGAG
ATTCTTCTCAGTGCCTGTTGGAGAACACCCAGTTTTAGAGGATGTCTTGGGAGGATCGCAAAGCAG
TCTGGTGAAGTGGAGAGCCACGCAGATCGCCTGTATGACTCTATCTTGGCCTCTCTGGACATGCTGGCT
GGCTGCACCCTTATCCCTGACAAAAAGCCATGGCGGCAATGGAGACCCATGTGAAGGGGTTTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:


ACCN: NM_015147

Insert Size: 2274 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_015147.2](#)

RefSeq Size: 5877 bp

RefSeq ORF: 2274 bp

Locus ID: 23177

UniProt ID: [Q76N32](#)

Cytogenetics: 2p14

MW: 81.1 kDa

Gene Summary: Involved in maintenance of centrosome cohesion, probably as part of a linker structure which prevents centrosome splitting (PubMed:18042621). Required for localization of CDK5RAP2 to the centrosome during interphase (PubMed:24554434, PubMed:25503564).[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (1) represents the longest transcript and encodes the longer isoform (a). Variants 1 and 2 both encode the same isoform (a).