

## Product datasheet for **SC318239**

### **KIRREL1 (NM\_018240) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	KIRREL1 (NM_018240) Human Untagged Clone
Tag:	Tag Free
Symbol:	KIRREL1
Synonyms:	KIRREL; NEPH1
Mammalian Cell Selection:	None
Vector:	<u><a href="#">pCMV6-XL4</a></u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >OriGene sequence for NM\_018240 edited

```

CCGCTCCTCTGGCTCGCTCGCAACTCCGGGCCCCAGCCGCGGGCGGGCGCACGGCGGGC
GGACAGCATGCTGAGCCTCCTCGTCTGGATCCTCACTCTCTCCGATACTTTCTCCCAAGG
GACCCAGACCCGCTTCAGCCAGGAGCCAGCTGACCAGACGGTGGTGGCTGGACAGCGGGC
CGTGCTCCCCTGTGTGCTGCTCAACTACTCTGGAATTGTGCAATGGACCAAGGACGGCT
GGCCCTGGGCATGGGCCAGGGCCTCAAAGCCTGGCCACGGTACCGGGTTGTGGGCTCCGC
AGACCTGGGCAGTACAACCTGGAGATCACAGATGCTGAGCTCTCTGACGACGCCTCTTA
CGAGTGCCAGGCCACGGAGGCCCTGCGCTCTCGGCGGGCCAAACTACCGTGTCTCAT
CCCCCAGAGGACACCAGGATTGACGGAGGCCCTGTGATTCTACTGCAGGCAGGCACCCC
CCACAACCTCACATGCCGGGCTTCAATGCGAAGCCTGCTGCCACCATCATCTGGTTCCG
GGACGGGACGCAGCAGGAGGGCGCTGTGGCCAGCACGGAATTGCTGAAGGATGGGAAGAG
GGAGACCACCGTGAAGCAACTGCTTATTAACCCACGGACCTGGACATAGGGCGTGTCTT
CACTTGCCGAAGCATGAACGAAGCCATCCCTAGTGGCAAGGAGACTTCCATCGAGCTGGA
TGTGCACCACCCTCTACAGTGACCCTGTCCATTGAGCCACAGACGGTGCAGGAGGGTGA
GGGTGTTGTCTTTACCTGCCAGGCCACAGCCAACCCCGAGATCTTGGGCTACAGTGGGC
CAAAGGGGGTTTTCTTGATTGAAGACGCCACGAGAGTCGTATGAGACAATGTGGATTA
TTCTTTTTTACGGAGCCTGTGTCTTGTGAGGTTCAACAAGTGGGAAGCACCATGT
CAGCACTTTAGTAAATGTCCACTTTGCTCCCCGATTGTAGTTGACCCCAAAACCCACAAC
CACAGACATTGGCTCTGATGTGACCCCTTACCTGTGTCTGGGTTGGGAATCCCCCCTCAC
TCTCACCTGGACAAAAGGACTCAAATATGGTCTGAGTAACAGCAACCAGCTGCTGCT
GAAGTCGGTGACTCAGGCAGACGCTGGCACCTACACCTGCCGGGCCATCGTGCCTCGAAT
CGGAGTGGCTGAGCGGGAGGTGCCGCTCTATGTGAACGGGCCCCCATCATCTCCAGTGA
GGCAGTGCAGTATGCTGTGAGGGGTGACGGTGGCAAGTGGAGTGTTCATTGGGAGCAC
ACCAACCCAGACCCGATAGCATGGGCCTGGAAGGAGAAGTCTTGGAGGTGGGGACCT
GGAACGCTATACAGTGGAGAGGACCAACTCAGGCAGTGGGGTGTATCCACGCTACCAT
CAACAATGTCATGGAGGCCGACTTTCAGACTCACTACAACCTGCACCGCCTGGAACAGCTT
CGGGCCAGGCACAGCCATCATCCAGCTGGAAGAGCGAGAGGTGTACCTGTGGGCATCAT
AGCTGGGGCCACCATCGGCGGAGCATCCTGCTCATCTTCTTCTCATCGCCTTGGTATT
CTTCTCTACCGGCGCCGAAAGGCAGTCGAAAGACGTGACCCTGAGGAAGCTGGATAT
CAAGGTGGAGACAGTGAACCGAGAGCCACTTACGATGCATTCTGACCGGGAGGATGACAC
CGCCAGCGTCTCCACAGCAACCCGGGTCATGAAGGCCATCTACTCGTCGTTTAAGGATGA
TGTGGATCTGAAGCAGGACCTGCGCTGCGACACCATCGACACCCGGGAGGAGTATGAGAT
GAAGGACCCCAACATGGCTACTACAACGTGCGTGCCCATGAAGACCGCCGCTTCCAG
GGCAGTGTCTATGCTGACTACCGTGGCCCTGGCCCTGCCCGCTTCGACGGCCGCCCTC
ATCCCGTCTCTCCCACTCCAGCGGCTATGCCAGCTCAACACCTATAGCCGGGGCCCTGC
CTCTGACTATGGCCCTGAGCCCACACCCCTGGCCCTGCTGCCCCAGCTGGCACTGACAC
AACCAGCCAGCTGTCTACGAGAAGTATGAGAAGTTCACTCCCATCCCTTCCCTGGGGC
AGCTGGGTACCCACCTACCGACTGGGCTACCCCAAGGCCCCACCTCTGGCCTGGAGCG
GACCCCATATGAGGCGTATGACCCATTGGCAAGTACGCCACAGCCACTCGATTCTCCTA
CACCTCCAGCACTCGGACTACGGCCAGCGATTCCAGCAGCGCATGCAGACTCACGTGTA
GGGGCCAGAGCCTGGCTGGGGCATCTCTGCGGGGAGAGGAGAAGGCTTTCACAGCTGTT
CCCTGATATTAGGGGATTGCTCATTGCTCCCTTCTCGGACCAGCCTTCTTCTCCAC
CATGGCAGGTGGGAGCAGGCTCTCCAGAAACACCCCGTCCCGAGGATGGTGTCTGTGC
ATGCCCCAGCCTCTGGGCTGCCCTTCCCTTCTTCTCGGGAGGATGTGTCTTCTGAC
CTGCACTTGTGCTGACCCTAGAATGGGGACAGGAAAGTGAAGGTTAGGGAAAGCAGAG
GGGGGCACTTTTTAGCATTCCCTTCTATCCCACCCCTCTGATCTCCATAAGTGGAAAT
GGGGGTACCCAGGATGGGCAGGCTTTGGCCTAGGGACATGAAGTATGGGAGTGGGTGGC
TGTGGCACAGACAGGTGGAACGGGATAGCCTGGCCAGTCCCTCTGTTGTCTGCATTCC
TGCCCTGGGTGCCTCTCTCTCCTCAGGGTACTGCAGAAGGGAGCGAACAGGG

```

**Restriction Sites:** Please inquire

<b>ACCN:</b>	NM_018240
<b>Insert Size:</b>	3000 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_018240.5</a> , <a href="#">NP_060710.3</a>
<b>RefSeq Size:</b>	3628 bp
<b>RefSeq ORF:</b>	2274 bp
<b>Locus ID:</b>	55243
<b>UniProt ID:</b>	<a href="#">Q96J84</a>
<b>Cytogenetics:</b>	1q23.1
<b>Protein Families:</b>	Transmembrane
<b>Gene Summary:</b>	<p>NEPH1 is a member of the nephrin-like protein family, which includes NEPH2 (MIM 607761) and NEPH3 (MIM 607762). The cytoplasmic domains of these proteins interact with the C terminus of podocin (NPHS2; MIM 604766), and the genes are expressed in kidney podocytes, cells involved in ensuring size- and charge-selective ultrafiltration (Sellin et al., 2003 [PubMed 12424224]).[supplied by OMIM, Mar 2008]</p> <p>Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>