

Product datasheet for **SC318248**

CDK11A (NM_033529) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CDK11A (NM_033529) Human Untagged Clone
Tag:	Tag Free
Symbol:	CDK11A
Synonyms:	CDC2L2; CDC2L3; CDK11-p46; CDK11-p58; CDK11-p110; p58GTA; PITSLRE
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_033529, the custom clone sequence may differ by one or more nucleotides

```

ATGGGTGATGAAAAGGACTCTTGAAAGTGAAAACCTTAGATGAAATCTTCAGGAAAAG
AAACGAAGGAAGGAACAAGAGGAGAAAGCAGAGATAAACGCTTAAAAAATCTGATGAC
CGGGATCCAAGCGGGATTCCCTTGAGGAGGGGAGCTGAGAGATCACTGCATGGAGATC
ACAATAAGGAAC TCCCCGTATAGAAGAGAAGACTCAATGGAAGACAGAGGAGAAGAAG
GATTCCTTTGGCCATCAAACCACCCCAAAATGTCTCGGAAAGAAAAAGTTCATCAGAA
AAAGATGAAAAGAGAAAAGAAAAATGGAAGCATGCTAGAGTGAAAGAAAAGAGAGCACGAA
CGTCGGAACGACATCGAGAAGAACAGGATAAAGCTCGCCGGAATGGGAAAGACAGAAG
AGAAGGGAAAATGGCAAGGGAGCATTCCAGGAGAGAAAAGGGACCGCTTGGAGCAGTTAGAA
AGGAAGCGGGAGCGGGAGCGCAAGATGCGGGAGCAGCAGAAGGAGCAGCGGGAGCAGAAG
GAGCGCGAGCGGGCGGGAGGAGCGGCCAAGGAGCGGGAGGCCCGCAGGGAAGTGTCT
GCACATCACCGAACGATGAGAGAGGACTACAGCGACAAAGTGAAAGCCAGCCACTGGAGT
CGCAGCCCGCTCGGCCCGCGGGAGCGGTTGAGTTGGGAGACGGCCGAAGCCAGTA
AAAGAAGAGAAAATGGAAGAAAGGGACCTGCTGTCCGACTTACAGGACATCAGCGACAGC
GAGAGGAAGACCAGCTCGGCCGAGTCTCGTCAGCGGAATCAGGCTCAGGTTCTGAGGAA
GAAGAGGAGGAGGAGGAAGAGGAGGAGGAGGAAAGGGAGCACCAAGTGAAGATCAGAGGAG
GAGGAGGAGGAAGAGGAAGAGGAGGAGGAGGAGACCGGCAGCAACTCTGAGGAGGCATCA
GAGCAGTCTGCCGAAGAAGTAAGTGAGGAAGAAATGAGTGAAGATGAAGAACCAGAAAAAT
GAAAACCACTCTTGGTTGTTCCAGAGTACCGTTGACCGGAGATTCCGGGGAGAGTGAA
GAAGCAGAGGAAGAAGTGGTGAGGGAACGCCGAGAGCAGCGCCCTGACAGAGGGCGAC
TATGTGCCCGACTCCCCTGCCCTGTTGCCATCGAGCTCAAGCAGGAGCTGCCAAGTAC
CTGCCGGCCCTGCAGGGCTGCCGAGCGTTCGAGGAGTTCAGTGCCTGAACAGGATCGAG
GAGGGCACCTATGGAGTGGTCTACAGAGCAAAAGACAAGAAAACAGATGAAATTGTGGCT
CTAAAGCGGCTGAAGATGGAGAAGGAGAAGGAGGGCTTCCCGATCACGTCCTGAGGGAG
ATCAACACCATCCTCAAGGCCAGCATCCCAACATTGTCACCGTTAGAGAGATTGTGGTG
GGCAGCAACATGGACAAGATCTACATCGTGATGAACTACGTGGAGCACGACCTCAAGAGC
CTGATGGAGACCATGAAACAGCCCTTCTGCCAGGGGAGGTGAAGACCCTGATGATCCAG
CTGCTGCGGGGGTGAACACCTGCACGACAATGGATCTGCACCGTGACCTCAAGACG
TCCAACCTGCTGCTGAGCCACGCCGCATCCTCAAGTGGGTGATTTTGGGCTGGCGGG
GAGTACGGATCCCCTCTGAAGGCCTACACCCCGTCTGTTGACCCAGTGGTACCCGCC
CCAGAGCTGCTGCTTGGTGCCAAGGAATACTCCACGGCCGTGGACATGTGGTCAGTGGGC
TGCATCTTCGGGGAGCTGCTGACTCAGAAGCCTCTGTTCCCCGGGAATTCGGAATCGAT
CAGATCAACAAAGTGTTCAAGGAGCTGGGGACCCCAAGTGAAGAAAATCTGGCCCGGCTAC
AGTGAGCTCCCAGTAGTCAAAAAGATGACCTTCAGCGAGCACCCCTACAACAACCTCCGC
AAGCGCTTCGGGGCTCTGCTCTCAGACCAGGGCTTCGACCTCATGAACAAGTTCCTGACC
TACTTCCCCGGGAGGAGGATCAGCGCTGAGGACGGCCTCAAGCATGAGTATTTCCGCGAG
ACCCCCCTCCCCATCGACCCCTCCATGTTCCCCACGTGGCCCGCAAGAGCGAGCAGCAG
CGTGTGAAGCGGGGCACCAAGCCGAGCCCTGAGGGAGGCCTGGGCTACAGCCAGCTG
GGTGACGACGACCTGAAGGAGACGGGCTTCCACCTTACCACCACGAACCGGGGGCTCT
GCCGCGGGCCCGGCTTACGCTCAAGTTC

```

Restriction Sites: Please inquire

ACCN: NM_033529

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_033529.2 , NP_277071.2
RefSeq Size:	2620 bp
RefSeq ORF:	2313 bp
Locus ID:	728642
UniProt ID:	Q9UQ88
Cytogenetics:	1p36.33
Protein Families:	Druggable Genome
Gene Summary:	<p>This gene encodes a member of the serine/threonine protein kinase family. Members of this kinase family are known to be essential for eukaryotic cell cycle control. Due to a segmental duplication, this gene shares very high sequence identity with a neighboring gene. These two genes are frequently deleted or altered in neuroblastoma. The protein kinase encoded by this gene can be cleaved by caspases and may play a role in cell apoptosis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2015]</p> <p>Transcript Variant: This variant (4) uses an alternate in-frame splice junction compared to variant 1. The resulting isoform (4) has the same N- and C-termini but is shorter compared to isoform 1.</p>