

## Product datasheet for **SC318254**

### Chromosome X open reading frame 6 (MAMLD1) (NM\_005491) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Chromosome X open reading frame 6 (MAMLD1) (NM_005491) Human Untagged Clone
Tag:	Tag Free
Symbol:	MAMLD1
Synonyms:	CG1; CXorf6; F18; HYSP2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


**Fully Sequenced ORF:** >SC318254 representing NM\_005491.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGATGACTGGAAAAGTCGGCTTGTAAATCAAGAGCATGCTTCCCATTTCGCCATGGTGGGAAATCGT
CAGGAGCCCAGAAAAGCTCCAGGAATCGGGAAGAAGCCCTCGTGGATGGAGGAAGAAGATTTATCTTTT
CTCTACAAGAGCAGCCCAGGAAGAAGCATCAGGGAAGTGTAAAGAGGAGACAAGAAGAAGACCCTTC
CAGTTTCCAGACATGGCTGATGGGGCTACCCTAATAAAATTAAGAGGCCTTGCCTTGAAGATGTCACC
CTTGCAATGGGCCAGGTGCTCATCTAGTACTGCTTGTGCAGAAGTGCAGGTCCCTCCATTGACAATA
AATCTAGCCCTGCGGCTATGGGAGTGGCTGGCCAGTCATTACTGCTGGAGAATAACCCTATGAATGGC
AACATCATGGGCTCACCATTTGTAGTACCACAGACTACAGAAGTGGGACTGAAAGGGCCCACTGTTCTT
TACTATGAGAAAAACAACAGCGTGCCGGCTGTAGACCAGGAGCTCAAGAGCTGTAGAGGAGCTCACC
AAAATTCAAGACCCTTCCAAAATGAGCTAGATCTTGAGAAGATACTGGGACGAAGCCAGAAGAGCCA
CTGGTTTTAGATCATCCCAGGCAACCCTAAGCACAACCTCCAAGCCTTCGGTTCAGATGTCACACTTG
GAGAGCCTGGCTTCCAGCAAGGAGTTTGTCTTCTAGTTGCAGCCAAGTTACTGGCATGTCACCTTCAGATC
CCATCCTCCTCCACAGGGATCAGCTATTCGATTCTTCCACCAGTAAGCAGATAGTGTACCCGAGTTCT
TCAATGGCACAGTCCAAGAGCCAGGTCCAGGCCATGCTCCCTGTGCGCTCTGCCCCCTTACCAGTGCCT
CAGTGGCATCACGCCACCAGCTGAAGGCGTTGGCAGCCAGCAAGCAGGGGTCTGCTACAAGCAGCAA
GGGCCACCCCAAGTGGTCTGGTCTGCCTCCTCCAGGACTCTCTCCACCTTACCGCCAGTGCCATCA
CCACACCACCACCGTGCCACTGCCACCACCACCACCCCAATCAGCCCCAGAGCCTCATGGTGTCC
TGCATGTCGTCCAATACCTTGTGGGTAGCACTCTCCGAGGCTCTCCAATGCCTTACTGCAAGCATG
ACGTCCAGCAGCAATGCTGCCCTGGGGCCCGCCATGCCCTATGCTCCTGAGAAGCTCCCCAGCCCTGCT
CTCACTCAACAGCCGAGTTTCGGCCCTCAGAGCTCCATTCTTGCCAACCTCATGTCTCTACCATAAA
ACCCCTCAAGGACACCTGATGTCTGCTCTGCTGCTGCCAGCAACCCTGGGCCGTCCCCACCCTATCGCCCA
GAGAAGCTCTCTAGCCCAGGCTTCCACAGCAGTCTTCCACCCACAGTGTCCCTGATCCGAAGCCTC
ACTCCCACCAGTAATCTTCTAAGCCAGCAACAGCAGCAGCAGCAGCAGCAGCAGCAAGCAATGTGATC
TTAAGCCATAAGCAGCAACTCATCCAAAACCCTGAGCATGATCATGCAGCAGGGGATGGCAAGCTCC
AGCCCAGGAGCCACGGAGCCATTTACTTTTGGCAACACCAAGCCCTTGTCCATTTTGTCTGAGCCG
GGTCCCCAGAAGATGCCCTCCATGCCTACCACCTTAGGCAGCCTTCCCTGCTCCACTACCTGCAGCAG
CCGACACCAACGCAGGCCTCCTCAGCCACTGCCTCCTCCACGGCCACTGCCACCTTGAGCTGCAGCAG
CAGCAGCAGCAACAGCAGCAGCAGCCTGACCATTCTTATTCTTCTGACAGCAGATGATGCAGCAACCC
CAGCGTTTTAGCGATCAGTGGCCTCAGATTCCATGCCTGCTCTGCCAGACAGGGCTGTTGCCATCTG
TTTGCATGGACTTCTGCAGCTAGCTCGGTGAAGCCCCAGCATCAACACGGGAAGCTTTTCACTAGCAGG
CAAGATCCTCAGCCTGGAGACGTGTACCGTCTAACATTACTCATGTAGACAAAGCCTGCAAGCTTGGG
GAAGCCAGGCACCCCAAGTCCAGCCTCGGGCGACAGCCCCGCTCCTGCCAGGCCCTGGGGAGTGAGTCC
TTCTGCCCAGCAGCTCCTTTGCTCATGAGCTGGCCGAGTCACTCCTCGTACAGCACCTCAGAGGCA
GGCCCTGGGGCAGCTGGGATCCGAAGGCCTGGAGGCAGGTGCCCGCTCCACTACTGCCTAGCTGCGAC
GCCACAGCAAGAGGAACAGAGATCAGGTCTTATGGCAATGACCCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**



**ACCN:** NM\_005491

**Insert Size:** 2325 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_005491.4](#)

**RefSeq Size:** 4607 bp

**RefSeq ORF:** 2325 bp

**Locus ID:** 10046

**UniProt ID:** [Q13495](#)

**Cytogenetics:** Xq28

**MW:** 83.2 kDa

**Gene Summary:**

This gene encodes a mastermind-like domain containing protein. This protein may function as a transcriptional co-activator. Mutations in this gene are the cause of X-linked hypospadias type 2. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Apr 2010]

Transcript Variant: This variant (2) has a shorter 5' UTR and has multiple differences in the coding region, compared to variant 1, one of which results in a translational frameshift. The resulting protein (isoform 2) has a distinct C-terminus and is shorter than isoform 1.