

Product datasheet for **SC318272**

COG4 (NM_015386) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	COG4 (NM_015386) Human Untagged Clone
Tag:	Tag Free
Symbol:	COG4
Synonyms:	CDG2J; COD1; SWILS
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC318272 representing NM_015386.
Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGGGACCAAGATGGCGGACCTTGATTGCGCTCCGAAGCTGTCAGGGGTGCAGCAGCCGTCTGAGGGG
GTGGGAGGTGGCGCTGCTCCGAAATCTCCGTGAGCTCATTGCTCCCTGACAGAGCTGCAGGAGCTG
GAGGCTGTATACGAACGGCTCTGCGGCCGAGGAGAAAGTGGTGGAGAGAGAGCTGGATGCTCTTTGGAA
CAGCAAAACACCAATTGAAAGTAAGATGGTCACTCTCCACCGAATGGGTCCTAATCTGCAGCTGATTGAG
GGAGATGCAAAGCAGCTGGTGGAAATGATCACCTTTACCTGCAACCTGGCTGAGAATGTGTCCAGCAAA
GTTGCTCAGCTTGACCTGGCCAAGAACCCTCTATCAGGCCATTGAGAGAGCTGATGACATCTTGAC
CTGAAGTTCTGCATGGATGGAGTTCAGACTGCTTTGAGGAGTGAAGATTATGAGCAGGCTGCAGCACAT
ACTCATCGCTACTTGTGCTGGACAAGTCGGTCTTGGAGTCTCAGCCGACAGGGCAAGAGGGGAGCATG
ATTGATGCCAACCTGAAATTGCTGCAGGAAGCTGAGCAACGTCTCAAAGCCATTGTGGCAGAGAAGTTT
GCCATTGCCACCAAGGAAGGTGATCTGCCCCAGGTGGAGCGCTTCTTCAAGATCTTCCCCTGCTGGGT
TTGCATGAGGAGGATTAAAGAAAGTTCTCGGAGTACCTTTGCAAGCAGGTGGCCAGTAAAGCTGAGGAG
AATCTGCTCATGGTGTGGGACAGACATGAGTATCGGAGAGCTGCAGTCACTTTGCAGATACACTT
ACTCTTCTGTTTGAAGGGATTGCCCGCATTGTGGAGACCCAGCCAATAGTGGAGACCTATTATGGG
CCAGGGGAGACTCTATACCTGATCAAAATCTGCGAGTGGAAATGTGACAGACAGGTGGAGAAGGTGGTA
GACAAGTTCATCAAGCAAAGGGACTACCACCAGCAGTTCGGCATGTTGAGAACAACCTGATGAGAAAT
TCTACAACAGAAAAATCGAACCAAGAGAACTGGACCCCATCTGACTGAGGTCAACCTGATGAATGCC
CGCAGTGAGCTATACTTACGCTTCTCAAGAAGAGGATTAGCTCTGATTTTGGAGTGGGAGACTCCATG
GCCTCAGAGGAAGTAAAGCAAGAGCACCAGAAGTGTCTGGACAACTCCTCAATAACTGCCTTTTGGAGC
TGTACCATGCAGGAGCTAATTGGCTTATATGTTACCATGGAGGAGTACTTCATGAGGGAGACTGTCAAT
AAGGCTGTGGCTCTGGACACCTATGAGAAGGGCCAGTGCATCCAGCATGGTGGATGATGTCTTCTAC
ATTGTTAAGAAGTGCATTGGGCGGGCTCTGTCCAGCTCCAGCATTGACTGTCTGTGCCATGATCAAC
CTCGCCACCACAGAGCTGGAGTCTGACTTCAGGGATGTTCTGTGTAATAAGCTGCGGATGGGCTTTCCT
GCCACCACCTTCCAGGACATCCAGCGGGGGTGACAAGTGCCGTGAACATCATGCACAGCAGCCTCCAG
CAAGGCAAATTTGACACAAAAGGCATCGAGAGTACTGACGAGGCGAAGATGTCCTTCTGGTGACTCTG
AACACGTGGAAGTCTGCAGTGAACATCTCCACTCTGAAGAAGACACTGGAGAGTACTGCACCAAG
CTTTTCCAGCCAGGCATTGGAGGGGAGCAGGCCAGGCCAAGTTTGCAGCTGCCTTTCTGACTTGCC
GCCGTGTCCAACAAATTCGAGACCTCTTGCAGGAAGGGCTGACGGAGCTCAACAGCACAGCCATCAAG
CCACAGGTGCAGCCTTGATCAACAGCTTTTTCTCCGTCTCCACAACATCGAGGAGGAAGAATTAAT
GACTATGAGGCCAACGACCCCTGGGTACAACAGTTCATCCTTAACCTGGAGCAGCAATGGCAGAGTTC
AAGGCCAGCCTGTCCCGGTCACTACGACAGCCTAACCGGCCTCATGACTAGCCTTGTGCCGTCGAG
TTGGAGAAAGTGGTGTGAAATCCACCTTTAACCGGCTGGGTGGTCTGCAGTTTGACAAGGAGCTGAGG
TCGCTCATTGCCTACCTTACCACGGTGACCACCTGGACCATCCGAGACAAGTTTGGCCGGCTCTCCCAG
ATGGCCACCATCCTCAATCTGGAGCGGGTGACCGAGATCCTCGATTACTGGGACCAATTCGGCCCA
TTGACGTGGCGCCTCACCCCTGCTGAAGTGCGCCAGGTGCTGGCCCTGCGGATAGACTTCCGCAGTGAA
GATATCAAGAGGCTGCGCCTGATAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_015386

Insert Size: 2370 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_015386.2](#)

RefSeq Size: 2838 bp

RefSeq ORF: 2370 bp

Locus ID: 25839

UniProt ID: [Q9H9E3](#)

Cytogenetics: 16q22.1

MW: 89.5 kDa

Gene Summary: The protein encoded by this gene is a component of an oligomeric protein complex involved in the structure and function of the Golgi apparatus. Defects in this gene may be a cause of congenital disorder of glycosylation type IIj. Two transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Aug 2010]
Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).