

Product datasheet for **SC318277**

ZNF606 (NM_025027) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF606 (NM_025027) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF606
Synonyms:	ZNF328
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC318277 representing NM_025027.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCAGCCATCAACCCGTGGGCTCCTGGGGTGCCCTTACGGACCAATCTTGGGGGATGACAGCTGTT
GACCCATGGGCTCCTGGGCTCTGTGTCTCAGTATCCTGCCTGGCACGTGGAGGGAAGCCTGGAGGAG
GGGAGGAGGGCCACTGGGCTCCAGCAGCCAGGTTTCAGGAACCAGTGACCTTCAAGGACGTGGCCGTG
GACTTCACCCAAGAAGAGTGGGGCAGCTGGACCTTGTTCAGAGGACCCTGTACCGTGATGTGATGCTG
GAGACCTATGGTCACCTGCTCTCTGTGGGAAATCAGATTGCCAAGCCTGAGGTCATCTCCCTGTTGGAG
CAAGGAGAAGAGCCGTGGTCAGTGGAGCAGGCATGTCTCAACGCACTTGTCCAGAATGGTGAGAAAT
CTTGAAAGCAAAGCATTGATCCAGCACAGAGCATTTTGAGGAAGAACAATCCCATGGCATGAAGTTG
GAAAGATATATGGGATGATCCTTGGTTCTCCAGTTAGAAGTTTGGGATGTAAGACCAATTAGAA
ATGTACCATGAACAGAGTACAGCTATGAGGCAGATGGTCTTCATGCAAAAGCAAGTACTATCCAG
AGAAGCTCTGAATTCTGTGGACTTGGGGCAGAGTTTAGCCAGAACTTAACTTTGTCCATCTCAGAGA
GTTTCTCAGATAGAACATTTCTATAAGCCTGATACACATGCTCAAAGTTGGAGATGTGACTCAGCCATA
ATGTATGCAGATAAGGTTACCTGTGAAAATAATGATTATGACAAAACCTGTTTATCAGTCCATTCAACCT
ATTTACCTGCAAGAATACAACTGGAGATAATCTTTCAAATGTACTGATGCTGTTAAATCTTTCAAT
CATATAATACATTTTGGTGATCATAAAGGAATTCACACAGGAGAAAACTCTATGAATATAAGGAATGC
CATCAAATCTTTAACCAGAGCCCATCATTTAATGAACACCCAAGGCTTCATGTTGGAGAAAAACAGTAT
AATTACAAAGAATATGAGAATATCTTTATTTCTCATCCTTTATGGAACATCAAAAAATGGTACTGTA
GAGAAAGCGTATAAATACAATGAATGGGAGAAAGTCTTTGGGTATGACTCTTTCCTTACTCAACATACA
AGCACTTACACTGCAGAGAAACCTATGACTACAATGAATGTGGGACGTCTTTCATCTGGAGCTCTTAC
CTTATTCAACATAAGAAACTCATACTGGAGAAAAACCCATGAATGTGATAAATGTGAAAAAGTTTTT
AGGAATCGCTCAGCCCTTACGAACATGAACGACTCACACTGGAATAAAACCCATGAATGTAATAAA
TGTGAAAAGCCTTCAGCTGGAATTCTCATCTTATTGTACATAAGAGAATTCATACAGGAGAAAAACCT
TATGTTTGTAAATGAGTGTGGGAAATCTTCAACTGGAACCTCATCTTATTGGACATCAGAGGACTCAT
ACAGGAGAGAAACCTTTTGAATGTACTGAATGTGGGAAATCATTGAGCTGGAGCTCCCATCTTATTGCC
CATATGAGAATGCATACTGGAGAGAAACCCTTTAAATGTGATGAATGTGAAAAAGCTTTTAGGGACTAC
TCAGCCCTTAGTAAACATGAAAGAACTCATTCTGGAGCAAAACCATATAAATGTACTGAATGTGAAAA
TCCTTCAGCTGGAGCTCCCATCTTATTGCCATCAGAGAACTCACACGGGAGAGAAACCATATAACTGT
CAGGAATGTGGCAAAGCATTGAGAGAACGCTCAGCCCTCACTAAACATGAGATAATTCATTCTGGAATT
AAGCCCTATGAATGTAATAAATGTGAAAAATCCTGTAGCCAGATGGCTCACCTTGTTAGACATCAAAGG
ACTCATACTGGAGAAAAACCCATGAATGCAATAAATGTGAAAAATCCTTCAGTCAGAGCTGTCACCTT
GTTGCTCATCGGAGAATTCACACTGGTGAGAAACCCATATAATGTAATCAGTGTGAAAGATCCTTTAAC
TGATGTTCTCACCTCATTGCACACCGGAGAACTCATCTGGAGAGAAACCATACAGGTGTAATGAATGT
GGGAAAGCATTTAATGAGAGTTCATCCCTTATTGTACACCTAAGAAACCATACTGGAGAAAGCCCTAC
AAATGTAATCATTGTGAAAAAGCATTTTGTAAGAATTCTCCCTTATTATTATCAGAGAAATGCATAGT
GGAGAGAAACGCTTTATATGCAGTGAATGTGAAAAAGCCTTATGTTGCTCACTCAGCCCTACTTCAACAC
CAGAGAAATCACAGTGAAGAGAACTGAATTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_025027
Insert Size: 2379 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_025027.3</u>
RefSeq Size:	4239 bp
RefSeq ORF:	2379 bp
Locus ID:	80095
UniProt ID:	<u>Q8WXB4</u>
Cytogenetics:	19q13.43
Domains:	zf-C2H2
Protein Families:	Druggable Genome, Transcription Factors
MW:	91.8 kDa
Gene Summary:	This gene encodes a zinc finger protein containing a Kruppel-associated box (KRAB) domain at its N-terminus, followed by contiguous C2H2 zinc finger motifs. The encoded protein is a nuclear protein that can act as a transcriptional repressor of growth factor-mediated signaling pathways in a reporter gene assay. This protein has been shown to interact with the SRY-box 9 gene product, and suppresses its transcriptional activity by inhibiting its DNA binding activity. Reduced expression of this gene promotes chondrocyte differentiation. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Dec 2016] Transcript Variant: This variant (1) encodes the longest isoform (1). Both variants 1 and 2 encode the same isoform (1).