

Product datasheet for **SC318281**

CASD1 (NM_022900) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CASD1 (NM_022900) Human Untagged Clone
Tag:	Tag Free
Symbol:	CASD1
Synonyms:	C7orf12; NBLA04196; SOAT
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC318281 representing NM_022900.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCGGCTCTGGCCTACAACCTGGCAAGCGGGAGATCAACCACTACTTCAGCGTGAGGAGCGCCAAG
GTGCTGGCGCTGGTGGCGGTGCTGCTGCTCGCAGCGTGCCACCTCGCCTCCCGCCGTACCGAGGCAAT
GATTCGTGTGAATACCTTCTCTCAAGTGGCAGATTTCTTGAGAGAAAAGTTGGCAACCTCACAGTTGT
ATGATGCATAAATACAAAATCAGTGAAGCAAAGAACTGCCTTGTAGATAAACATATTGCATTTATTGGA
GATTCCAGAATTCGTCAATTGTTTTATTCTTTTGTAAAAATAATTAATCCCAATTCAAAGAAGAAGGA
AATAAGCATGAAAACATTCCTTTTGAAGACAAGACTGCATCAGTTAAAGTGGATTTTCTGTGGCCTCCT
GAAGTTAATGGTCTATGAAACAGTGTATCAAAGTGTGGACTGAGGATCCATTGCAAAGCCACATGTG
ATTGTAGCAGGAGCTGCCACATGGTCCATCAAGATTCACAATGGTAGCAGTGAAGCGTTTCTCAATAT
AAAATGAACATCACCTCCATAGCACCACTTTAGAAAAATTGGCAAAGACTAGTGATGTTTATTGGGTC
TTACAAGATCCTGTTTATGAAGATCTATTAAGTAAAATAGGAAGATGATCACTAATGAGAAGATAGAT
GCTTACAATGAAGCTGCAGTCAATTTGAATAGTAGCACCAGAAATCTAAATCAAATGTTAAGATG
TTCAGTGTTCCAAATTAATTGCTCAAGAAACCATCATGGAATCTTTGGATGGCTTACATCTTCCTGAA
TCGAGCAGAGAACTACTGCAATGATTCTTATGAATGTGTATTGCAATAAGATTTTGAAGCCTGTAGAT
GGGTCCTGTTGTCAACCTCGGCCCTCTGTACTCTCATACAGAAGCTAGCTGCTTGTTTTTTCACTTTA
TCTATTATCGGATATTTAATTTTTTACATAATTCATCGTAATGCTCATCGGAAGAATAAGCCGTGTACT
GATTTGGAAAGTGAGAGAGAAAAGAAAAATATTATCAATACCCCTGTGTCTTATTAGAAATACTTTTA
CAATCTTTCTGCAAACCTGGCCTGATTATGGCATAATTTCTATATGTGTGACCGTGCAAATCTGTTTCATG
AAGGAAAACAAAATTTTATACACATTCATCTTTCTTTTATTCCAATTATCTACATTTTGGTTTTGGGAGTA
TTTTATAATGAAAATACTAAAGAGACTAAAGTATTAATAGAGAACAAACAGACGAATGGAAAGGCTGG
ATGCAACTTGTGATTTTATTTATCACATTTCTGGAGCAAGTACATTTTTGCCTGTATACATGCACATT
CGAGTTCTGGTGTGCATATTTATTTAGACAGGGTATGGGCATTTCTCATACTTTTGGATAAAAAGGA
GATTTTGGAACTATAGAGTATGTCAGGTTTTATTTCTGCTCAATTTCTGGTAGTGGTGTATGTATA
GTAATGGATCGACCTTATCAATTCATTACTTTGTCCCCTTGGTCACTGTATGGTTCATGGTCATATAT
GTTACTTTAGCACTATGGCCACAATAATCCAAAAAAGCAAACGGAATGTTTCTGGCATTTTGGC
TTACTGTTGAACTAGGCTTTTTGCTGTTATTCATATGTTTTTGGCATATCTCAGGGTGCATTTGAG
AAGATCTTTTCTCTTTGGCCATTGTCCAAGTGTTTGAACTGAAAGGGAATGTATATGAATGGTGGTTC
AGATGGAGGTTAGACCGTTATGTAGTTTTCCACGGAATGCTGTTTCTTTATTTATCTGGCTTTGCAG
AAGCGTCAAATACTTTCTGAAGGAAAGGGTGAACCTCTTTTTTCAAACAAAATTTCAAATTTTCTGTTG
TTTATTTCAAGTAGTTTCTTTCTTGACCTATTCCATCTGGGCTAGCAGTTGTAAAAACAAAGCAGAGTGC
AATGAACTCCATCCGTCTGTTTCTGTGGTACAGATTTTAGCCTTCATCCTAATAAGAAACATCCCTGGA
TATGCCGTTTCAGTTTACAGTTCATTTTTGCTTGGTTTGGAAAAATTTCAATTAGAGCTATTTATTGTC
CAGTATACATATGGCTGGCAGCGGACACAAGGGGTATCTTGGTACTGATACCTGAAACCTATGCTC
AACATCATTGTCAGCACTTTTCAATTTGTTTGTGTGGCAGTAAATTTCTCAGATCACTAATGATCTT
GCACAGATTATTTCTAAAGATAACTCATCTCTTTGAAAAGGTTGGCATGTATAGCTGCATTTTTT
TGTGGACTCCTCATCTTATCATCCATTCAGATAAATCAAACATTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_022900

Insert Size: 2394 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_022900.4](#)

RefSeq Size: 3898 bp

RefSeq ORF: 2394 bp

Locus ID: 64921

UniProt ID: [Q96PB1](#)

Cytogenetics: 7q21.3

Protein Families:	Transmembrane
MW:	91.7 kDa
Gene Summary:	O-acetyltransferase that catalyzes 9-O-acetylation of sialic acids (PubMed:20947662, PubMed:26169044). Sialic acids are sugars at the reducing end of glycoproteins and glycolipids, and are involved in various processes such as cell-cell interactions, host-pathogen recognition (PubMed:20947662, PubMed:26169044).[UniProtKB/Swiss-Prot Function]