

Product datasheet for **SC318297**

PPP1R9B (NM_032595) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PPP1R9B (NM_032595) Human Untagged Clone
Tag:	Tag Free
Symbol:	PPP1R9B
Synonyms:	PPP1R6; PPP1R9; SPINO; Spn
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)


Fully Sequenced ORF: >SC318297 representing NM_032595.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGATGAAGACGGAGCCACGGGGCCCGGGGTCCCCTCCGAGCGCCTCCCCGACCCGAGCGCCTAC
GAGGCGGGCATCCAGGCGCTGAAGCCGCCGACGCGCCCGGGCCGACGAGGCACCAAGGGGGCCAC
CACAAAGAAATATGGCTCCAACGTCCACCCGATCAAAAGTATGTTCTGCAGATGGGCACGACGGCGGG
CCCTCGGGCGAGGCGGGCGGCGCGGGCCTGGCCGAGGCCACGGGCGTCCGAGCGGGCTGCGC
CTGTGCTGCCGCGGGCCAGCAGCCTGAACGAGAAGCTGGACCACAGCGCCTGTGAAGCTGGGCACC
AGCGTGTGCGAGCGGTGAGCGCTTCGACTCCAAGCCGCGCCCTCCGCGCAGCCTGCGCCGCCCG
CACCCGCGTCCCGGCTGCAGGAGACGCGGAAGCTGTTCAACGGAGCGCCACGCGCCGAGGCGGC
GACAAGGAGGCCGCGCGCGGGCTGTGAGGAGGAGCGCGCCGCCCTGCAGGACCGGAAGCTGGAC
GTCGTGGTGGCTTCAACGGCAGCACCGAGGCGTGGACAAGCTGGACGCTGACGCCGTGCCCCACG
GTCAGCCAGCTCAGCGCGTCTTCGAGAAGGCCGACTCGAGGACCGGCCCTCCACCGGGGCCGGGCTC
CCAGGGCCGACGGGGTTCGCCAGGTCAACTCGAAGCTGGTCAAGCGGTCCCGGGTGTTCAGCCC
CCGCCCGCGCCGCCGCCCGCCCGTCCGGGGATGCCCCGGCCGAGAAAGAGCGATGCCCCGAGGGCAG
CAGCCCCCGCAGCACCGAGTGGCCCTGCCCGGCCGCCCAAGCCCCGGGAGGTGCGCAAGATTAAG
CCGGTGGAGGTGGAGGAGCGGGGAGTCCGAGGCCGAGTCCGGCGCCCGGGGAGGTGATCCAGGCCGAG
GTTACGGTCCACGCGCCCTGGAGAATGGCAGCACCGTGGCAACTGCAGCCAGCCCCGCGCCGAGGAG
CCAAAGGCCCAAGCGGCCCGGAGAAGGAGGCGGGCGGGTACGCCGCCAGAGAGGGGGTGGCAAT
GGCCGGGCCCGGACGTGGCCCTGAGGAGGTAGATGAATCCAAGAAGGAGGACTTCTCGGAGCGGAC
TTGGTGGACGTGAGCGCTACAGTGGGCTCGGGGAGGACTCTGCGGGCAGTCCCTGGAGGAGGAC
GAAGACGACGAGGAGGATGGGGAGCCCCCTACGAGCCGAGTCCGGGTGCGTGGAGATCCCGGGCTG
TCGGAGGAGGAGGCCAGCCCCGAGCCGGAAGATCCATTTAGCACGGCGCCATCCAAGTGTTCAGC
ACTTACTCCAACGAGGATTACGATCGTCGCAACGAGGATGTGGATCCCATGGCAGCCTCTGCTGAGTAC
GAGCTGGAGAAGCGTGTGGAGAGGTTGGAGCTGTTCCCTGTGGAGCTGGAGAAGGACTCCGAGGGCTG
GGCATCAGCATCATCGGCATGGCGCCGGGGCAGACATGGGCCTGGAGAAGCTGGGTATCTTCGTAAG
ACCGTGACGGAGGTTGGTGGGCCCATCGGGATGGCAGGATCCAGGTGAATGATCTCCTGGTGGAGGTG
GATGGAACAAGTCTGGTGGGAGTACCCAGAGCTTCGCGCGTCTGTGCTCCGGAACACCAAGGGCCGA
GTGCGGTTTATGATTGGCCGGGAGCGGCCGGGAGAGCAGAGCGAAGTGGCCAGCTAATTCAGCAGACT
TTGGAACAGGAGCGATGGCAGCGGGAGATGATGGAGCAGAGATACGCCAGTATGGGGAGGATGACGAG
GAGACGGGAGGATATGCCACTGACGAGGATGAGGAGCTGAGCCCCACGTTCCCGGGTGGTGGATGGCC
ATCGAGGTGTTTGGAGTACGCGGAGAACGAGGATGCACTGTCCCCTGTGGACATGGAGCCCGAGAAGCTG
GTGCACAAGTTCAAGGAGCTCCAGATCAAGCATGCGGTCACTGAGGCAGAGATCCAGCAGCTGAAAAGA
AAGCTGCAGAGCCTGGAGCAGGAGAAGGGGCGCTGGCGGGTGGAGAAGGCGCAGTTGGAGCAGAGTGTG
GAGGAGAACAAGGAGCGCATGGAGAACTGGAAGGCTACTGGGGTGGAGCCAGAGCCTGTGCCAGGCT
GTGGACGAGCACCTGCGGGAGACTCAGGCGCAGTACCAGGCCCTGGAGCGCAAGTACAGCAAGGCCAAG
CGCCTCATCAAGGACTACCAGCAGAAGGAGATCGAGTTCCTGAAAAAGGAGACTGCACAGCGTCCGGTT
CTGGAGGAGTCCGAGCTGGCCAGAAAGGAGGAGATGGACAAGTCCCTGGACAAGATCTCAGAAGTGGAA
GGAAACTTGCAAACTGAGGAATTCGAATTTACTTA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_032595

Insert Size: 2454 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq:	<u>NM_032595.4</u>
RefSeq Size:	4072 bp
RefSeq ORF:	2454 bp
Locus ID:	84687
UniProt ID:	<u>Q96SB3</u>
Cytogenetics:	17q21.33
Protein Families:	Druggable Genome, Phosphatase
MW:	89.3 kDa
Gene Summary:	This gene encodes a scaffold protein that functions as a regulatory subunit of protein phosphatase 1a. Expression of this gene is particularly high in dendritic spines, suggesting that the encoded protein may play a role in receiving signals from the central nervous system. The encoded protein has putative tumor suppressor function and decreased expression has been observed in tumors. [provided by RefSeq, Feb 2014]