

Product datasheet for **SC318298**

PLA2G4D (NM_178034) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PLA2G4D (NM_178034) Human Untagged Clone
Tag:	Tag Free
Symbol:	PLA2G4D
Synonyms:	cPLA2delta
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC318298 representing NM_178034.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGAGAGCCTGTCACCTGGGGGACCACCTGGCCACCCTACCAGGGGGAGGCCTCTACCTGCTGGCAG
CTCACAGTGAGGGTCTGGAGGCGCGAACCTGCGCTGGGCTGACCTGTTGAGTGAGGCCGACCCTTAC
GTGATCCTACAGCTGTCGACCGCAGCTGGAATGAAGTTTAAGACCAAGACGCTCACCGACACCAGTCAT
CCTGTGTGGAATGAGGCCTTCCGTTTCTTATCCAAGTCAGGTCAAGAATGTTCTGAGGCTTAGCATC
TATGATGAGGACTCAGTCACGGAGGATGACATCTGCTTCAAGTTCTCTATGACATCTCAGAAGCTCTC
CCTGGCAAGCTGCTCCGAAAACCTTCTCCAGAGTCCCAGGGAGAGGAGGAGCTGGATGTGGAGTTC
CTGATGGAAGAAACGTCAGATCGCCAGAAAACCTCATACCAACAAAGTCATTGTGCCCGAGAGCTG
TCATGCCTGGATGTGCATCTGGACAGCACAGGGAGCACCGCTGTGGTTCAGATCAGGACAAGCTGGAG
CTGGAGCTGGTGTGAAGGGTCTATGAGGACACACAGACATCCTTCTGGGCACAGCCTCTGCCTTC
CGTTCCACTACATGGCAGCCCTAGAGACAGAGCTGAGCGGGCGCCTGAGGAGCTCCAGAAGCAATGGC
TGGAATGGGGACAACCTAGCTGGGTACCTCACTGTGCCCTGAGGCCCTTGACCATTGGGAAGGAGGTG
ACTATGGATGTTCTGCTCCAAATGCCCCAGGAGTGAGGCTGCAGCTCAAGGCAGAGGGCTGCCCTGAG
GAGCTGGCCGTGCACCTGGGCTTCAATCTCTGTGCAGAGGAGCAGGCCTTCTGAGCAGGAGGAAGCAG
GTGGTGGCCAAGGCCCTGAAGCAGGCCCTGCAGCTGGACAGAGACCTGCAGGAGGATGAGGTACCCGTT
GTGGGCATCATGGCCACAGGAGGAGGTGCCCGGCCATGACCTCACTCTACGGCCACCTATTGGCCTTG
CAGAAGCTGGGCTCCTAGACTGTGTGACCTACTTCAGTGGCATCTCTGGCTCTACGTGGACAATGGCC
CACCTGTACGGGGACCCTGAGTGGTCGCAGAGGGACCTGGAGGGACCTATCAGATACGCCCGGAGCAC
CGGGCAAGGCAAGCTGGAGGTCTTTTCCCGAGAGCCCTGGGAGCTACCGCCGGGAGCTGGAGCTG
CGGGCTGAGCAGGGCACCCACGACCTTTGTGGACCTGTGGGCGCTAGTGTGGAGTCCATGCTGCAC
GGCCAGGTGATGGATCAGAAGCTGTGAGACAGAGAGCCGCCCTGGAACGGGGTCAAGAACCTCTGCC
CTCTACTTGAGCCTCAATGTCAAAGAGAACAATCTGGAGACACTGGACTTCAAGGAGTGGTTGAGTTC
TCCCCCTATGAGGTCGGTTTCTGAAGTACGGGCTTCTGTCCTCTGAGCTCTTCGGCTCCGAGTTC
TTCATGGGACGGCTGATGAGGAGGATCCCGAGCCCGGATCTGCTTTCTGGAAGCCATCTGAGCAAC
ATTTTCTCCCTGAACCTGCTGGATGCCTGGTATGACCTCACCAGTTCTGGGGAGTCTGGAACAGCAC
ATCAAGGACAAGACCAGGAGCTTAGAGAAGGAGCCCTGACCACCTCGGGACCTCCTCGCGCTGGAG
GCCTCGTGGCTGCAGCCAGGCACGGCGCTGGCCAGGCATTTAAAGGCTTCTGACAGGCAGGCCCTC
CACCAGCGCAGCCCAACTTCTCCAGGGCCTCCAGCTGCACCAGGACTACTGTAGCCACAAAGACTTC
TCCACCTGGGCAGACTACCAGCTTGACTCCATGCCAGCCAGCTGACCCCAAGGAGCCCGGCTCTGC
CTGGTGGACGCCGCTACTTCAACACCAGCTCTCCCTCCATGTTCCGGCCAGGCCGAGGCTGGAC
CTCATCTCTCCTTCGACTACTCCCTATCTGCGCCCTCGAGGCACTGCAGCAGACGGAGCTGACTGC
CGGGCCCGGGGGTGCCTTCCCCGGGTGGAACCCAGCCCTCAGGACCAGCACCAGCCAAGGGAATGC
CACCTCTTCTCAGACCCCGCTGCCCGAGGCCCGATCCTGCTGCACTTCCCGTGGTCAATGCCTCC
TCAAGGACCACTCAGCCCCGGTGTCCAGCGCAGCCCCGAGAGCTCCAGGGTGGCCAAGTGGATCTC
ACCGGGCCACCTGCCCTACACCCTGTCCAACATGACCTACAAGGAGGAAGACTTCGAGCGCCTGCTG
CGGCTCAGTGACTACAACGTGCAGACCAGCCAGGTGCCATCCTGCAGGCCCTGAGGACCCGCTGAAG
CACCGGACTCTAGAGGCGAGGCCTCCAAGGACACAGACCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_178034

Insert Size: 2457 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_178034.3](#)

RefSeq Size: 3584 bp

RefSeq ORF: 2457 bp

Locus ID: 283748

UniProt ID: [Q86XP0](#)

Cytogenetics: 15q15.1

Protein Pathways: Fc gamma R-mediated phagocytosis

MW: 92 kDa

Gene Summary: The phospholipase A2 enzyme family, including PLA2G4D, catalyze the hydrolysis of glycerophospholipids at the sn-2 position and then liberate free fatty acids and lysophospholipids (Chiba et al., 2004 [PubMed 14709560]).[supplied by OMIM, Jun 2009]