

## Product datasheet for **SC318316**

### **BRIT1 (MCPH1) (NM\_024596) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	BRIT1 (MCPH1) (NM_024596) Human Untagged Clone
Tag:	Tag Free
Symbol:	BRIT1
Synonyms:	BRIT1; MCT
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC318316 representing NM\_024596.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGGCGGCCCCCATCTGAAAGATGTAGTGGCCTATGTTGAAGTGTGGTCATCCAATGGAACAGAAAAAT
TATTTCAAAGACATTTACAACACAGCTTGTGGATATGGGGCAAAGGTTTCAAAAACTTTTAACAAACAA
GTAACCTACGTTATCTTCAAAGATGGCTACCAGAGCACTTGGGACAAAGCTCAGAAGAGAGGCGTAAAG
CTCGTTTCGTTGCTCTGGGTGGA AAAATGCAGGACAGCTGGAGCACACATTGATGAATCATTGTTCCCT
GCAGCTAATATGAATGAACACTTATCAAGCCTAATTA AAAAAAACGTAATGTATGCAGCCCAAAGAT
TTAATTTTTAAAACACCAGAAAAATGATAAGAGATTT CAGAAGAAAATTTGAGAAAAATGGCTAAAGAGCTA
CAAAGGCAAAAAACAAATCTAGATGATGATGTACCTATTCTCTTATTTGAATCTAATGGTTCATTAATA
TATACTCCACAATTGAAATTAATAGTAGTCACCACAGCGCAATGGAGAAGAGATTACAAGAGATGAAG
GAGAAAAGGGAAAATCTTCCCCACCTCTTCCCAAATGATT CAGCAGTCTCATGATAATCCAAGTAAC
TCTCTGTGTGAAGCACCTTTGAACATTTACGTGATACTTTGTGTTCAGATGAATACTTTGCTGGTGGC
TTACACTCATCTTTTGATGATCTTTGTGAAAAC TCAAGGATGTGAAAATCAGGAAAAGGAAGTTGGAAGGA
TCCATTAATGACATTA AAGTGTGTGTATTCTTCACTTGTATTGAAAGCAAATAATATTCATTCA
TCACCATCTTCACTCACCTCGATAAATCAAGTCTCAGAAAATTTCTGAGTAATCTTTCAAAGGAAGAA
ATAAACTTGCAAAGAAAATATTGCAGGTAAGTAGTCACCCCTGACCAAAAGCAGGCTGCAGGTATGTCT
CAGGAGACGTTTGAAGAGAAGTATCGTTTGTCTCTACCTTATCTTCAACAAAAGGCCACCTTTTGATA
CATTCAAGACCCAGGAGTTCCTCAGTAAAGAGAAAAAGAGTATCATATGGCTCCCATTACCTCCGAG
GAAAAATGCAAGAGAAAGAGGAGCACCAGGAGATCTATCATGCCGAGGCTGCAGCTGTGAGGTGCGAA
GACAGGTCGAGCACGTGGCGGACCTGCCCTGGAGGCTTTAGCTGTGGGAGTCTTCATATGATGAC
TATTTTTACCTGATAATCTTAAGGAAAAGGTATTCAGAGAATCTTCTCCTGAATCTCAGCTGCCATCA
AGCCCTGCTCAGTTGAGCTGCAGAAGTCTTTCTAAGAAGGAGAGAACAAGCATATTTGAAATGTCTGAT
TTTTCTGCGTTGGCAAAAAACCCAGAACAGTTGACATTACCAATTTACAGCAAAAACCATCTCCAGT
CCTCGGAAAACCTGAAAATGGTGAAGGCCGTGCAACTTCGAGTTGCGTGACTTCTGCCCTGAAGAAGCC
CTAAGGTGTTGTAGACAGGCTGGGAAAGAAGACGCATGCCCAGAGGGAATGGCTTTTCTACACCATT
GAGGACCTGCTCTTCAAAGGACATGATGATGATTTAACTCCTTTGGAAGGAAGCCTTGAGAAAATG
AAAGAAGCGGTTGGTCTGAAAAGCACACAGAACAAGGTACCCTTCAAATATCAAACCTCTCTGAA
GGCGAAGCCAGAGTGAACATGAGCCATGTTTTATAGTTGACTGTAACATGGAGAGCTACAGAAGAG
AAGGAAAACCTTACCCGGAGGATACAGTGGAAAGTGTAAAAATAGACCAACAAGGCATGATGTTTTAGAT
GACTCATGTGACGGCTTAAAGGACCTCATCAAACCTCATGAGGAATTGAAGAAAAGTGGGAGAGGCAAA
AAGCCAACAAGAACATTAGTCATGACAAGCATGCCATCTGAAAAGCAGAAATGTCGTATCCAGGTTGTG
GATAAATTGAAAGGCTTTCAATTGCACCAGACGTCTGTGAGACCAGACTCACGTGCTTTCCGGGAAG
CCACTTCGCACCCTGAATGTGCTGCTGGGAATTCGCGTGGCTGCTGGGTTCTCTTATGATTGGGTG
CTATGGTCTTTAGAATTGGGTCAGTGGATTTCTGAGGAGCCGTTCAACTGTCTCACCCTCCCTGCA
GCTCCCCTGTGCCAAGCGAGTGCCACTGTCTGCAGGGCCGTACCGGGAACCCCTTTTGCCGACCAG
CCAGCGATGTTTGTCTCGCTGCCAGCAGCCCCCAGTGGCCAAGCTCTGTGAAC TAGTCCACCTGTGC
GGAGGCCGGTACGCCAAGTCCCCGCCAGGCCAGCATCGTCATCGGGCCCTACAGCGGAAAGAAGAAA
GCCACAGTCAAGTATCTGTCTGAGAAAATGGGTCTTAGATTCCATCACCCAGCACAAAGTCTGTGCCCT
GAAAACCTACCTATTGTCACAATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_024596  
**Insert Size:** 2508 bp

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_024596.4</a>
<b>RefSeq Size:</b>	8056 bp
<b>RefSeq ORF:</b>	2508 bp
<b>Locus ID:</b>	79648
<b>UniProt ID:</b>	<a href="#">Q8NEM0</a>
<b>Cytogenetics:</b>	8p23.1
<b>Domains:</b>	BRCT
<b>MW:</b>	92.8 kDa
<b>Gene Summary:</b>	<p>This gene encodes a DNA damage response protein. The encoded protein may play a role in G2/M checkpoint arrest via maintenance of inhibitory phosphorylation of cyclin-dependent kinase 1. Mutations in this gene have been associated with primary autosomal recessive microcephaly 1 and premature chromosome condensation syndrome. Alternatively spliced transcript variants have been described. [provided by RefSeq, Feb 2010]</p> <p>Transcript Variant: This variant (1, also known as MCPH1-FL) encodes isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>