

## Product datasheet for **SC318317**

### **AKNAD1 (NM\_152763) Human Untagged Clone**

#### **Product data:**

|                           |                                         |
|---------------------------|-----------------------------------------|
| Product Type:             | Expression Plasmids                     |
| Product Name:             | AKNAD1 (NM_152763) Human Untagged Clone |
| Tag:                      | Tag Free                                |
| Symbol:                   | AKNAD1                                  |
| Synonyms:                 | C1orf62                                 |
| Mammalian Cell Selection: | Neomycin                                |
| Vector:                   | pCMV6-Entry (PS100001)                  |
| E. coli Selection:        | Kanamycin (25 ug/mL)                    |



[View online »](#)

**Fully Sequenced ORF:** >SC318317 representing NM\_152763.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGATGAGGCTGATTTTTCAGAACACACGACTTATAAGCAGGAGGATTTGCCTTATGATGGGGACCTC
TCTCAGATTAAGATAGGCAATGATTACAGTTTTACCTCAAAAAAGGATGGCCTTGAAGTCTTAAATCAA
ATTATTTTCATAGCAGATGACCCTCAAGAGAAGGCTATGCATAGTGAGACTTGTGGAAATACAGCTGTG
ACCATACCCCTGGGTAAAATTACTGAAAATGCTGCCAACAAAAAGACGAGAAAGAGAAACAATGCACT
GCAGCTCTTCATATTCCAGCAAATGAAGGAGACGCTTCTAAGTCAAGCATTTCTGATATTTTACTTCAT
CATCTTTCCAAAGAGCCATTCTTAAGAGGTCAAGGCATTGATTGTGAAACCCCTCCAGAGATCTCAAAT
GCCGACAGTTTGAAGAGGAAGCTATTATTAAGTATTATTTTCATGTTATAATAAGAATTCTTGCCCA
AAAGAACAACCCAGAACTCACTGACCAACTCAACCCGAAAAGGGATGGTGAACACAGCAATAAGCCT
GGTCTGCCACCACGACAGAGGAAAATACCTCTGATTTAGAAGGGCCAGTGGCTGCTGGAGATAGCAGC
CATCAAGAAAATGTGAATGTTCTAACTAAAACCAAGGGTCCAGGTGATAAACAAAAAGTTATCAAGGG
CAGTCACCCAGAAACAGCAGACTGAAAAAGCAAATTCAGGCAACACGTTCAAATACGGCCAAGGTCAA
GTTCAATTACCAGCTCCCTGATTTCTCTAAGATTGCTCCCAAAGTAAAAATTCCTAAAAATAGATAATT
AATAAACCACTTGCAATAGCTAAACAAGCCAGCTTTTCTTCCAAGTCGAGAGATAAACCCACTCTTGTC
CAAGATAGTCTAGAAACACGCCTGAGTCAAATGTGTTGAAAAACAACATCAAGAGCAGAAAGGAAA
ATCACTGAACCTTCACAACAAATCCAGATGGAGCCCATAGTACATATCCACCAAGAACTTCTCACAGGA
ATAGAATCTGAGGCAAGTCTCTAAGTTGTCAACACCTCTCAGAAAGGCACTTCTCAAGTTCTTCT
TACATATTTCAAAGATATCCCAAGGGAACAGATGTGTGAGAAGTTGAAAGAACAGACTGATCAACTG
AAGACTAAAGTACAAGAATTTTCAAAGAATAAAACAGGACTCTCCTTACCATTGCAAGACAAGAAG
CTGGTCTGGAGAACTGCAGGACACCTTGAAGTGTGAGCAGAACTTCTGGCCACCAAGGACAAG
CATCTGACTTTGCAGCAGCAAGTCCACAAGCACGAATCAACAATCGTTGGTGACTTTGATCCAGAAAGA
AAAGTGAAGGTGAAATCTTCAAATTGGAGATGCTGCTTGAAGATGTTAAAGAGAAAAATGGATGAAAGC
AAATATACTTCAGCTCCTTCCCTTCTGTGAGTTCTCCAGTTACCCTGGATGACTTGGCTCCACCTTC
TCTTCACTCTCAAATGAGATTCCCAAGGAGCACCCCGGCCACCCTTCTGGGCCTCGAGGCTCAGGTGGG
AGTGAGGTAACAGGGACACCCCAAGGAGGGCCGAGGAGGCCCCAACGAGGAGCTCTGTGAGCTGGCC
CCGACAGCTTACCTAAACGGTCATTATGGAGATGCAGCTGCCAGAACAGCCAGACCAAGTGGCCATG
AGGCTGTCTTCTAACTCTGGGAGGATCCCAACGGCACTCCAAGAAGGCAGGATTGTGCAGAGATGACG
GCACCCAGTCCGAGCTGTGCCTTCTGTGCGAGGCTCCTTGAATGGAAGCAAAACGTGGAGAAAAAGGGC
CACGGAAGGATCAACTGTGGAAGATTTTCAATTGTCTTCATGAAAAGGCACCACACTCAGATTTCACT
CCCAATTCTGATACAGGACACAGCTTCTGTTCTGATTCTGGCACTGAAATGCAGAGTAACAAATGTCAG
GACTGTGGCACTAAGATTCTACCTCCCGAAGAGCCTGCAGGAAAGAACCAACTAAAGAATTTTATTAT
AGATACAACACTCCAGGACAGAATTACTCAAATCATAGCAAAAGAGGTGCCTTTGTCCAGCCCCATTCT
TTAGATGAAAGTAAAACTCTTACCCTCTTTTTTAAACCCAAACGGATCTGTTCTCAGAGAGTGAAT
TCAAAATCCTTTAAAGGTGAACATGAGCCACACCAGGAAAAAAAACCTTCAGGCTTTCATGACCTAC
AGCTCAGACCCCTGCAACACCCCTACCCCATTTCTACTCCTGCAGGATTTCTGGAAGCAAAATCCTTATGT
GACTTTGATAGCACTGAAGAAATAAAATCTGAGATTTTAACTCGGCTTTGGATCATGCCCTAAGGACA
GCAACCATTTTGAAGAACTACAGATCAAATGATTAACGATTGCAGAAGACCTTGCTAAAGCACAG
AGGTGGAGGAATCGACTGAAATAGAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** Sgfl-MluI  
**ACCN:** NM\_152763  
**Insert Size:** 2511 bp

|                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|-------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).                                                                                                                                                                                                   |
| <b>OTI Annotation:</b>        | This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.                                                                                                                                                                                                                                                                               |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                                                             |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>                          |
| <b>RefSeq:</b>                | <u><a href="#">NM_152763.4</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>RefSeq Size:</b>           | 3063 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>RefSeq ORF:</b>            | 2511 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Locus ID:</b>              | 254268                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>UniProt ID:</b>            | <u><a href="#">Q5T1N1</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Cytogenetics:</b>          | 1p13.3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>MW:</b>                    | 92.9 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Gene Summary:</b>          | <p>This gene encodes a protein which contains a domain found in an AT-hook-containing transcription factor. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2012]</p> <p>Transcript Variant: This variant (1) encodes the functional protein. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p> |