

Product datasheet for **SC318320**

SMPD4 (NM_017751) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SMPD4 (NM_017751) Human Untagged Clone
Tag:	Tag Free
Symbol:	SMPD4
Synonyms:	NEDMABA; NEDMEBA; NET13; NSMASE-3; NSMASE3; SKNY
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_017751, the custom clone sequence may differ by one or more nucleotides

```

ATGACGACTTTTCGGCGCCGTGGCGGAATGGCGGCTTCCATCTCTGAGGCGAGCGACGCTA
TGGATCCCACAGTGGTTTGTAAAGAGGCCATTTCAACTCTCCACTGGAGGCTGCTATG
GCGTTCCTCACCTGCAGCAGCCAGCTTTCTACTGGCTAGCCTGAAAGCTGACTCTATA
AATAAGCCCTTTGCACAGCAGTCCCAAGACTTGGTTAAAGTCATTGAGGACTTCCAGCA
AAGGAGCTGCACACCATTTCCCATGGCTGGTAGAAAGCATTTTTGGCAGCCTAGATGGT
GTCTCGTTGGCTGGAACCTCCGCTGCTTACAGGGGCGCGTGAATCCTGTGGAGTACAGC
ATCGTGATGGAATTTCTCGACCCTGGTGGCCCAATGATGAAGTTGGTTTATAAGCTTCAA
GCTGAAGACTATAAGTTCGACTTTCCTGTCTCCTACTTGCCTGGTCTGTGAAGGCGTCC
ATCCAGGAGTGCATCCTCCCTGACAGTCTCTGTACCACAACAAGGTCCAGTTCACCCCT
ACTGGGGGCTTGGTCTGAACCTTGGCCCTGAATCCGTTGAGTATTACATATTCTTCTTT
GCCTTGAGCCTCATCACTCAGAAGCCACTTCTGTGTCCCTCCACGTCCGTACTTCAGAC
TGTGCCTATTTTCATCCTGGTGGACAGGTACCTGTATGGTTTCTGCCACCGAAGGCAGT
GTGCCCCACCCTCTCCTCCAGCCAGGGGGGACCAGCCCTCACCACCTCCCAGGACA
CCAGCCATACCCTTTGCTTCTATGGCCTCCACCACACTAGCCTCTAAAGCGACACATC
TCTCATCAGACGCTGTGAATGCAGACCCCGCCTCCCACGAGATCTGGAGGTCAGAACT
CTGCTCCAGGTTTTTGTGAAATGTGGCTTCATCACTATTCTTGGAGATGTATCAAAAA
ATGCAGTCCCCTCATGCCAAGGAGTCGTTACAGCCTACTGAGGAGCATGTGTTGGTGGTG
CGCCTGTCTGAAGCACCTGCACGCCTTTGCCAACAGCCTGAAGCCAGAGCAGGCTCA
CCCTCCGCCCACTCCCACGCCACCAGCCCCCTGGAGGAGTTCAAACGGGCTGCTGTCCCG
AGGTTTCGTCAGCAGAACTCTACCTCTTCTTGCAGCATTGCTTTGGCCACTGGCCCTG
GACGCTCGTTCAGAGCTGTCTGGAGATGTGGCTGAGTACCTGCAGCCGTGGCGGTAC
GGCCTGACAAGCAGGCTCCGGGACGCACTCCAGCCCGGTGTGTGCGGAGAAATGG
GCACCTTTTGTCCAGGAGAACCTGCTGATGTACCCAAGTTGTTTGTGGGCTTTCTGAAC
CGCGCGCTCCGCACAGACCTGGTCCAGCCCAAGCACGCGCTCATGGTGTCCGAGTGGCC
AAAGTCTTTGCCAGCCCAACCTGGCTGAGATGATTAGAAAAGGTGAGCAGTATTCTCTG
GAGCCAGAGTGGTATCCCCACCGCCAGCACCAGTCTTACGGCCCCACATTCCTACT
GGGAGCTTCTGTACCCTGGCCACCAGCGGTCACTGATGCCTCCTTCAAGGTGAAGAGC
CACGTCTACAGCCTGGAGGGCCAGGACTGCAAGTACACCCCGATGTTTGGGCCGAGGCC
CGCACCTGGTCTGCGCTCGCTCAGCTCATCACACAGGCCAAACACACAGCCAAGTCC
ATCTCCGACCAAGTGTGCGGAGAGCCCGGTGGCCACTCCTTCTCATGGCTGGGCTTT
AGCTCCATGGACACCAATGGCTCCTACACAGCCAACGACCTGGACGAGATGGGGCAAGAC
AGTGTCCGGAAGACAGATGAATACCTGGAGAAGGCCCTGGAGTACCTGCGCCAGATATTC
CGGCTCAGCGAAGCGCAGCTCAGGCAGTTCACACTCGCCTTGGGCACCACCCAGGATGAG
AATGGAAGAAAGCAACTCCCCGACTGCATCGTGGGTGAGGACGGACTCATCCTTACGCC
CTGGGGCGGTACCAGATCATCAATGGGCTGCGAAGTTTGAATTTGAGTACCAGGGGGAC
CCGGAGCTGCAGCCATCCGGAGCTATGAGATCGCCAGCTTGGTCCGCACACTCTTTAGG
CTGTCGTCGCCATCAACCACAGATTTGCAGGACAGATGGCGGCTCTGTGTTCCCGGAT
GACTTCTCGGCAGCTTCTGTGCTACCACTCACAGAACCTGGGCTGGCCAGCAGGCAC
CTGCTGAGCCCTGTGGGCGGAGGCAGGTGGCCGGCCACACCCGCGGCCCCAGGCTCAGC
CTGCGCTTCTGGGCAGTTACCGGACGCTGGTCTCGTGTGCTGGCCTTCTTGTGGCC
TCTCTGTTCTGCGTGGGGCCCTCCCATGCACGCTGCTGCTACCCTGGGCTATGCTCCT
TACGCCTCTGCCATGACTGCTGACCGAGCGGGGAAGCTGCACCAGCCC

```

Restriction Sites: Please inquire
ACCN: NM_017751
Insert Size: 3833 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_017751.2</u> , <u>NP_060221.2</u>
RefSeq Size:	3833 bp
RefSeq ORF:	2514 bp
Locus ID:	55627
UniProt ID:	<u>Q9NXE4</u>
Cytogenetics:	2q21.1
Protein Families:	Transmembrane
Protein Pathways:	Metabolic pathways, Sphingolipid metabolism

Gene Summary:

The protein encoded by this gene is a sphingomyelinase that catalyzes the hydrolysis of membrane sphingomyelin to form phosphorylcholine and ceramide. This gene is activated by DNA damage, cellular stress, and tumor necrosis factor, but it is downregulated by wild-type p53. The encoded protein localizes to the endoplasmic reticulum and Golgi network.

[provided by RefSeq, Mar 2017]

Transcript Variant: This variant (1) lacks an alternate in-frame exon in the central coding region, compared to variant 2, resulting in an isoform (1) that is shorter than isoform 2. CCDS Note: This CCDS representation uses the 5'-most in-frame start codon. While this start codon is well-conserved, it should be noted that other than human and chimp, all other species have indels that follow this start codon. Other primates and rodents have frameshifting indels, and thus cannot use this start codon. However, some mammalian species, including cow, pig and dog, contain additional compensating indels that bring them back in frame; these species could therefore use the start codon represented in this CCDS. An alternative downstream start codon, which is more widely conserved and has a stronger Kozak signal, also exists. It is possible that leaky scanning by ribosomes would allow the downstream start codon to be used, at least some of the time. The use of the downstream start codon would result in a protein that is 39 aa shorter at the N-terminus. There is no experimental evidence showing which start codon is preferentially used in vivo.