

Product datasheet for **SC318353**

Protein Kinase D2 (PRKD2) (NM_016457) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Protein Kinase D2 (PRKD2) (NM_016457) Human Untagged Clone
Tag:	Tag Free
Symbol:	PRKD2
Synonyms:	HSPC187; nPKC-D2; PKD2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

Fully Sequenced ORF: >SC318353 representing NM_016457.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGCCACCGCCCCCTTTATCCCGCCGGGCTCCCTGGCTCTCCCGGGCCGGGGTCTCCTCCGCCCCC
GGCGGCTAGAGCTGCAGTCGCCGCCACCGTACTGCCAGATCCCGGCCCGGGTTCGGGGTCTCC
TTTCACATCCAGATCGGGCTGACCCGCGAGTTCGTGCTGTTGCCCGCCGCTCCGAGCTGGCTCATGTG
AAGCAGCTGGCCTGTTCCATCGTGGACCAGAAGTTCCTGAGTGTGGCTTCTACGGCCTTACGACAAG
ATCCTGCTTTTCAAACATGACCCACGTCGGCCAACTCCTGCAGCTGGTGCCTCGTCCGGAGACATC
CAGGAGGGCGACCTGGTGGAGGTGGTGTGTCGGCCTCGGCCACCTTCGAGGACTCCAGATCCGCCCG
CACGCCCTCACGGTGCACCTCTATCGGGCGCTGCCTTCTGTGATCACTGCGGGGAGATGCTCTTCGGC
CTAGTGCGCCAGGGCCTCAAGTGCATGGCTGCGGGCTGAACTACCACAAGCGCTGTGCCTTCAGCATC
CCCAACAAGTGTAGTGGGGCCGCAAACGGCGCTGTCATCCACGTCTCTGGCCAGTGGCCACTCGGTG
CGCCTCGGCACCTCCGAGTCCCTGCCCTGCACGGCTGAAGAGCTGAGCCGTAGCACCACCGAACTCTG
CCTCGCCGTCCCCGTATCCTTCTCTCTCTCTGCTCATCGTATACGGGCCGCCCAATTGAGCTG
GACAAGATGCTGCTCTCAAGGTCAAGGTGCCGCACACCTTCTCATCCACAGCTATACAGGCCACC
GTTTGCCAGGCTTGAAGAACTCCTCAAGGGCCTTCCGGCAGGGCCTGCAATGCAAAGACTGCAAG
TTAACTGTACAAACGCTGCGCCACCCGCTCCCTAATGACTGCCTGGGGGAGGCCCTTATCAATGGA
GATGTGCCGATGGAGGAGGCCACCGATTCAGCGAGGCTGACAAGAGCGCCCTCATGGATGAGTCAGAG
GACTCCGGTGTATCCTGGCTCCCACTCAGAGAATGCCTCCACGCCAGTGAGGAGGAGGAAGGCGAG
GGAGGCAAGGCCAGAGCTCCCTGGGGTACATCCCCTAATGAGGGTGGTGAATCGGTGCGACACAG
ACGGGAAATCCAGCACACGCTGCGGGAGGGTGGGTGGTTTATTACAGCAACAAGGACACGCTGAGA
AAGCGGCACTATTGGCGCTGGACTGCAAGTGTATCACGCTCTCCAGAACAACAGACCAACAGATAC
TATAAGGAAATCCGCTGTGAGAATCCTCACGGTGGAGTCCGCCAGAACTTCAGCCTTGTGCCGCCG
GGCACCAACCCACACTGCTTTGAGATCGTCACTGCCAATGCCACCTACTTCGTGGGCGAGATGCCTGGC
GGGACTCCGGGTGGGCAAGTGGGCAGGGGGCTGAGGCCGCCCGGGGCTGGGAGACAGCCATCCGCCAG
GCCCTGATGCCGTATCCTTCAGGACGCACCCAGCGCCCGAGGCCACGCGCCCCACAGACAAGCTTCT
CTGAGCATCTCTGTGTCCAACAGTCAGATCCAAGAGAATGTGGACATTGCCACTGTCTACCAGATCTTC
CCTGACGAAGTGTGGGCTCAGGGCAGTTTGGAGTGGTCTATGGAGGAAAACACCGGAAGACAGGCCGG
GACGTGGCAGTTAAGGTCAATTGACAAATGCGCTTCCCTACCAAGCAGGAGAGCCAGCTCCGGAATGAA
GTGGCCATTCTGCAGAGCCTGCGGCATCCCGGGATCGTGAACCTGGAGTGCATGTTTCGAGACGCCTGAG
AAAGTGTGGTGGTGGATGGAGAAGTGCATGGGGACATGTTGGAGATGATCCTGTCCAGTGAAGGGC
CGGCTGCCTGAGCGCCTCACCAAGTTCCTCATACCCAGATCCTGGTGGCTTTGAGACACCTTCACTTC
AAGAACATTGTCCACTGTGACTTGAAACAGAAAACGTGTTGCTGGCATCAGCAGACCCATTTCTCAG
GTGAAGCTGTGTGACTTTGGCTTTGCTCGCATCATCGGCGAGAAGTCTTCCGCCGCTCAGTGGTGGGC
ACGCCGGCCTACCTGGCACCCGAGGTGCTGCTCAACCAGGGCTACAACCCTCGCTGGACATGTGGTCA
GTGGGCGTGATCATGTACGTGAGCTCAGCGCACCTTCCCTTCAACGAGGATGAGGACATCAATGAC
CTCATCAACAACCTGCTGAGGTGAAGATGCGCAACGCTACAGCGTGGACAAATCTCTCAGCCACCCC
TGGTTACAGGAGTACCAGACGTGGCTGGACCTCCGAGAGCTGGAGGGGAAGATGGGAGAGCGATACATC
ACGCATGAGAGTGACGACGCGCTGGGAGCAGTTTGCAGCAGAGCATCCGCTGCCTGGGTCTGGGCTG
CCCACGGACAGGGATCTCGGTGGGGCTGTCCACCACAGGACCACGACATGCAGGGGCTGGCGGAGCGC
ATCAGTGTCTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:


ACCN: NM_016457

Insert Size: 2637 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_016457.4](#)

RefSeq Size: 3617 bp

RefSeq ORF: 2637 bp

Locus ID: 25865

UniProt ID: [Q9BZL6](#)

Cytogenetics: 19q13.32

Protein Families: Druggable Genome, Protein Kinase

MW: 96.7 kDa

Gene Summary: The protein encoded by this gene belongs to the protein kinase D (PKD) family of serine/threonine protein kinases. This kinase can be activated by phorbol esters as well as by gastrin via the cholecystokinin B receptor (CCKBR) in gastric cancer cells. It can bind to diacylglycerol (DAG) in the trans-Golgi network (TGN) and may regulate basolateral membrane protein exit from TGN. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]
Transcript Variant: This variant (1) represents the longest transcript. Variants 1, 2, and 3 encode the same isoform (A).