

## Product datasheet for **SC318372**

### CHD1L (NM\_004284) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CHD1L (NM_004284) Human Untagged Clone
Tag:	Tag Free
Symbol:	CHD1L
Synonyms:	ALC1; CHDL
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)




[View online »](#)

**Fully Sequenced ORF:** >SC318372 representing NM\_004284.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGCGCGCGGGCGCTACTAGCCGCGGGGCAAGCGCCTGGCTTCTTACTGCGGCTTCATACTGAG
GGCCGAGCCGAGCGCGGGTGCAGGAGCAGGACTTACGGCAGTGGGGCTGACAGGGATTACACCTA
CGCTCTTACCAGCTGGAGGGAGTAACTGGCTCGCCAGCGCTTCCATTGTCAGAATGGCTGATCCTG
GGAGATGAGATGGGCTGGGGAAGACCTGCCAGACTATTGCTCTTCTATTATTTGGCAGGAAGATTA
AATGATGAAGGGCCATTTCTGATTCTTTGTCCCTTGCTGTTTTGAGCAACTGGAAGAAGAAATGCAG
AGATTTGCTCCAGGTCTTCTGTGTAACATATGCAGGCGACAAGGAGGAAAGAGCCTGCCTTCAGCAA
GACCTGAAACAGGAGTCACGTTTTCATGTGCTACTGACTACCTATGAGATTTGCTGAAAGATGCATCA
TTTCTAAAATCATTCCCTTGGAGTGTCTTGTGTGGATGAAGCTCACAGGTTGAAAACCAAAGCTCC
CTGCTGCATAAGACCTTGTGAGAGTTCAGTAGTCTTCCAGTCTCCTGTTGACCGGAACCTCCATCCAG
AACAGCCTCCAAGAGCTCTACTCCCTCCTCAGTTTTGTGGAGCCTGATCTTTTTCCAAGGAAGAGGTG
GGAGATTTTATCAACGCTACCAGGATATTGAGAAAGAATCTGAGTCAGCAAGTGAAGTGCACAACTC
TTGCAGCCATTTCTGCTGAGGCGAGTGAAGCTGAGGTAGCTACAGAGCTTCCCAAGAAGACAGAAGTA
GTGATATACCATGGCATGTCAGCATTGCAGAAGAAATACTACAAGGCCATTTTGATGAAAGACCTAGAT
GCATTTGAAAAAGAGACGGCAAAGAAAGTTAACTACAGAACATTTTGTCCAGCTTCGAAAAGTGTGTG
GATCACCCATATTTGTTTGTGAGTGTGGAGCCGGAGCCTTTTGAAGTTGGAGACCACCTGACTGAGGCT
AGTGGGAAGCTTACCTGCTGGATAAGCTACTAGCATTCTGTATTCTGGGGCCATCGGGTTTTACTT
TTCTCCAAATGACCCAGATGTTGGATATTCTCCAAGACTATATGGATTACAGAGGCTACAGCTATGAG
CGTGTGGATGGTCTGTGAGAGGAGAAGAGACACTTGGCCATTAAGAAGCTTTGGACAGCAGCCCAT
TTGTTTTTCTCCTGAGTACTAGGGCAGGTGGAGTTGGCATGAACTTAACAGCAGCAGATACTGTGATT
TTTGTGACAGTGACTTTAATCCTCAGAATGACTTGAAGCAGCTGCCAGGGCTCATCGATTGGCCAA
AACAAAGTCTGTTAAAGTTATTCGGCTGATTGGTCGAGACACTGTGGAAGAAATAGTCTATAGGAAAGCA
GCCTCCAAACTGCAGCTACCAACATGATCATAGAAGGAGGCCATTTTACTCTGGGAGCCAGAAAACCC
GCTGCCGATGCTGACCTCCAGTTGAGTGAGATACTCAAATTTGGTTTGGATAAACTGCTGGCCTCTGAG
GGGAGCACCATGGATGAAATAGACCTGGAGTCCATCCTGGGAGAAAACAAAAGATGCCAGTGGGTCTCT
GATGCCTTGCTGCAGCAGAAGGAGGGAGCAGAGATCAAGAGGAAGGAAAAATCATATGTACTTATTT
GAAGGTAAGATTATTCTAAAGAGCCCAGTAAGGAAGACAGAAAATCATTGAACTGGTAAACCTT
CAGAAAACCCTTTTGGAGAAAGCTAGTCAAGAGGGCCGATCACTCCGAAATAAAGGCAGTGTCTCATC
CCAGGCCTTGTGGAGGGATCTACCAAAAGGAAGCGGGTCTGAGTCCAGAAGAGCTGGAGGACAGACAG
AAGAAAAGACAAGAAGCAGCTGCCAAGAGAAGGAGACTCATAGAGGAGAAGAAGAGGCAAAAAGGAAGAG
GCTGAACATAAGAAAAAGATGGCCTGGTGGGAATCCAACAATTACCAGTCTTCTGCCTGCCCTCTGAG
GAGAGCGAGCCAGAGGACCTTGAGAAATGGGGAAGAGAGCTCTGCTGAGCTGGATTACCAAGACCCAGAT
GCTACTTCCCTCAAGTACGTTAGTGGTATGTCACCCACCCTCAGGCTGGGGCCGAGGATGCTCTCATT
GTGCACTGCGTAGATGACTCTGGCCACTGGGCAGAGGTGGTTTTATTTACAGCTCTGAAAAGCCATCC
GCTGAGCCAAGAAAAATATAGAGCTGGCTGGGAAAATGAAAGACCTGAGTTTGGGAGGTGTCCTTTTA
TTTCTGTTGATGATAAAGAATCAAGAAACAAAGGGCAAGATTTGTTGGCCTTGATTGTGGCTCAGCAT
CGTGATCGTTCCAATGTCTGTCTGGCATTAAAGATGGCAGCCCTAGAAGAGGGCCTGAAGAAGATATTT
TTAGCAGCAAAAAGAAGAAAGCAAGTGTTCATCTCCACGATTGGACATGCCACGAAAGGTTTTAAC
TGGTATGGTACTGAGCGACTTATTCGAAACATCTGGCTGCAAGAGGCATCCCAACTTACATATATTAT
TTTCTAGAAGCAAGTCTGCTGTCTTCAATCACAGTCTTCACTTCTCCTCAAGACAGCTGGTGCCT
TAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

**Restriction Sites:** Sgfl-Mlul

**Plasmid Map:**


**ACCN:** NM\_004284

**Insert Size:** 2694 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_004284.4](#)

**RefSeq Size:** 3036 bp

**RefSeq ORF:** 2694 bp

**Locus ID:** 9557

**UniProt ID:** [Q86WJ1](#)

**Cytogenetics:** 1q21.1

**Domains:** SNF2\_N, DEAD

**MW:** 101 kDa

**Gene Summary:** This gene encodes a DNA helicase protein involved in DNA repair. The protein converts ATP to add poly(ADP-ribose) as it regulates chromatin relaxation following DNA damage. Overexpression of this gene has been linked to several types of cancers. [provided by RefSeq, Feb 2017]

Transcript Variant: This variant (1) encodes the longest isoform (1). Sequence Note: The RefSeq transcript and protein were derived from genomic sequence to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on alignments.