

Product datasheet for **SC318412**

FBXO10 (NM_012166) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FBXO10 (NM_012166) Human Untagged Clone
Tag:	Tag Free
Symbol:	FBXO10
Synonyms:	FBX10; PRMT11
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF:

>SC318412 representing NM_012166.

Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGGCTGGTGGCCTCCCTTGGAGCTGTGGCGCATGATCTTAGCCTACTTGACCTTCCCGACCTG
GGCCGCTGCAGCCTGGTATGCAGGGCCTGGTATGAAGTATCCTCAGTCTCGACAGCACCCGCTGGCGG
CAGCTGTGTCTGGGTTGCACCGAGTGCCGCCATCCCAATTGGCCCAACCAGCCAGATGTGGAGCCTGAG
TCTTGAGAGAAGCCTTCAAGCAGCATTACCTTGCATCCAAGACATGGACCAAGAATGCCTTGGACTTG
GAGTCTTCCATCTGCTTTTCTCTATTCCGCCGAGGAGGGAACGACGTACCCTGAGTGTGGGCCAGGC
CGTGAGTTTGACAGCCTGGGCAGTGCCTTGGCCATGGCCAGCCTGTATGACCGAATTGTGCTCTTCCCA
GGTGTGTACGAAGAGCAAGGTGAAATCATCTTGAAGGTGCCTGTGGAGATTGTAGGGCAGGGGAAGTTG
GGTGAAGTGGCCCTGCTGGCCAGCATTGATCAGCACTGCTCAACCACACGCTGTGCAACCTCGTCTTC
ACGCCAGCCTGGTTCTCACCCATCATGTATAAGACAACATCAGGTCACGTCCAGTTTGACAACTGCAAC
TTTGAGAACGGGCACATCCAGGTCCATGGCCCGGGTACTTGCCAAGTGAAGTTCTGTACCTTCAAAAAC
ACCCATATCTTCTGCACAACGTGCCCTGTGTCTGCTGAAAACGTGAATTTGTGGGCAGTGA AAAAC
AACTCTGTGACTGTTGAGGGTCACCCATCTGCAGATAAGAACTGGGCCTACAAGTATCTACTAGGGCTT
ATCAAGTCCTCACCCACTTTTCTCCCCACAGAGGACTCTGACTTTTAAATGTCCCTGGACCTAGAGAGC
CGGGACCAAGGCTGGAGCCCAAAGACCTGTGACATTGTTATCGAGGGCAGCCAGAGCCCTACCAGCCCA
GCCTCTAGCTCCCAAAGCCAGGCTCCAAGGCTGGCTCACAGGAGGCAGAGGTGGGTAGTGATGGTGAA
AGGGTGGCCAGACCCCGGACAGCAGCGATGGAGGCTGAGTCCAGCGGTGAGGATGAAGATGAGGAC
CAGCTGATGTACAGACTATCCTACCAAGTGCAGGGCCACGCCCTGTATTGGGGGGCTCATTCTGGGC
CCACCTCTACCAGGAGCATCCATTACAGTGCAGCTGCCTAGTGCTGAAGTCACTGCAGCAGGAGCTG
CAGAAGGATAAGGAGGCCATGGCACTGGCCAACCTCCGTGCAGGGCTGCCTCATCCGCAAGTGCCTCTTC
CGGGACGGGAAGGGAGGCGTCTTCGTCTGCTCCACGGCAGAGCCAAGATGGAAGGAAACATCTTCGG
AACCTGACTTACGCAGTGCGGTGTATACATAATAGCAAGATCATCATGCTCAGGAACGACATTTACCGC
TGCCGAGCGTCAGGCATCTTCTTCGCTTGGAGGGCGGTGGCTTGATTGCCGGCAACAACATTTACCAC
AATGCAGAGGCTGGTGTAGACATCCGAAAAAGTCCAACCCACTCATACTGTGTAACCAGATCCACCAT
GGCCTTCGCTCTGGCATTGTCTGCTTGGCAATGGGAAAGGCATCATCCGGAACAATCAAATCTTTTCC
AATAAGGAGGCTGGCATTACATCCTGTACCACGGAAACCCGTTGTGAGCGGGAACACATCTTCAAG
GGCCGTGCTGCCGCATAGCAGTGAATGAGAACGGCAAGGCCTCATCACAGAAAATGTATCCGTGAG
AATCAGTGGGAGGTGTGGACATCCGCCGTGGAGGGATCCCCGTTCTCAGGAGTAACCTCATCTGCTTT
GGCTATTAGATGGTGTGGTGTGGGAGACGAAGGCAAGGCCTCATAGAAGGAAATACCATCTACGCT
AACAAGGGCTGTGGTGTGGATGATGTCGTCCAGCCTCCCCATGTACCCAGCAACCACGTGAGTAC
AATGGCCTGTATGGAGTGGCAGTATTTAGCCAGAAGGATGGCTCCAGCGAGTTACCTCGAGGCCACAGG
GCTCAAGAGAAGTTACGCGAGGATGGGGACGCCATCCTCTGGGAGACAGAGCTGGAGAAGGAGGACGAC
CCACTGCGCCGGCCCATCACCATAGCTCTTGTTGAGTCTAACAGTATTAATCACAAATGGAGCCTCAGGA
CTCTATGTCCAGAGCAGCGAGGCACTGCATGTCATACCAATGTGATCCACGCAATGGGACAGAGGC
ATTACTGTGGCCAGAGCAGCAACCCACCCGAGTGGCCAACAACAGCATCTCTGCAACCGGCAAGT
GGGGTCAAGTTGAGGCCAGTGCAAAAGTGGAGCTCCGGGGCAATGGTATCTATGACAACAGAGGCCAC
GGCATTATACCAAGGGCGACAGCACCATCGTCATTGAAAACGATATCATTGGCAACCGGGGACGCGG
CTGCAGCTGCTGCCAGGTCCGACACTAAAGTAATAAAGAACCGGATCCACTCGTTCCGGGCCTACGGC
ATCGCCGTGCGGGGCCGTGCCAAGGCCCTGGTGCAGGAAAACATCATCTCCAGGGCAAAACAGTAAG
ACCATCTTTCAGCAGATCTCAAACAACCGAGAATGCATCATGAAAACAACAAGTTCTTGGTCTTCAAG
AAAAAGTCTGATACGTGGCGCTGGTGAACCCACCAGCACGGCCCCACCTGAAAATTCTCTCAGACGT
CCCTCGGCAGCCACAATGGGCAGAAGGTGACAGCCATGGCAACGAGGATCACAGCCCGGTGGAAGGT
GGTTACCACAGCAACCGCAGTGTCTTCTGCACCATCCTGTA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:

Sgfl-MluI

ACCN:	NM_012166
Insert Size:	2871 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_012166.2</u>
RefSeq Size:	4592 bp
RefSeq ORF:	2871 bp
Locus ID:	26267
UniProt ID:	<u>Q9UK96</u>
Cytogenetics:	9p13.2
Protein Families:	Druggable Genome
MW:	105.2 kDa
Gene Summary:	Members of the F-box protein family, such as FBXO10, are characterized by an approximately 40-amino acid F-box motif. SCF complexes, formed by SKP1 (MIM 601434), cullin (see CUL1; MIM 603134), and F-box proteins, act as protein-ubiquitin ligases. F-box proteins interact with SKP1 through the F box, and they interact with ubiquitination targets through other protein interaction domains (Jin et al., 2004 [PubMed 15520277]).[supplied by OMIM, Mar 2008]