

Product datasheet for **SC318417**

VPS35L (NM_020314) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	VPS35L (NM_020314) Human Untagged Clone
Tag:	Tag Free
Symbol:	VPS35L
Synonyms:	C16orf62; EC97
Vector:	<u>pCMV6 series</u>



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_020314, the custom clone sequence may differ by one or more nucleotides

```

ATGGCCGTCCTTCTGGCACTCCAGGAATAGGAACTACAAAGCTGAATTTGCATCATGC
CGACTGGAGGCTGTACCATTTGGAGTTTGGGGACTATCACCTCTGAAACCCATAACTGTC
ACAGAGTCAAAGACAAAGAAAGTGAACCGGAAAGGAAGCACTTCTCCACGTCTCTCTCC
TCCTCCAGCTCCGTGGTGGACCCGCTGAGCAGCGTCTCGATGGGACTGACCCCTCTCC
ATGTTTGCAGCCACTGCTGACCCCGCAGCCTTGGCAGCTGCCATGGACAGCTCCAGAAGG
AAACGTGATAGAGATGATAACTCCGTTGTAGGATCGGATTTTGAGCCTTGACCAACAAA
CGGGGAGAAATCCTTGCCCGGTACACCACTACCGAAAAGCTGTCTATTAATCTGTTTATG
GGATCTGAAAAAGGCAAAGCTGGGACTGCCACATTGGCAATGTCAGAGAAGGTGCGGACC
CGGCTGGAGGAGCTGGATGACTTTGAGGAGGGTCCCAAAGGAGCTGTTGAACTTGACT
CAGCAGGATTACGTGAACCGCATAGAGGAGCTCAACCAATCGCTGAAGGATGCCTGGGCC
TCAGACCAGAAAGTGAAGGCTCTAAAAATAGTCATCCAGTGTTCAAAGCTTCTTTCAGAC
ACCAAGTGTATTACAGTTCTACCAAGCAAATTTGCTTATCACCGACATACTTGATACA
TTTGAAAGCTCGTGACGAGCGCATCTTTCCATGTGTGGATAGCCGAGCGCTCTTA
CCAGATCACTTTTCTCCAGAGAATGCAAATGACACGGCCAAGGAAACATGCCTAAATTGG
TTTTTCAAGATTGCCTCCATCAGGGAACCTATTCCAAGATTTTACGTGGAGGCATCCATC
CTGAAATGTAACAAATTCCTCTCCAAACGGGAATTTAGAGTGCCTGCCCGGTTGACA
TGCATGATCAGAGGGATCGGAGACCCACTAGTGTCCGTGTATGCCCGTGCCTACCTGTGC
CGGGTGGGAATGGAAGTGGCCCCACATCTCAAAGAAACCTAAATAAGAACTTTTTTAC
TTCCTCCTTACGTTCAAACAGATTATGGGGATACGGTCCAGAACAGCTGGTGGTCCAA
GGAGTGGAGCTCCCATCTTACCTCCCTTGTACCCGCTGCCATGGACTGGATCTTCCAG
TGCATCTCCTACCATGCCCCGAGGCTCTGCTGACCGAGATGATGAAAGGTGTAAGAAA
CTAGGAAACAATGCCTTGCTGTTGAATTCTGTGATGTCTGCCTTCCGGGCTGAGTTTCATC
GCCACAAGGTCTATGGATTTTATTGGCATGATTAAAGAGTGTGATGAATCTGGTTTCCCC
AAGCATCTTCTTTTTCGATCACTGGGATTAACCTTGGCCTTGGCTGATCCTCTGAGAGT
GACCGACTTCAGATTCTCAACGAAGCTTGGAAGTCATCACTAAGCTGAAGAACCACAG
GACTACATTAATTGTGCCGAAGTGTGGGTGGAATACACCTGCAAGCATTTACGAAACGA
GAGGTGAATACCGTTTTTGGCAGATGTCATCAAGCACATGACTCCAGATCGTGCATTTGAA
GATTCTACCCCCAGCTTCAGTTAATAATTAAGAAAGTTATTGCCCACTTCCATGACTTC
TCAGTTCTTTTCTCAGTGAAAAATTTCTGCCGTTTCTGGACATGTTCCAAAAAGAGAGT
GTGCGGGTGGAGTTTGCAAATGCATCATGGACGCCTTTATCAAGCATCAACAAGAGCCC
ACCAAGGACCCGGTCATCTTGAATGCCCTTTTGCATGTTTGAAGACCATGCATGACTCT
GTGAATGCACTCACTCTTGAGGATGAGAAAAGAATGCTGTGATATTTGATTAATGGATTT
ATAAAATGGTTTCCTTTGGCCGTGATTTTGAACAACAGCTGAGTTTTTATGTTGAGTCC
AGGTCGATGTTTTGCAATCTGGAGCCTGTTCTTGTGCAGTTGATTCATAGTGTGAACCGG
TTGGCAATGGAGACAAGAAAAGTAATGAAAGGAAATCATTCCAGAAAGACAGCTGCATTT
GTCCGGGCCTGTGTTGCCTACTGCTTCATCACCATCCCCTCCCTGGCGGGCATCTTACA
CGTCTCAATCTCTACCTGCATTCTGGTCAGGTGGCCTTGGCCAACCAAGTGCCTCTCCAA
GCTGATGCTTTTTTCAAAGCCGCTATAAGCCTTGTTCGGGAAGTTCCAAAGATGATTAAT
ATTGATGGGAAGATGCGGCCATCGAATCGTTCCTTCTGGAATTCCTCTGCAATTTCTTT
TCTACTTTATTAATAGTTCGGATCATCCTGAACATGGGGTCTGTTTCTTGTTCGAGAG
CTTCTCAACGTGATCCAGGACTACACCTGGGAGGACAACAGCGATGAGAAAATCCGCATC
TACACCTGCGTCTGCATCTCTCTCCGCCATGAGCCAGGAGACGTACCTTTACCACATA
GACAAAGTGGACTCCAACGACAGCCTCTACGGGGGAGACTCCAAGTTCTTGGCAGAAAAC
AACAAGCTGTGTGAGACGGTGTGGCTCAGATCCTAGAGCATCTGAAAACCTGGCCAAG
GACGAGGCCCTGAAGCGCCAGAGCTCGTTGGGCCTTTCCTTCTTAAACAGCATCTTGGCC
CATGGGGACCTACGCAACAACAAGCTCAACCAGCTCTCCGTCAACCTGTGGCACCTGGCA
CAGAGGCACGGCTGTGCAGACACCAGGACCATGGTGAAAACGCTAGAATACATCAAGAAG
CAAAGCAAACAACCAGACATGACTCATCTGACGGAGCTGGCCCTCAGACTCCCTCTGCAA
ACAAGGACC

```

Restriction Sites:	Please inquire
ACCN:	NM_020314
Insert Size:	3894 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_020314.4</u> , <u>NP_064710.3</u>
RefSeq Size:	3894 bp
RefSeq ORF:	3894 bp
Locus ID:	57020
UniProt ID:	<u>Q7Z3J2</u>
Cytogenetics:	16p12.3
Gene Summary:	Acts as component of the retriever complex. The retriever complex is a heterotrimeric complex related to retromer cargo-selective complex (CSC) and essential for retromer-independent retrieval and recycling of numerous cargos such as integrin alpha-5/beta-1 (ITGA5:ITGB1) (PubMed:28892079). The recruitment of the retriever complex to the endosomal membrane involves CCC and WASH complexes (PubMed:28892079). In the endosomes, drives the retrieval and recycling of NxxY-motif-containing cargo proteins by coupling to SNX17, a cargo essential for the homeostatic maintenance of numerous cell surface proteins associated with processes that include cell migration, cell adhesion, nutrient supply and cell signaling (PubMed:28892079). Involved in copper-dependent ATP7A trafficking between the trans-Golgi network and vesicles in the cell periphery; the function is proposed to depend on its association with the CCC complex and cooperation with the WASH complex on early endosomes. Seems not to be required for CCC complex stability (PubMed:25355947). [UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) encodes the longer isoform (1).