

Product datasheet for **SC318807**

FOLR2 (NM_001113536) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FOLR2 (NM_001113536) Human Untagged Clone
Tag:	Tag Free
Symbol:	FOLR2
Synonyms:	BETA-HFR; FBP; FBP/PL-1; FOLR1; FR-BETA; FR-P3; FRbeta
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001113536, the custom clone sequence may differ by one or more nucleotides ATGGTCTGGAATGGATGCCACTTCTGCTGCTTCTGGTCTGTGTAGCCACCATGTGCAGT GCCCAGGACAGGACTGATCTCCTCAATGTCTGTATGGATGCCAAGCACCACAAGCAAAG CCAGGTCTGAGGACAAGCTGCATGACCAATGCAGTCCCTGGAAGAAGAATGCCTGCTGC ACAGCCAGCACCAGCCAGGAGCTGCACAAGGACACCTCCCGCCTGTACAACCTTAACTGG GACCACTGCGGCAAGATGGAGCCCCTGCAAGCGCCACTTCATCCAGGACACCTGTCTC TATGAGTGCTCACCAACCTGGGGCCCTGGATCCAGCAGGTGAATCAGAGCTGGCGCAA GAACGCTTCTGGATGTGCCCTTATGCAAAGAGGACTGTCAGCGCTGGTGGGAGGATTGT CACACCTCCCACAGTGAAGAGCAACTGGCACAGAGGATGGGACTGGACCTCAGGAGTT AACAAGTGCCAGCTGGGGCTCTCTGCCGACCTTTGAGTCTACTTCCCCTCCAGCT GCCCTTTGTGAAGGCCTCTGGAGTCACTCATAAAGGTGAGCAACTACAGCCGAGGGAGC GGCCGCTGCATCCAGATGTGGTTTGATTGAGCCAGGGCAACCCCAACGAGGAAGTGGCG AGGTTCTATGCTGCAGCCATGCATGTGAATGCTGGTGAGATGCTTCATGGGACTGGGGGT CTCCTGCTCAGTGGCCCTGATGCTGCAACTCTGGCTCCTTGGC
Restriction Sites:	Please inquire
ACCN:	NM_001113536
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001113536.1, NP_001107008.1</u>
RefSeq Size:	1127 bp
RefSeq ORF:	768 bp
Locus ID:	2350
UniProt ID:	<u>P14207</u>
Cytogenetics:	11q13.4
Protein Families:	Druggable Genome, Secreted Protein
Gene Summary:	<p>The protein encoded by this gene is a member of the folate receptor (FOLR) family, and these genes exist in a cluster on chromosome 11. Members of this gene family have a high affinity for folic acid and for several reduced folic acid derivatives, and they mediate delivery of 5-methyltetrahydrofolate to the interior of cells. This protein has a 68% and 79% sequence homology with the FOLR1 and FOLR3 proteins, respectively. Although this protein was originally thought to be specific to placenta, it can also exist in other tissues, and it may play a role in the transport of methotrexate in synovial macrophages in rheumatoid arthritis patients. Multiple transcript variants that encode the same protein have been found for this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (4) differs in the 5' UTR, compared to variant 1. All variants (1-4) encode the same protein.</p>