

Product datasheet for SC318819

ARIH2OS (NM_001123040) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ARIH2OS (NM_001123040) Human Untagged Clone
Tag:	Tag Free
Symbol:	ARIH2OS
Synonyms:	C3orf71
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC318819 representing NM_001123040. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
 ATGCTCGGACAACGTGCCGGCGACGGGGAGCGCCGGGCTCCCGGGCGACGGCGAAGGCGGAGTCCCG
 GCGCGGCCAGGGAGGCGCGCGGAGAGGCCCCCAGCGGCCAGCCAAGGTAAACAAGGCCGTGACGTGC
 GCCGCGCACTTACCGGGAGCTGCGGCCTCGCGGCCGTGAGCCGAACAAGCCAGACCGGGTCAGGCCG
 GGCCAGAGGGACCGATTGGGGCGAAGCGGCAGCGGAGGCGGGCGGCGACGCCGCTCAAGCCCGCGCT
 GCTTCTCCAGAAGAGTCGTCCACAGCTCCGGAAGTGCTTGGCGCCGTTGCGTCACTTCCGGATCGG
 GGTGACCCACGGTCGCTCGGGTCGCGACAGGCTCCCGGCTAGAGGGCCTGTTTAGCGCCGCTCCTTG
 AAATTTAGCGCTCTGACCCAGAGTCTGACCAGGGTACGGCAGGCGCCGACCGCGTCTGGAGCCACTATT
 CGCCTACCAGCGTCTCCCGTCGAGATGTTTTTAACCGCGGTTTCTACCGGCTTCTCATTTCACTGT
 TTGTACTCCGGGATCGGGCACGGAGAAGACATCCTGGCGTCAGTGGAGCAGATAACCATTGTTTCTCGG
 CCGCTATCTGGTCAGAGGGGAGCTGGGCCGGGAAATTCGCTACACCCCGAGGCGGTGCGAGGGTGGC
 CCCAGAGCCGCAACCACACAGGCTTTCGCTTCCCTGCCGAGGCTCGTTCGCCGCGCAGTTCTCCGA
 CTTACGGTCACCGTGCAAGATTGCATCTTAAGTGCCTTGCTTGCAGTTTCTTTTCACAGTATAGGAGTT
 GTCATCATGACTTCCAGTTACCTTCTGGGACCGGTTGTCAAATGA
 ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001123040
Insert Size:	873 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001123040.1</u>
RefSeq Size:	1598 bp
RefSeq ORF:	873 bp
Locus ID:	646450
Cytogenetics:	3p21.31
MW:	30.4 kDa