

Product datasheet for **SC318843**

ESE1 (ELF3) (NM_001114309) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	ESE1 (ELF3) (NM_001114309) Human Untagged Clone
Tag:	Tag Free
Symbol:	ELF3
Synonyms:	EPR-1; ERT; ESE-1; ESX
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC318843 representing NM_001114309. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTGAAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTGCAACCTGTGAGATTAGCAACATTTTTAGCAACTACTTCAGTGCATGTACAGCTCGGAGGAC
TCCACCCTGGCCTCTGTTCCCCCTGCTGCCACCTTTGGGGCCGATGACTTGGTACTGACCCTGAGCAAC
CCCCAGATGTCTGGAGGGTACAGAGAAGGCCAGCTGGTTGGGGGAACAGCCCCAGTTCTGGTCTGAAG
ACGCAGTTCTGGACTGGATCAGCTACCAAGTGGAGAAGAACAAGTACGACGCAAGCGCCATTGACTTC
TCACGATGTGACATGGATGGCGCCACCCTCTGCAATTGTGCCCTTGAGGAGCTGCGTCTGGTCTTTGGG
CCTCTGGGGACCAACTCCATGCCAGCTGCGAGACCTCACTTCCAGCTCTTCTGATGAGCTCAGTTGG
ATCATTGAGCTGCTGGAGAAGGATGGCATGGCCTTCCAGGAGGCCCTAGACCCAGGGCCCTTTGACCAG
GGCAGCCCTTTGCCAGGAGCTGCTGGACGACGGTCAGCAAGCCAGCCCTACCACCCCGCAGCTGT
GGCGCAGGAGCCCTCCCTGGCAGCTCTGACGTCTCCACCGCAGGGACTGGTGCTTCTCGGAGCTCC
CACTCCTCAGACTCCGGTGAAGTGACGTGGACCTGGAATCCCACTGATGGCAAGCTTCCCCAGCGAT
GGTTTTCGTGACTGCAAGAAGGGGGATCCCAAGCACGGGAAGCGGAAACGAGGCCCGCCCGAAAGCTG
AGCAAAGAGTACTGGGACTGTCTCGAGGGCAAGAAGAGCAAGCACGCGCCAGAGGCACCCACCTGTGG
GAGTTCATCCGGGACATCCTCATCCACCCGGAGCTCAACGAGGGCCTCATGAAGTGGGAGAATCGGCAT
GAAGGCGTCTTCAAGTTCCTGCGCTCCGAGGCTGTGGCCCAACTATGGGGCCAAAAGAAAAGAACAGC
AACATGACCTACGAGAAGCTGAGCCGGGCATGAGGTACTACTACAAACGGGAGATCCTGGAACGGGTG
GATGGCCGGCAGCTCGTCTACAAGTTTGGCAAAAACCAAGCGGCTGGAAGGAGGAAGAGGTTCTCCAG
AGTCGGAAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_001114309

Insert Size: 1116 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001114309.1](#)

RefSeq Size: 3263 bp

RefSeq ORF: 1116 bp

Locus ID: 1999

UniProt ID: [P78545](#)

Cytogenetics: 1q32.1

Protein Families: Transcription Factors

MW: 41.5 kDa

Gene Summary: Transcriptional activator that binds and transactivates ETS sequences containing the consensus nucleotide core sequence GGA[AT]. Acts synergistically with POU2F3 to transactivate the SPRR2A promoter and with RUNX1 to transactivate the ANGPT1 promoter. Also transactivates collagenase, CCL20, CLND7, FLG, KRT8, NOS2, PTGS2, SPRR2B, TGFB2 and TGM3 promoters. Represses KRT4 promoter activity. Involved in mediating vascular inflammation. May play an important role in epithelial cell differentiation and tumorigenesis. May be a critical downstream effector of the ERBB2 signaling pathway. May be associated with mammary gland development and involution. Plays an important role in the regulation of transcription with TATA-less promoters in preimplantation embryos, which is essential in preimplantation development (By similarity).[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1 and 2 encode the same protein.