

## Product datasheet for **SC318866**

### TBC1D14 (NM\_001113363) Human Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** TBC1D14 (NM\_001113363) Human Untagged Clone  
**Tag:** Tag Free  
**Symbol:** TBC1D14  
**Vector:** pCMV6 series  
**Fully Sequenced ORF:** >NCBI ORF sequence for NM\_001113363, the custom clone sequence may differ by one or more nucleotides

```
ATGGAATATGAAGACAAGGCTGGAAGACCTAGCAAGCCACCCTCTCCAAAGCAGAATGTG
AGGAAGAATCTTGACTTTGAACCACTTTCCACCACCGCACTCATCCTCGAGGACAGACCA
GCAAATCTCCCAGCAAACCAGCTGAAGAAGCTCAGAAGCACAGACAGCAGTATGAAGAA
ATGGTGGTTTCAGGCCAAAAGCGAGAGCTGAAAGAAGCCAGCGAAGGAAGAAGCAGCTG
GAAGAAAGATGCAGAGTCGAGGAAAGCATTGGAACGCTGTGCTCACCTGGAATAATGAG
ATCTTACCTAACTGGGAAACAATGTGGTGTCTAGAAAAGTTCGAGATTTATGGTGGCAG
GGAATCCCTCAAAGTGTGAGAGGCAAAGTCTGGAGCTTAGCCATTGGCAACGAGTTAAAT
ATCACCCACGAGCTCTTTGACATCTGTCTTGCCCGAGCCAAGGAGAGGTGGCGGTCCCTT
AGCACAGGAGGCTCTGAAGTGGAGAACGAAGATGCTGGTTTTTTCAGCAGCAGACAGAGAA
GCCAGTCTGGAGCTTATTAACCTGGACATTTCTAGAACATTTCTAATCTCTGCATTTTC
CAGCAAGTGGTCCATATCATGACATGTTGCACAGTATTTGGGCGCTTATACTTGTTAC
CGGCCAGATGTGGGTTATGTCCAGGCCATGTCCTTCATAGCAGCAGTGTGATCTTGAAC
TTAGATACTGCAGATGCCTTTATTGCCTTTTCTAACCTTCTGAATAAACCTGTCAAATG
GCGTTTTTTAGAGTGGACCATGGCCTTATGTTGACTTATTTGCTGCATTTGAAGTATTC
TTTGAAGAAAATTTGCCGAAATTTTGGCGATTTCAAGAAGAACAACCTAACTCCAGAT
ATCTACCTAATTGATTGGATCTTTACCTTATATAGTAAATCTCTGCCCTCGACCTGGCC
TGTCGATCTGGGACGTGTTCTGTCGCGATGGGGAAGAGTTCTGTTCCGCACGGCCCTG
GGCATCTGAAGCTGTTGAGGACATCCTGACCAAGATGGACTTCATTACATGGCCAG
TTCTGACCCGGCTGCCCGAGGACCTGCCCGCCGAGGAGCTGTTGCTCCATCGCCACG
ATCCAGATGCAGAGCCGAAACAAGAAGTGGGCTCAGGTACTGACTGCATTGCAGAAAGAC
AGCCGGGAAATGGAGAAGGGAAGTCCGTCCTCCGACAC
```

**Restriction Sites:** Please inquire  
**ACCN:** NM\_001113363  
**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001113363.1</a></u> , <u><a href="#">NP_001106834.1</a></u>
<b>RefSeq Size:</b>	4153 bp
<b>RefSeq ORF:</b>	1242 bp
<b>Locus ID:</b>	57533
<b>UniProt ID:</b>	<u><a href="#">Q9P2M4</a></u>
<b>Cytogenetics:</b>	4p16.1
<b>Gene Summary:</b>	<p>Negative regulator of starvation-induced autophagosome formation.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (3) lacks part of the 5' coding region, contains an alternate 5'-terminal exon, and initiates translation at an alternate start codon, compared to variant 1. The encoded isoform (b) has a shorter N-terminus, compared to isoform a.</p>