

Product datasheet for **SC318878**

HIP55 (DBNL) (NM_001122956) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	HIP55 (DBNL) (NM_001122956) Human Untagged Clone
Tag:	Tag Free
Symbol:	HIP55
Synonyms:	ABP1; HIP-55; HIP55; SH3P7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC318878 representing NM_001122956. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCGCGAACCTGAGCCGGAACGGGCCAGCGTGCAAGAGGCCTACGTGCGGGTGGTACCGAGAAG
TCCCGACCGACTGGGCTCTCTTACCTATGAAGGCAACAGCAATGACATCCGCGTGGCTGGCACAGGG
GAGGGTGGCCTGGAGGAGATGGTGGAGGAGCTCAACAGCGGGAAGGTGATGTACGCTTTCGAGAGTG
AAGGACCCCAACTCTGGACTGCCAAATTTGTCTCATCAACTGGACAGGCGAGGGCGTGAACGATGTG
CGGAAGGGAGCCTGTGCCAGCCACGTGACACCATGGCCAGCTTCTGAAGGGGGCCCATGTGACCATC
AACGCACGGGCCGAGGAGGATGTGGAGCCTGAGTGCATCATGGAGAAGGTGGCCAAGGCTTCAGGTGCC
AACTACAGCTTTCACAAGGAGAGTGGCCGCTTCCAGGACGTGGGACCCCAAGGCCCAAGTGGGCTCTGTG
TACCAGAAGACCAATGCCGTGTCTGAGATTAAGGGTTGGTAAAGACAGCTTCTGGGCCAAAGCAGAG
AAGGAGGAGGAGAACCGTCGGCTGGAGGAAAAGCGGCGGGCCGAGGAGGCACAGCGGCAGCTGGAGCAG
GAGCGCCGGGAGCGTGAGCTGCGTGAGGCTGCACGCCGGGAGCAGCGCTATCAGGAGCAGGGTGGCGAG
GCCAGCCCCAGAGCAGGACGTGGGAGCAGCAGCAAGAAGTGGTTTCAAGGAACCGAAATGAGCAGGGG
TCAACATGTGCTTCCCTCCAGGAGTCTGCCGTGCACCCGAGGGAGATTTTCAAGCAGAAGGAGAGGGCC
ATGTCCACCACCTCCATCTCCAGTCTCAGCCTGGCAAGCTGAGGAGCCCTTCTGCAGAAGCAGCTC
ACCCAACAGAGACCCACTTTGGCAGAGAGCCAGCTGCTGCCATCTCAAGGCCAGGGCAGATCTCCCT
GCTGAGGAGCCGCGCCAGCACTCCTCATGTCTGGTGCAGGCAGAAGAGGAGGCTGTGTATGAGGAA
CCTCCAGAGCAGGAGACCTTCTACGAGCAGCCCCACTGGTGCAGCAGCAAGGTGCTGGCTCTGAGCAC
ATTGACCACACATTCAGGGCCAGGGGCTCAGTGGGCAAGGGCTCTGTGCCCTGCCCTGTACGACTAC
CAGGCAGCCGACGACACAGAGATCTCCTTTGACCCCGAGAACCCTCATCACGGGCATCGAGGTGATCGAC
GAAGGCTGGTGGCGTGGCTATGGCCGGATGGCCATTTTGGCATGTTCCCTGCCAACTACGTGGAGCTC
ATTGAGTGA
ACGGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001122956
Insert Size:	1320 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001122956.1
RefSeq Size:	2237 bp
RefSeq ORF:	1320 bp
Locus ID:	28988
UniProt ID:	Q9UJU6
Cytogenetics:	7p13
MW:	49 kDa
Gene Summary:	<p>Adapter protein that binds F-actin and DNM1, and thereby plays a role in receptor-mediated endocytosis. Plays a role in the reorganization of the actin cytoskeleton, formation of cell projections, such as neurites, in neuron morphogenesis and synapse formation via its interaction with WASL and COBL. Does not bind G-actin and promote actin polymerization by itself. Required for the formation of organized podosome rosettes (By similarity). May act as a common effector of antigen receptor-signaling pathways in leukocytes. Acts as a key component of the immunological synapse that regulates T-cell activation by bridging TCRs and the actin cytoskeleton to gene activation and endocytic processes.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (3) represents the longest transcript and encodes the longest isoform (c).</p>